Porphyrine HPLC Kit

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

Porphyrins HPLC Kit

For the determination of Porphyrins (free acids) in urine

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1. INTENDED USE	14
2. SUMMARY AND EXPLANATION OF THE TEST	14
3. PRINCIPLE OF THE TEST	15
4. MATERIAL SUPPLIED	16
5. MATERIAL REQUIRED BUT NOT SUPPLIED	16
6. PREPARATION AND STORAGE OF REAGENTS	16
7. PRECAUTIONS	17
8. SPECIMEN COLLECTION AND PREPARATION	17
9. ASSAY PROCEDURE	17
PROCEDURAL NOTES	17
SAMPLE AND STANDARD PREPARATION	18
CHROMATOGRAPHIC CONDITIONS	18
10. TREATMENT OF THE COLUMN	18
11. RESULTS	19
CALCULATION	19
TYPICAL CHROMATOGRAM	19
12. LIMITATIONS	20
13. QUALITY CONTROL	20
EXPECTED VALUES	20
CONTROLS	20
14. PERFORMANCE CHARACTERISTICS	21
PRECISION AND REPRODUCIBILITY	21
LINEARITY	21
DETECTION LIMIT	21
15. DISPOSAL	21
16. TROUBLESHOOTING	22
17. REFERENCES	23
18. GENERAL NOTES ON THE TEST AND TEST PROCEDURE	23

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
Porphyryns HPLC Kit

For the determination of Porphyryns (free acids) in urine

Valid from 14.06.2007

REF

KC 2601

CAL
CTRL 1
CTRL 2

CE

IVD

1. INTENDED USE

The *Immundiagnostik* Assay is intended for the quantitative determination of porphyrins in urine. This Assay is designed for *in vitro* diagnostic use only.

2. SUMMARY AND EXPLANATION OF THE TEST

The porphyrins are precursors of heme, which is a cofactor in hemoglobin, myoglobin and the cytochromes. They play an important role in the oxygen metabolism. Heme biosynthetic activity is located quantitatively in bone marrow and liver. The series of reactions leading to heme synthesis begins with the condensation of succinyl coenzyme A (CoA) and glycine and ends with the insertion of an iron atom into a molecule of protoporphyrin IX.

There are several genetic defects of heme biosynthesis described. In all of them an increase of porphyrins is recognized.

The laboratory diagnosis and classification of porphyrin disorders are based primarily on the measurement of the excretion of porphyrins and porphyrin precursors, mainly in urine. Uro-, copro-, penta- and tricarboxyporphyrin are increased in autosomal dominant acute porphyrin disorders. In chronic porphyrin disorders, including *Porphyria cutanea tarda*, uro- and heptaporphyrin are elevated. In chronic lead intoxication the coproporphyrins are slightly increased (0.5-2 $\mu\text{mol}/24$ h urine), acute lead poisoning is characterized by extremely high excretion of total porphyrins (up to 15 $\mu\text{mol}/24$ h urine).

Indications

- Hereditary hepatic porphyria (*porphyria variegata*)
- *Porphyria cutanea tarda* in chronic liver disease and alcoholic-liver-syndrom
- Hepatic porphyria in prostate and liver tumors
- Acute and chronic lead poisoning

3. PRINCIPLE OF THE TEST

The first step in the determination of Porphyrins includes an easy sample preparation. The pH value of the sample and the calibrator is adjusted below 2.5 by addition of 20 µl of 25% HCl. After centrifugation, the supernatant is injected into the HPLC system.

The separation *via* HPLC follows a gradient method, using a „reversed phase“ column at 30 °C; one run lasts 25 minutes. The quantification is performed with the delivered calibrator. The concentration is calculated *via* integration of the peak areas or peak heights.

Summary

The HPLC technique provides an easy, fast and precise method for quantitative determination of porphyrins. The kit contains all reagents necessary for sample preparation and separation in ready-to-use form except the column.

Like in the case of many other parameters, the advantage of the HPLC analysis is the simultaneous handling of many components in a single test. A complete HPLC system allows even laboratories without any experience in „high performance liquid chromatography“ a quick and precise use of the technique for clinical chemistry routine analyses. Mostly, a one-point calibration of the test system is sufficient - unlike immuno assays with up to 6 calibrators per test. With short test series, the one-point calibration is much more economic than the 6-point calibration for immuno assays. In addition, it is possible to automate the sample application and the result calculation, so that even a higher number of samples can be handled nearly without control.

4. MATERIAL SUPPLIED

Cat. No	Content	Kit Components	Quantity
KC 2601LA	MOPHAA	Eluent A	1000 ml
KC 2601LB	MOPHAB	Eluent B	2 x 1000 ml
KC 2601KA	CAL	Calibrator, lyophilized 1 ml	10 vials
KC 2601HC	HCL	HCl (25%)	10 ml
KC 2601RE	RECSOL	Reconstitution solution	10 ml
KC 2601KO	CTRL1 CTRL2	Control 1 and 2; 1 ml lyophilized (concentration, see product data sheet)	2 x 3 vials

HPLC column (KC 2601RP) as well as individual components can be ordered separately from Immundiagnostik. Please ask for the price list of the individual components.

5. MATERIAL REQUIRED BUT NOT SUPPLIED

- 1.5 ml reaction tubes (Eppendorf)
- Centrifuge
- Various pipettes
- HPLC gradient pump with Fluorescence-detector
- Reversed phase C₁₈-column
- Vortex mixer

6. PREPARATION AND STORAGE OF REAGENTS

- Reconstitute the **calibrator** (CAL) and **controls** (CTRL1, CTRL2) in **1 ml** reconstitution solution (RECSOL). The concentration of the porphyrins might have minor changes from lot to lot. The concentrations of the different porphyrins are given in the attached product specification sheet. After reconstitution the calibrator should not be stored longer than 24 h at 4°C.
- Store **Calibrator** (CAL) and **controls** (CTRL1, CTRL2) at – 20 °C. All other test reagents are stable at 20-25 °C, up to the date of expiry stated on the label.

7. PRECAUTIONS

- For *in vitro* diagnostic use only.
- The 25% HCl (HCL) solution contains acid. Although diluted, it still must be handled with care. It can cause burns and should be handled with gloves, eye protection, and appropriate protective clothing. Any spill should be wiped out immediately with copious quantities of water. Do not breathe vapor and avoid inhalation.
- Reagents should not be used beyond the expiration date shown on kit label.

8. SPECIMEN COLLECTION AND PREPARATION

Urine is suited for this test system.

The porphyrins are light- and temperature sensitive. Therefore, protect samples from light and cool immediately after collection.

9. ASSAY PROCEDURE

Procedural notes

- Quality control guidelines should be observed.
- Incubation time, incubation temperature and pipetting volumes of the components are defined by the producer. Any variation of the test procedure, which is not coordinated with the producer, may influence the results of the test. Immundiagnostik AG can therefore not be held responsible for any damage resulting from wrong use.
- The assay should always be performed according the enclosed manual.

Sample and standard preparation

- Mix **1 ml** of a cool and dark stored 24-h-urine (without additives) or calibrator (CAL) with 20 μ l 25% HCl (HCL) (pH should be below 2,5).
- Centrifuge for 5 min at 10.000 x g.
- Inject **100 μ l** in the HPLC-system.

Chromatographic conditions

Column material:	Bischoff ProntoSil 120-5-C18 ace EPS; 5 μ m
Column dimension:	125 mm x 4 mm
Flow rate:	0,75 ml/min
Fluorescence detection:	Excitation: 400 nm Emission: 620 nm Sensitivity high
Temperature:	30 °C
Injection volume	100 μ l
Run time	25 min
Gradient:	0 min (38% B / 62 % A) 1,5 min (38% B / 62 % A) 9,5 min (80% B / 20 % A) 11,5 min (80% B / 20 % A) 11,6 min (95% B / 5 % A) 15 min (95% B / 5 % A) 15,1 min (38% B / 62 % A) 24 min (38% B / 62 % A)

10. TREATMENT OF THE COLUMN

After the analysis, the column should be flushed with 30 ml aqua bidest. (1 ml/min) and stored in 50% methanol in aqua bidest (ca. 30 ml, flow 0.7 ml/min). Before use, the system should be equilibrated with ca. 30 ml eluent.


11. RESULTS

Calculation

Samples are evaluated with the supplied calibrator by comparison of the peak areas.

$$\text{Concentration sample} = \frac{\text{Peak height sample} \times \text{Concentration of the calibrator}}{\text{Peak height calibrator}}$$

Typical chromatogram

Column: ProntoSil 120-5-C18 ace-EPS

Order of peaks (from the left to the right)

1. 8-carboxyl Porphyrin (Uro-)
2. 7-carboxyl Porphyrin (Hepta-)
3. 6-carboxyl Porphyrin (Hexa-)
4. 5-carboxyl Porphyrin (Penta)
5. 4-carboxyl Porphyrin (Copro) I
6. 4-carboxyl Porphyrin (Copro) III
7. Mesoporphyrin IX

12. LIMITATIONS

Do not use blood samples for analysis.

13. QUALITY CONTROL

Expected values

The normal range is related to 24 h urine.

	Mean ($\mu\text{g}/24\text{ h}$)	Area ($\mu\text{g}/24\text{ h}$)
8-carboxyl Porphyrin (Uro-)	7.0	0-33
7-carboxyl Porphyrin (Hepta-)	< detection limit	0-10
6-carboxyl Porphyrin (Hexa-)	< detection limit	0-7
5-carboxyl Porphyrin (Penta)	< detection limit	0-5
4-carboxyl Porphyrin (Copro)	44.8	0-120
Coprop. I 17 – 31 % (Ref. 3)		
Coprop. III 69 – 83 % (Ref. 3)		

It is recommended that each laboratory should establish its own normal range. Above mentioned values are only for orientation and may vary from other published data.

Controls

Control samples should be analyzed with each run. Results, generated from the analysis of control samples, should be evaluated for acceptability using appropriate statistical methods. The results for the patient samples may not be valid, if within the same assay one or more values of the quality control sample are outside the acceptable limits.

14. PERFORMANCE CHARACTERISTICS

Precision and reproducibility

Intra-Assay CV:

Uroporphyrin	1,9% (104,6µg/l)	[n=10]
Heptaporphyrin	3,8% (44,0µg/l)	[n=10]
Hexaporphyrin	3,7% (46,5µg/l)	[n=10]
Pentaporphyrin	3,5% (58,4µg/l)	[n=10]
Coproporphyrin I	3,3% (68,1µg/l)	[n=10]
Corpoporphyrin III	3,8% (58,26µg/l)	[n=10]

Inter-Assay VK:

Uroporphyrin	6,6% (52,5µg/l)	[n=5]
Heptaporphyrin	6,0% (30,1µg/l)	[n=5]
Hexaporphyrin	5,0% (30,8µg/l)	[n=5]
Pentaporphyrin	5,1% (35,9µg/l)	[n=5]
Coproporphyrin I	7,5% (24,0µg/l)	[n=5]
Corpoporphyrin III	4,7% (65,9µg/l)	[n=5]

Linearity

up to 1 mg/l

Detection limit

2,5 µg/l

15. DISPOSAL

The mobile phase (MOPHA) must be disposed as non-halogenated solvent. The 25% HCl solution could be neutralized with NaOH and if the pH value is neutral it can be disposed as salt solution.

Important: Reaction will produce heat, be carefull!. Please refer to the appropriate national guidelines.

16. TROUBLESHOOTING

Problem	Possible reason	Solution
No signal	No or defect connection to evaluation system	Check signal cord and connection
	Detector lamp is altered	Change lamp
No peaks	Injector is congested	Check Injector
Double peaks	Dead volume in fittings and / or column	Renew fittings and / or column
Contaminating peaks	Injector dirty	Clean injector
	Contamination at the head of the column	Change direction of the column and rinse for 30 min at low flow rate (0.2 ml/min) with mobile phase
	Air in the system	Degas pump
	Auto sampler vials contaminated	Use new vials or clean them with methanol
Broad peaks, tailing	Precolumn / column exhausted	Use new precolumn / column
Variable retention times	Drift in temperature	Use a column oven
	Pump delivers imprecise	Check pump, degas the system
	System is not in steady state yet	Rinse system mobile phase for 15 min
Baseline is drifting	Detector lamp did not reach working temperature yet	Wait
	Detector lamp is too old	Renew lamp
	System is not in steady state yet	Rinse system mobile phase for 15 min
	Pump delivers imprecise	Check pump, degas the system
Baseline is not smooth	Pump delivers imprecise	Check pump, degas the system
	Detector flow cell is dirty	Clean flow cell

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18. GENERAL NOTES ON THE TEST AND TEST PROCEDURE

- This assay was produced and put on the market according to the IVD guidelines of 98/79/EC.
- The test components contain organic solvents. Contact with skin or mucous membranes has to be avoided.
- All reagents in the test package are for research use only.
- The reagents should not be used after the date of expiry stated on the label.
- Do not interchange different lot numbers of any kit component within the same assay.
- Quality control guidelines should be observed.
- Incubation time, incubation temperature and pipetting volumes of the components are defined by the producer. Any variation of the test procedure, which is not coordinated with the producer, may influence the results of the test. Immundiagnostik AG can therefore not be held responsible for any damage resulting from wrong use.