# Vitamin B<sub>6</sub> HPLC Kit

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Kalibrator


# Vitamin B<sub>6</sub> HPLC Kit

For the determination of Vitamin B<sub>6</sub> in plasma, serum and  
EDTA-whole blood

Valid from 02.03.2010

KC 2150

CAL
CTRL 1
CTRL 2

**Immundiagnostik AG**  
Stubenwald-Allee 8a  
64625 Bensheim, Germany

# Content

<b>1. INTENDED USE</b> .....	<b>15</b>
<b>2. SUMMARY AND EXPLANATION OF THE TEST</b> .....	<b>15</b>
<b>3. PRINCIPLE OF THE TEST</b> .....	<b>15</b>
<b>4. MATERIAL SUPPLIED</b> .....	<b>16</b>
<b>5. MATERIAL REQUIRED BUT NOT SUPPLIED</b> .....	<b>17</b>
<b>6. PREPARATION AND STORAGE OF REAGENTS</b> .....	<b>17</b>
<b>7. PRECAUTIONS</b> .....	<b>17</b>
<b>8. SPECIMEN COLLECTION AND PREPARATION</b> .....	<b>18</b>
<b>9. ASSAY PROCEDURE</b> .....	<b>18</b>
Procedural notes .....	18
Sample and standard preparation.....	19
Chromatographic conditions .....	19
<b>10. TREATMENT OF THE COLUMN</b> .....	<b>20</b>
<b>11. RESULTS</b> .....	<b>20</b>
Calculation .....	20
Typical chromatogram .....	20
<b>12. QUALITY CONTROL</b> .....	<b>21</b>
Expected values .....	21
Controls .....	21
<b>13. PERFORMANCE CHARACTERISTICS</b> .....	<b>21</b>
Precision and reproducibility .....	21
Linearity.....	21
Detection limit: .....	22
Recovery: .....	22
<b>14. DISPOSAL</b> .....	<b>22</b>
<b>15. TROUBLESHOOTING</b> .....	<b>22</b>
<b>16. REFERENCES</b> .....	<b>23</b>
<b>17. GENERAL NOTES ON THE TEST AND TEST PROCEDURE</b> .....	<b>24</b>

## 1. INTENDED USE

The Immundiagnostik Assay is intended for the quantitative determination of vitamin B<sub>6</sub> in plasma, serum and EDTA-whole blood. This Assay is designed for in vitro diagnostic use only.

## 2. SUMMARY AND EXPLANATION OF THE TEST

Vitamin B<sub>6</sub> is the term for pyridoxin, pyridoxal and pyridoxamin and the appropriate phosphate products. All forms can be transformed into the active form pyridoxal-5-phosphate (here described as „vitamin B<sub>6</sub>“).

Vitamin B<sub>6</sub> functions as a coenzyme and is essential for more than 50 reactions in the protein metabolism thereby synthesizing, transforming or degrading amino acids. Vitamin B<sub>6</sub> supports the resorption of amino acids and their transport into the cells. Furtheron vitamin B<sub>6</sub> contributes to the synthesis of neuro transmitters and amine products (histamin).

Due to the fact that vitamin B<sub>6</sub> contributes to a variety of different reactions lack of vitamin B<sub>6</sub> results in various clinical pictures.

Applications:

- determination of vitamin B<sub>6</sub> status
- homocysteinaemia
- skin diseases
- movement disorders
- Anaemia, depression

## 3. PRINCIPLE OF THE TEST

The first step in the determination of Vitamin B<sub>6</sub> includes the sample preparation with additional derivatisation. During the precipitation higher molecular substances are removed. After centrifugation the supernatant is used for derivatisation (20 min at 60°C) thereby transforming the vitamin B<sub>6</sub> into a fluorescent product. The sample is cooled, centrifuged and injected into the HPLC system.

The separation via HPLC follows an isocratic method at 30°C using a „reversed phase“ column; one run lasts 10 minutes. The quantification is performed with the delivered plasma calibrator; the concentration is calculated via integration of the peak heights.

### Summary

The application of vitamin B<sub>6</sub> for HPLC makes it possible to determine the vitamin in an easy, fast and precise method. The kits includes all reagents in ready to use form for preparation and separation of the samples with exception of the columns.

As many other parameters the advantage of HPLC analytics is the simultaneous handling of many analytes in a single test. The HPLC complete system enables even laboratories without experience in „high performance liquid chromatography“ to use this technique for clinical chemical routines quickly and precisely. Mostly a one-point calibration is sufficient for calibrating the test system - unlike immuno assays with up to 6 calibrators per test. It is possible to automate the sample application and calculation of the results so that even higher number of samples can be handled nearly without control. (With short test series the one-point calibration is much more economic than 6-point calibration for immuno assays).

## 4. MATERIAL SUPPLIED

Cat.No	Label	Kit Components	Quantity
KC2100LM	MOPHA	Mobile phase	1000 ml
KC2100KA	CAL	Calibrator, 4 ml lyophilized	1 vial
KC2150FR	PREC	Precipitating reagent	30 ml
KC2100RE	RECSOL	Reconstitution solution	10 ml
KC2100DL	DER	Derivatization solution)	3 x 8,5 ml
KC2100KO	CTRL 1 CTRL 2	Control 1 and 2; 250 µl lyophilized	2 x 3 vials

HPLC column (KC 2100 RP) as well as individual components can be ordered separately from Immundiagnostik. Please ask for the price list of the individual components.

## 5. MATERIAL REQUIRED BUT NOT SUPPLIED

- 1.5 ml reaction tubes (Eppendorf)
- Centrifuge
- Various pipettes
- HPLC with Fluorescence-detector
- Reversed phase C18-column
- Heatable shaker or water bath
- Vortex mixer

## 6. PREPARATION AND STORAGE OF REAGENTS

- Reconstitute the **CAL** (calibrator) in **4 ml RECSOL** (reconstitution solution). Take aliquots and store them at -20°C. Avoid repeated freeze-thaw cycles. The concentration of vitamin B<sub>6</sub> might have minor changes from lot to lot.
- Reconstitute **CTRL 1, CTRL 2** (control 1 and 2) in **250 µl RECSOL** (reconstitution solution).
- **All other test reagents are stable at 2-8 °C**, the **CAL** (calibrator) and **CTRL 1, CTRL 2** (control 1 and 2) **at -20 °C**, up to the date of expiry stated on the label.

## 7. PRECAUTIONS

- Human materials used in kit components were tested for hepatitis-B (HBsAg), HCV, HIV-1 and HIV-2 antibodies and found to be negative. However, for safety reasons, all kit components should be treated as potentially infectious.
- The **PREC** (precipitating reagent) contains acid. Even diluted, it still must be handled with care. It can cause acid burns and should be handled with gloves, eye protection, and appropriate protective clothing. Any spills should be wiped out immediately with copious quantities of water.
- Reagents should not be used beyond the expiration date shown on the kit label.
- The **DER** (derivatization solution) contains KCN and should be pipetted in a fume hood.

## 8. SPECIMEN COLLECTION AND PREPARATION

EDTA-plasma and serum or EDTA-whole blood are suitable for use in the assay.

Vitamin B<sub>6</sub> is light- and temperature sensitive; therefore samples have to be protected from light and cooled and centrifuged immediately.

The samples are stable in the dark at 2-8°C for 1 week. For longer storage samples should be frozen at -20°C.

## 9. ASSAY PROCEDURE

### Procedural notes

- Quality control guidelines should be observed.
- Incubation time, incubation temperature and pipetting volumes of the components are defined by the producer. Any variation of the test procedure, which is not coordinated with the producer, may influence the results of the test. Immundiagnostik AG can therefore not be held responsible for any damage resulting from wrong use.
- The assay should always be performed according the enclosed manual.

## Sample and standard preparation

Pipette into 1.5 ml reaction tubes:	
1.	<b>200 µl CAL</b> (calibrator), serum sample or <b>CTRL 1</b> and <b>CTRL 2</b> (control 1 and 2) +
2.	<b>300 µl PREC</b> (precipitating reagent) Vortex for at least 30 seconds.
3.	Mix well. Leave the tubes for <b>10 minutes at 2-8°C</b> and centrifuge afterwards at <b>10.000g for 2 minutes</b> .
4.	<b>250 µl supernatant</b> +
5.	<b>250 µl DER</b> (derivatization solution)
6.	Mix well. Incubate for <b>20 minutes at 60°C on a shaker</b> or in a waterbath; cool to <b>2-8°C</b> and centrifuge at <b>10.000g for 5 minutes</b> . Take supernatant. It is stable <b>in the dark for 5 days at 2-8°C</b> .
7.	Inject <b>50 µl of the supernatant</b> into the HPLC-system.

## Chromatographic conditions

<b>Column material:</b>	Bischoff Prontosil Eurobond, 5 µm Merck Lichrospher 5 µm MZ- Inertsil ODS-2; 5 µm		
<b>Column dimension:</b>	125 mm x 4 mm		
<b>Flow rate:</b>	1-1.5 ml/min		
<b>Detection :</b>	Fluorescence:	Excitation	320 nm
		Emission	415 nm
<b>Injection volume:</b>	20 µl		
<b>Running time:</b>	7 min (Dialysis patients 15 minutes)		
<b>Temperature:</b>	30 °C		

## 10. TREATMENT OF THE COLUMN


After the analysis the column should be flushed with 30 ml aqua bidest (1 ml/min) and stored in 50% methanol in aqua bidest (ca. 30 ml, flow 0.5 ml/min). Before use, the system should be equilibrated with ca. 30 ml Eluent.

## 11. RESULTS

### Calculation

$$\text{Concentration sample} = \frac{\text{Peak height sample} \times \text{Concentration of the calibrator}}{\text{Peak height calibrator}}$$

### Typical chromatogram

## 12. QUALITY CONTROL

### Expected values

Serum and Plasma

**4,3 - 17,5 ng/ml (n = 90)**

Whole blood

**8,7 - 27,2 ng/ml\***

\*Reference: den Ottolander, Diagnostisch Kompas, 1997 Ziekenfondsraad, Amstelveen, Holland

We recommend, that each laboratory should develop their own normal range. The values mentioned above are only for orientation and can deviate from other published data.

Control samples or serum pools should be analyzed with each run of calibrators and patient samples. Results generated from the analysis of the control samples should be evaluated for acceptability using appropriate statistical methods. In assays in which one or more of the quality control sample values lie outside the acceptable limits, the results for the patient sample may not be valid.

### Controls

Control samples or serum pools should be analyzed with each run of calibrators and patient samples. Results generated from the analysis of the control samples should be evaluated for acceptability using appropriate statistical methods. In assays in which one or more of the quality control sample values lie outside the acceptable limits, the results for the patient sample may not be valid.

## 13. PERFORMANCE CHARACTERISTICS

### Precision and reproducibility

**Intra-Assay CV:** 2.3 % (28.4 ng/ml) [n = 12]

**Inter-Assay CV:** 5.4 % (12.9 ng/ml) [n = 12]

### Linearity:

up to 250 ng/ml

Detection limit:

1.04 ng/ml

Recovery:

84,1 %

## 14. DISPOSAL

The DER (derivatisation solution) must be oxidized with hydrogen peroxide, the pH value adjusted to 6-8, and disposed as aqueous salt solution. The MOPHA (mobile phase) and the PREC (precipitation solution) must be neutralized with NaOH to neutral pH and disposed as salt solution.

**Important:** Reaction will produce heat, be careful! Please refer to the appropriate national guidelines.

## 15. TROUBLESHOOTING

Problem	Possible reasons	Solution
No signal	No or defect connection to evaluation system.	Check signal cord and connection.
	Detector lamp is altered	Change lamp
No peaks	Injector is congested	Check Injector
Double peaks	Dead volume in fittings and / or column	Renew fittings and / or column
Contaminating	Injector dirty	Clean injector
peaks	Contamination at the head of the column	Change direction of the column and rinse for 30 min at low flow rate (0.2 ml/min) with mobile phase
	Air in the system	Degas pump
	Autosampler vials contaminated	Use new vials or clean them with methanol

Problem	Possible reasons	Solution
Broad peaks, tailing	Precolumn / column exhausted	Use new precolumn / column
Variable retention times	Drift in temperature	Use a column oven
	Pump delivers imprecise	Check pump, degas the system
	System is not in steady state yet	Rinse system mobile phase for 15 min
Baseline is drifting	Detector lamp did not reach working temperature yet	Wait
	Detector lamp is too old	Renew lamp
	System is not in steady state yet	Rinse system mobile phase for 15 min
	Pump delivers imprecise	Check pump, degas the system
Baseline is not calm	Pump delivers imprecise	Check pump, degas the system
	Detector flow cell is dirty	Clean flow cell

## 16. REFERENCES

Ambrosch A. et al. (2000). Relation between homocysteinaemia and diabetic neuropathy in patients with Type 2 diabetes mellitus. *Diabetic Med* 18; 185-192.

Dierkes J. et al. (2001). Vitamin supplementation can markedly reduce the homocysteine elevation induced by fenofibrate. *Atherosclerosis* 158; 161-164.

Dierkes et al. (2001). Homocysteine lowering effect of different multivitamin preparations in patients with end-stage renal disease. *J Renal Nut* 11; 67-72.

## 17. GENERAL NOTES ON THE TEST AND TEST PROCEDURE

- This assay was produced and put on the market according to the IVD guidelines of 98/79/EC.
- The test components contain organic solvents. Contact with skin or mucous membranes has to be avoided.
- All reagents in the kit package are for in vitro diagnostic use only.
- Reagents should not be used beyond the expiration date shown on the kit label.
- Do not interchange different lot numbers of any kit component within the same assay.
- Guidelines for medical laboratories should be observed.
- Incubation time, incubation temperature and pipetting volumes of the components are defined by the producer. Any variation of the test procedure, which is not coordinated with the producer, may influence the results of the test. Immundiagnostik AG can therefore not be held responsible for any damage resulting from wrong use.

### Used Symbols:

	Store at		Catalog Number
	In Vitro Diagnostic Device		No. of tests
	Manufacturer		Use by
	Chargennummer		Controls
	Calibrator		

**Immundiagnostik AG**

Stubenwald-Allee 8a  
D-64625 Bensheim

Tel.: +49(0) 62 51/70 19 00

Fax: +49(0) 62 51/84 94 30

[info@immundiagnostik.com](mailto:info@immundiagnostik.com)

[www.immundiagnostik.com](http://www.immundiagnostik.com)