25-OH-Vitamin D direct ELISA Kit

Zur in-vitro-Bestimmung von 25-OH-Vitamin D
in humanem Serum

25-OH-Vitamin D direct ELISA Kit

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Hersteller

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
25-OH-Vitamin D direct ELISA Kit

For the determination of 25-OH-Vitamin D
in human serum

Valid from 12.02.2009

K 2109

EU patent #EP1097132

Australian patent # 763458

Content

1. INTENDED USE	17
2. SUMMARY AND EXPLANATION OF THE TEST	17
3. MATERIAL SUPPLIED	18
4. MATERIAL REQUIRED BUT NOT SUPPLIED	18
5. PREPARATION AND STORAGE OF REAGENTS	19
6. SPECIMEN COLLECTION AND PREPARATION	19
7. ASSAY PROCEDURE	20
Principle of the test	20
Test procedure	20
8. RESULTS	22
9. LIMITATIONS	22
10. QUALITY CONTROL	22
Expected values	22
11. PERFORMANCE CHARACTERISTICS	23
Precision and reproducibility	23
Specificity	24
Sensitivity	24
Linearity	24
Correlation data	25
12. PRECAUTIONS	26
13. TECHNICAL HINTS	27
14. GENERAL NOTES ON THE TEST AND TEST PROCEDURE	27
15. REFERENCES	28

1. INTENDED USE

The Immundiagnostik ELISA is intended for the quantitative determination of the 25-OH-Vitamin D in serum. For in vitro diagnostic use only.

2. SUMMARY AND EXPLANATION OF THE TEST

Vitamin D is a steroid hormone involved in the intestinal absorption of calcium and the regulation of calcium homeostasis. There are two different forms of Vitamin D, named D₃ and D₂, which are very similar in structure. The D₂ is a synthetic product, which is predominantly absorbed by fortified food.

Physiological Vitamin D₃ levels result not only from dietary uptake but can also be produced from a cholesterol precursor, 7-dehydrocholesterol, in the skin during sun exposure. In the liver, the vitamin is hydroxylated to 25-hydroxyvitamin D (25-OH-Vitamin D), the major circulating metabolite of Vitamin D. Although 1,25-(OH)₂ Vitamin D portrays the biological active form of Vitamin D, which is synthesized in the kidney, it is widely accepted that the measurement of circulating 25-OH-Vitamin D provides better information with respect to patients Vitamin D status and allows its use in diagnose hypovitaminosis (1,2).

The concentration of 25-OH-Vitamin D decreases with age and a deficiency is common among elderly persons.

Clinical applications of 25-OH-Vitamin D measurements are the diagnosis and therapy control of postmenopausal osteoporosis, rickets, osteomalacia, renal osteodystrophy, pregnancy, neonatal hypocalcemia and hyperparathyroidism. In addition, a prevalence of subclinical Vitamin D deficiency has been discussed in different European countries.

Vitamin D intoxication mostly occurs during a large intake of pharmaceutical preparations of Vitamin D and may lead to hypercalcemia, hypercalcuria and nephrocalcinosis in susceptible infants.

Frequency of vitamin D deficiency

Mild vitamin D deficiency caused by too low exposure to UV-B rays (sunlight deficiency) or decreased ability of cutaneous vitamin D synthesis (older persons) is very common in Europe. In Germany and France, approximately one third of the population over 50 years of age suffer from vitamin D deficiency during winter. Immigrants with dark skin present a special risk group. Serious vitamin D deficiencies can be detected in 3% of all geriatric patients. Vitamin D deficiencies can be observed even in young persons with gastrointestinal illnesses (liver function defects, malabsorption) or accelerated metabolism (from drugs such as antiepileptics).

In Germany, low values are very often detected in the general population. However, most of the cases are associated with a pre-existing secondary hyperparathyroidism. 25-hydroxyvitamin D values in the range 12-30 ng/ml, 30-75 nmol/l respectively, are considered as "decreased values" by the most authors (sub-clinical vitamin D deficiency).

3. MATERIAL SUPPLIED

Catalogue No.	Content	Kit Components	Quantity
K 2109MTP	PLATE	One holder with precoated strips	12 x 8 wells
K 2109WP	WASHBUF	ELISA wash concentrate 20x	50 mL
K 2109RS	RECSOL	Reconstitution solution	20 mL
K 2109RR	RELREAG	Releasing reagent	1 vial
K 2109AK	AB	Anti 25-OH-Vitamin D3 antibody ready to use	18 mL
K 2109ST	STD	Standards, ready to use	6 vials 300 µL each
K 2109KO	CTRL	Controls, ready for use (see specification for range)	2 vials 300 µL each
K 2109K	CONJ	Conjugate, peroxidase labeled, ready to use	22 mL
K 2109TMB	SUB	TMB substrate (Tetramethylbenzidine)	2 x 15 mL
K 2109 SD	SAMDIL	Sample Dilution Buffer	35 mL
K 2109AC	STOP	ELISA stop solution, ready to use	1 x 15 mL

4. MATERIAL REQUIRED BUT NOT SUPPLIED

- Bidistilled water (aqua bidest.)
- Deep freezer -20 °C
- Precision pipettors calibrated to deliver 10–1000 µl
- Horizontal microtiter plate shaker
- Multi-channel dispenser or repeating dispenser
- Vortex-Mixer
- Water bath or heating block
- Standard laboratory glass or plastic vials, cups, etc. (one time products)
- Microtiter plate reader 450 nm (reference wave length 620 or 690 nm)
- Refridgerator with defined 8–10 °C

5. PREPARATION AND STORAGE OF REAGENTS

- To run assay more than once, ensure that reagents are stored at conditions stated on the label. The kit can be used up to 3 times within the expiry date stated on the label.
- The **ELISA WASHBUF** (wash buffer concentrate) should be diluted with aqua bidest. **1:20** before use (50 ml WASHBUF + 950 ml aqua bidest.), mix well. Crystals could occur due to high salt concentration in the stock solutions. The crystals must be redissolved at room temperature or at 37°C in a water bath before dilution of the buffer solutions. The **buffer concentrate** is stable at **2-8°C** until the expiry date stated on the label. **Diluted buffer solution** can be stored in a closed flask at **2-8°C for two weeks**.
- **Pre-heat RECSOL (reconstitution solution) for 5 min at 37°C in a water bath before use.**
- Reconstitute **RELREAG** (releasing reagent) in **16 ml RECSOL** (reconstitution solution), mix gently (do not vortex), and leave at room temperature for 5 min. After use, aliquot, freeze remaining releasing reagent and store at -20 °C. Bring thawed releasing reagent **to room temperature before use** (e.g. incubate for 5 minutes at 37 °C in a water bath).
- Bring **AB** (antibody) at room temperature at least one hour before use.
- All other test reagents are ready for use. The test reagents are stable up to the date of expiry (see label of test package) when stored at **2 –8 °C**. **Note:** Unused strips are covered and stored at 2-8° C. The covered strips should be used within 4 weeks.

6. SPECIMEN COLLECTION AND PREPARATION

1. Vitamin D is an inert substance. The samples can be stored at room temperature. However, sample storage at 2-8°C is recommended when the analysis is performed within 24 h after collection. Otherwise, the samples must be stored at -20°C until analyzed. Avoid repeated freeze-thaw cycles.
2. Prefer **serum** as sample matrix; whole blood is not suitable.
3. Indicated incubation times and temperatures must be strictly observed.
4. Mix samples well before use.

7. ASSAY PROCEDURE

Principle of the test

The assay utilizes of a competitive ELISA technique with a selected monoclonal antibody recognizing 25-hydroxy vitamin D. For a reliable determination of 25-OH-Vitamin D, it is necessary to release it from the 25-OH-Vitamin D-DBP-complex.

Standards, controls and patient samples which are assayed for 25-hydroxy vitamin D are incubated with the releasing reagent. The pre-incubated solutions are then transferred to the microplate coated with 25-hydroxy vitamin D, and an anti-25-hydroxy vitamin D antibody is added. During an over night incubation step, 25-hydroxy vitamin D in the sample and a fixed amount of 25-hydroxy vitamin D bound to the microtiter well compete for the binding of the antibody. Then a peroxidase-conjugated antibody is added into each microplate well. A complex of 25-hydroxy vitamin D - anti-25-hydroxy vitamin D antibody – peroxidase conjugate is formed. Tetramethylbenzidine (TMB) is used as a peroxidase substrate. Finally, an acidic stop solution is added to terminate the reaction, whereby the color changes from blue to yellow. The intensity of the yellow color is inversely proportional to the concentration of 25-hydroxy vitamin D. A dose response curve of the absorbance unit (optical density, OD at 450 nm) vs. concentration is generated using the values obtained from the standard. 25-hydroxy vitamin D in the samples is determined from this curve.

Test procedure

- | | |
|----|--|
| 1. | Prior to use in the assay allow all reagents and samples to come to room temperature (18 - 26 °C) and mix gentle, avoid foam formation |
| 2. | Mark the positions of STD (Standards)/SAMPLE/CTRL (Control) on a protocol sheet |
| 3. | Reconstitute RELREAG (Releasing reagent) (see P.5., Page 19) |
| 4. | Label V-tubes (e.g. 1.5 ml Eppendorf-tubes) |
| 5. | Pipette 30 µl of STD (Standard)/ SAMPLE/CTRL (Control) respectively, into the corresponding tube |
| 6. | Add 300 µl of RELREAG (releasing reagent) into each tube, vortex shortly |
| 7. | Incubate for 1 hour at 37 °C in a water bath or heating block (do not use an incubator) |
| 8. | Open tubes carefully and add 600 µL SAMDIL (sample dilution buffer). Close the tubes and vortex carefully. |

9.	Take microtiter strips out of the kit. Unused strips are covered and stored at 2-8° C. The covered strips should be used within 4 weeeeks
10.	Transfer 50 µl of STD (Standard)/ SAMPLE/CTRL (Control) from the V-tubes in duplicate to respective well
11.	Add 150 µl of AB (anti 25-OH-Vitamin D ₃ antibody) into each well
12.	Cover the plate tightly and incubate over night (min. 18 – max. 22 hours) at 8-10 °C in the dark
13.	Aspirate and wash the wells 5x with 250 µl of diluted wash buffer. The use of 8-channel pipette is recommended. Remove remaining wash buffer by hitting the plate against paper towel after the last wash For TECAN and Dynex instruments a programming protocol can be requested from Immundiagnostik AG
14.	Add 200 µl CONJ (Conjugate) into each well
15.	Cover the plate tightly and incubate for 1 hour at room temperature while shaking
16.	Aspirate and wash the wells 5x with 250 µl of diluted wash buffer. The use of 8-channel pipette is recommended. Remove remaining wash buffer by hitting the plate against paper towel after the last wash
17.	Add 200 µl of SUB (Substrate) into each well
18.	Incubate for 15 - 20 minutes at room temperature (18-26°C) in the dark
19.	Add 50 µl of STOP (Stop solution) into each well
20.	Determine the absorption with an ELISA reader at 450 nm against 620 nm (or 690 nm) as a reference. If no reference wave length is available, read only at 450 nm. If the extinction of the highest standard exceeds the range of the photometer, absorption must be measured immediately at 405 nm against 620 nm (690 nm) as a reference

8. RESULTS

The following algorithms can be used alternatively to calculate the results. We recommend the use of the „4-Parameter-algorithm“.

1. 4-parameter-algorithm

It is recommended a linear ordinate for optical density and a logarithmic abscissa for concentration. When using a logarithmic abscissa, the zero calibrator must be specified with a value less than 1 (e. g. 0.01).

2. Point-to-point-calculation

We recommend a linear ordinate for optical density and a linear abscissa for concentration.

The plausibility of the pairs of values should be examined before the automatic evaluation of the results. If this option is not available with the used program, a control of the paired values should be done manually.

9. LIMITATIONS

Samples with a 25-OH Vitamin D level greater than the highest calibrator should be diluted with RECSOL (reconstitution solution) and re-assayed.

Whole blood is not suitable.

10. QUALITY CONTROL

Control samples should be analyzed with each run. Results, generated from the analysis of control samples, should be evaluated for acceptability using appropriate statistical methods. The results for the patient samples may not be valid, if within the same assay one or more values of the quality control sample are outside the acceptable limits.

Expected values

Normal ranges for 25-OH-Vitamin D₃

1 ng/ml = 2.5 nmol/l

1 nmol/l = 0.4 ng/ml

Information from ASBMR 2006

Deficiency (seriously deficient)	< 12 ng/ml	resp. < 30 nmol/l
Insufficiency (deficient)	12 - 30 ng/ml	resp. 30 - 75 nmol/l
Sufficiency (adequately supplied)	> 30 ng/ml	resp. > 75 nmol/l

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Note

The vitamin D production in the skin is high variable and depends on the season- and daily time, degree of latitude, age, sun protection etc.

The normal ranges depend on the method used (e. g. vitamin-D-release from the vitamin D binding protein, DBP) and serve only as orientation.

Literature references

Visser M, Deeg DJ, Puts MT, Seidell JC, Lips P. (2006) Low serum concentrations of 25-hydroxyvitamin D in older persons and the risk of nursing home admission. *Am J Clin Nutr.* Sep;84(3):616-22; quiz 671-2

Grant WB, Holick MF. Benefits and requirements of vitamin D for optimal health: a review. (2005) *Altern Med Rev.* Jun;10(2):94-111. Review

Wicherts IS, van Schoor NM, Boeke AJ, Visser M, Deeg DJ, Smit J, Knol DL, Lips P. (2007) Vitamin D status predicts physical performance and its decline in older persons. *J Clin Endocrinol Metab.* Jun;92(6):2058-65

11. PERFORMANCE CHARACTERISTICS

Precision and reproducibility

Intra-Assay (n = 20)		
Sample	25-OH-Vitamin D [nmol/L]	CV [%]
1	41.4	10.0

Inter-Assay (n = 20)		
Sample	25-OH-Vitamin D [nmol/L]	CV [%]
1	39.0	8.0

Specificity

The specificity was tested by measuring the cross-reactivity against a range of compounds having structural similarity to 25-OH-Vitamin D. The specificity is calculated in per cent, based on the cross-reactivity with the anti-25-OH-Vitamin D antibody:

25-OH-Vitamin D ₃	100.0 %
25-OH-Vitamin D ₂	67.8 %
24, 25-OH-Vitamin D ₃	≥ 100.0 %
Vitamin D ₂ (Ergocalciferol)	0.3 %

Sensitivity

The sensitivity was set as B₀ - 2SD. The zero-standard was measured 20 times.

Sensitivity n=20

Sample	Vitamin D mean value [OD]	Standard variation	Detection limit [nmol/L]
1	1.31	0.05	2.0




Linearity


Two patient's serum samples were diluted with RECSOL (reconstitution solution) and analyzed. The results in nmol/l are shown below:

Sample	Dilution	Observed [nmol/L]	Expected [nmol/L]
A	undiluted	35.9	36.0
	1 : 2	18.4	18.0
	1 : 4	9.9	9.0
B	undiluted	118.0	120.0
	1 : 2	61.8	60.0
	1 : 4	27.8	30.0

Correlation data

Good correlation with HPLC / RIA / LC-MS / ELISA:

* mainly patients with Vitamin-D deficiencies were tested.

12. PRECAUTIONS

- For in vitro diagnostic use only.
- The quality control guidelines should be observed.
- Human material used in the kit components was tested and found to be negative for HIV, Hepatitis B and Hepatitis C. However, for safety reasons, all kit components should be treated as potentially infectious.
- Reagents of the kit package contain sodium azide and thimerosal as bactericides. Sodium azide and thimerosal are toxic. The substrates for the enzymatic color reactions are described to be also toxic and carcinogenic. Contact with skin or mucous membranes has to be avoided.
- Stop solution consists of sulfuric acid, a strong acid. Although diluted, it still must be handled with care. It can cause burns and should be handled with gloves, eye protection, and appropriate protective clothing. Any spill should be wiped out immediately with copious quantities of water.

13. TECHNICAL HINTS

- Do not interchange different lot numbers of any kit component within the same assay.
- Reagents should not be used beyond the expiration date shown on the kit label.
- Substrate solution should remain colourless until use.
- To ensure accurate results, proper adhesion of plate sealers during incubation steps is necessary.
- Avoid foaming when mixing reagents.
- The assay should always be performed according the enclosed manual.
- Incubation time, incubation temperature and pipetting volumes of the components are defined by the producer. Any variation of the test procedure, which is not coordinated with the producer, may influence the results of the test. Immundiagnostik AG can therefore not be held responsible for any damage resulting from wrong use.

14. GENERAL NOTES ON THE TEST AND TEST PROCEDURE

- This assay was produced and put on the market according to the IVD guidelines of 98/79/EC.
- All reagents in the kit package are for in vitro diagnostic use only.
- The guidelines for medical laboratories should be observed.
- Incubation time, incubation temperature and pipetting volumes of the components are defined by the producer. Any variation of the test procedure, which is not coordinated with the producer, may influence the results of the test. Immundiagnostik AG can therefore not be held responsible for any damage resulting from wrong use.
- Warranty claims and complaints in respect of deficiencies must be lodged within 14 days of receipt of the product. The product shall be send to Immundiagnostik together with the complaint in writing.

15. REFERENCES

1. Scharla SH et al. (1996) Exp Clin Endocrinol Diabetes 104:289-292
2. Scharla SH et al. (1998) Osteoporos Int 8 (Supplement 2):S7-S12
3. Chapuy MC et al. (1996) J Clin Endocrinol Metab 81:1129-33
4. Oster P et al. (1983) Akt Gerontol 13:221-2
5. Scharla S (1997) Schattauer Verlag, Stuttgart, Seiten 217-242
6. Offermann G (1978) Dtsch med Wschr 103:1387-1388

Used symbols:

Temperature limitation

Catalogue Number

In Vitro Diagnostic Medical Device

Contains sufficient for <n> tests

Manufacturer

Use by

Lot number

Immundiagnostik AG

Stubenwald-Allee 8a
D-64625 Bensheim

Tel.: +49(0) 62 51/70 19 00

Fax: +49(0) 62 51/84 94 30

info@immundiagnostik.com

www.immundiagnostik.com