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# Pankreatische Amylase ELISA Kit

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
# Pancreatic amylase ELISA Kit

*For the in vitro determination of human pancreatic amylase in  
stool*

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|   |            |
|---|------------|
| Table of content  | Seite/Page |
|   | 2          |
| <b>1. INTENDED USE</b>                                  | <b>13</b>  |
| <b>2. SUMMARY AND EXPLANATION OF THE TEST</b>           | <b>13</b>  |
| <b>3. PRINCIPLE OF THE TEST</b>                         | <b>13</b>  |
| <b>4. MATERIAL SUPPLIED</b>                             | <b>14</b>  |
| <b>5. MATERIAL REQUIRED BUT NOT SUPPLIED</b>            | <b>14</b>  |
| <b>6. PREPARATION AND STORAGE OF REAGENTS</b>           | <b>15</b>  |
| <b>7. PRECAUTIONS</b>                                   | <b>15</b>  |
| <b>8. SPECIMEN COLLECTION AND PREPARATION</b>           | <b>16</b>  |
| FAECES  | 16         |
| <b>9. ASSAY PROCEDURE</b>                               | <b>16</b>  |
| PROCEDURAL NOTES  | 16         |
| TEST PROCEDURE  | 16         |
| <b>10. RESULTS</b>                                      | <b>18</b>  |
| TYPICAL CALIBRATION CURVE                               | 18         |
| <b>11. LIMITATIONS</b>                                  | <b>19</b>  |
| <b>12. QUALITY CONTROL</b>                              | <b>19</b>  |
| EXPECTED VALUES   | 19         |
| <b>13. PERFORMANCE CHARACTERISTICS</b>                  | <b>20</b>  |
| SENSITIVITY   | 20         |
| <b>14. REFERENCES</b>                                   | <b>20</b>  |
| <b>15. GENERAL NOTES ON THE TEST AND TEST PROCEDURE</b> | <b>20</b>  |

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
# Pancreatic amylase ELISA Kit

*For the in vitro determination of human pancreatic amylase in stool*

Valid from 06.03.2008

## 1. INTENDED USE

The *Immundiagnostik* Assay is intended for the quantitative determination of **Pancreatic Amylase** in stool. For *in vitro* diagnostic use only.

## 2. SUMMARY AND EXPLANATION OF THE TEST

Like lipase and elastase the **Pancreatic Amylase** is synthesized in the pancreas. In case of a chronic pancreatitis, **Pancreatic Amylase** is decreased. Some forms of disease (alcoholism, accidental trauma) are associated with a pancreatic insufficiency. In the routine laboratory **Pancreatic Amylase** has proved itself to be a genuine alternative to Elastase-I in stool diagnosis.

### Indication

- Chronic pancreatitis
- Exocrine pancreas insufficiency

## 3. PRINCIPLE OF THE TEST

In a first incubation step, the Pancreatic Amylase in the samples is bound to monoclonal mouse antibodies (in excess), which are immobilized to the surface of the microtiter wells. To remove all unbound substances, a washing step is carried out. In a second incubation step, an anti-Pancreatic Amylase (POD-monoclonal antibody) antibody is added. After another washing step, to remove all unbound substances, the substrate, tetramethylbenzidine (TMB) is added. An acidic stop solution is then added to stop the reaction. The colour converts from blue to yellow. The intensity of the yellow color is directly proportional to the concentration of Pancreatic Amylase in the sample. A dose response curve of the absorbance unit (optical density, OD) vs. concentration is generated, using results obtained from the calibrators. Pancreatic Amylase, present in the patient samples, is determined directly from this curve.

The combination of two specific antibodies in the Pancreatic Amylase ELISA drastically reduces the possibility of wrong-negatives results and offers a secure diagnostic system to the user.

#### 4. MATERIAL SUPPLIED

| Cat. No   | Label   | Kit Components  | Quantity       |
|-----------|---------|---|----------------|
| K 6410MTP | PLATE   | One holder with precoated strips                              | 12 x 8         |
| K 6410WP  | WASHBUF | ELISA wash buffer concentrate 10x                             | 2 x 100 ml     |
| K 6410K   | CONJ    | Conjugate (mouse anti-pancreatic amylase, peroxidase-labeled) | 1 x 30 $\mu$ l |
| K 6410ST  | STD     | Calibrators, lyophilized<br>(0; 440; 1750; 7000; 28000 mU/l)  | 5 vials        |
| K 6410 KO | CTRL    | Control, lyophilized  | 1 vial         |
| K 6410TMB | SUB     | TMB substrate<br>(Tetramethylbenzidine)                       | 1 x 15 ml      |
| K 6410AC  | STOP    | ELISA stop solution, ready to use                             | 1 x 15 ml      |

On request we will send you 3 calibrator sets free of charge. Any further calibrator sets will be charged.

#### 5. MATERIAL REQUIRED BUT NOT SUPPLIED

- Bidistilled water (aqua bidest.)
- Laboratory balance
- Precision pipettors and disposable tips to deliver 10-1000  $\mu$ l
- Foil to cover the microtiter plate
- Horizontal microtiter plate shaker
- A multi-channel dispenser or repeating dispenser
- Centrifuge capable of 3000 x g
- Vortex-Mixer
- Standard laboratory glass or plastic vials, cups, etc.
- Microtiter plate reader at 450 nm  
(reference wave length 620 or 690 nm)

## 6. PREPARATION AND STORAGE OF REAGENTS

- To run the assay more than one time, make sure that the reagents are stored at the conditions stated on the label. **Prepare just the appropriate amount necessary for the assay.** The kit can be used up to 4 times within the expiry date stated on the label.
- The **ELISA WASHBUF** (wash buffer concentrate) should be diluted with aqua bidest. **1:10** before use (100 ml WASHBUF + 900 ml a. bidest.), mix well. Crystals could occur due to high salt concentration in the stock solutions. The crystals must be redissolved at 37°C in a water bath before dilution of the buffer solutions. The **WASHBUF** is stable at **2-8°C** until the expiry date stated on the label. Diluted **buffer solution** could be stored in a closed flask at **2-8°C for one month.**
- The **lyophilized STD** (standards) and **CTRL** (control) are stable at **2-8°C** until the expiry date stated on the label. The **STD** (standards) and **CTRL** (control) must be reconstituted with **250 µl** aqua bidest. Allow the vial content to dissolve for **10 minutes** and mix thoroughly by gentle inversion to insure complete reconstitution.
- The **CONJ** (POD antibody) must be diluted **1:1000** in wash buffer (10 µl CONJ + 10 ml wash buffer). The antibody is stable at 2 -4 °C until expiry date given on the label. **Diluted antibody solution is not stable and could not be stored.**
- All other test reagents are ready to use. The test reagents are stable at **2-8°C** until the expiry date stated on the label.

## 7. PRECAUTIONS

- For *in vitro* diagnostic use only.
- Human material used in the kit components was tested and found to be negative for HIV, Hepatitis B and Hepatitis C. However, for safety reasons, all kit components should be treated as potentially infectious.
- Reagents of the kit package contain sodium azide or thimerosal as bactericides. Sodium azide and thimerosal are toxic. The substrates for the enzymatic color reactions are toxic and carcinogenic. Avoid contact with skin or mucous membranes.

- Stop solution is composed of sulfuric acid, which is a strong acid. Even diluted, it still must be handled with care. It can cause acid burns and should be handled with gloves, eye protection, and appropriate protective clothing. Any spill should be wiped out immediately with copious quantities of water.
- Reagents should not be used beyond the expiration date shown on kit label.

## 8. SPECIMEN COLLECTION AND PREPARATION

### *Faeces*

Give about **100 mg** of the sample (please note the real weight for the calculation) to **5 ml** of the wash buffer and mix. Centrifuge the diluted sample for 10 min at 3000 rpm.

Afterwards the supernatant is diluted **1:40** in wash buffer (25 µl supernatant is added to 975 µl wash buffer). **100 µl** of this solution is used for the assay.

We recommend to weight out the stool samples for each run new. Supernatant is not stable and can't be stored.

Stool samples can be stored at  $-20^{\circ}\text{C}$  for 4 weeks. Avoid thawing and freezing cycles.

Immundiagnostik recommends the use of the sample tubes from Roche Diagnostics for sample preparation.

## 9. ASSAY PROCEDURE

### *Procedural notes*

- Do not interchange different lot numbers of any kit component within the same assay.
- The quality control guidelines should be observed.
- Incubation time, incubation temperature and pipetting volumes of the different components have been defined by the producer. Any variations of the test procedure, without consulting the producer with the producer, may influence the test results. Immundiagnostik AG can therefore not be held responsible for any damage.
- The assay should always be performed according the enclosed manual.

### *Test procedure*


Wash the PLATE (precoated microtiter plate) 5 x with 250 µl ELISA diluted wash buffer. After the final washing step remove residual buffer by tapping the plate on absorbent paper. Carry out the tests in duplicate.

1. Add **100 µl STD** (standards), **CTRL** (control) and **SAMPLE** (prediluted patient samples) in duplicate into the respective well.
2. Incubate for **1 hour** shaking on a horizontal mixer at room temperature.
3. Decant the content of the plate and wash the wells **5 x with 250µl** diluted ELISA wash buffer. After the final washing step remove residual buffer by tapping the plate on absorbent paper.
4. Add **100 µl of diluted CONJ** (peroxidase-labeled anti-pancreatic amylase antibody).
5. Incubate for **1 hour** shaking on a horizontal mixer at room temperature.
6. Decant the content of the plate and wash the wells **5 x with 250µl** diluted ELISA wash buffer. After the final washing step remove residual buffer by tapping the plate on absorbent paper.
7. Add **100 µl of SUB** (TMB substrate).
8. Incubate for **5-15 minutes** at room temperature.
9. Add **50 µl of STOP** (stop solution) and mix shortly.
10. Determine absorption immediately with an ELISA reader at **450 nm** against 620 nm (or 690 nm) as a reference. If no reference wavelength is available, read only at 450 nm. If the extinction of the highest standard exceeds the range of the photometer, absorption must be measured immediately at 405 nm against 620 nm as a reference.

## 10. RESULTS

We recommend the 4-parameter-algorithm for evaluation of the test results. Alternatively, a point-to-point-calculation or a spline-algorithm can be used. Before the automatic evaluation of the results, the plausibility of the pairs of values should be examined. If this option is not available with the used program, a control of the paired values should be done manually.

### *Typical calibration curve*

|                      |      |       |       |       |       |
|----------------------|------|-------|-------|-------|-------|
| Concentration [mU/l] | 0    | 440   | 1750  | 7000  | 28000 |
| OD Mean value        | 0.05 | 0.195 | 0.651 | 1.812 | 2.544 |

These data are for demonstration only and cannot be used instead of data obtained from the actual assay.

## Faeces

For the Pancreatic Amylase concentration of faeces samples, calculate as described in the following example:

weight: 80 mg (1ml Stool = 1g) = 0,08 ml

dilution step 1: 5ml / 0,08ml = 62,5

dilution step 2: 40

dilution factor: 62,5 x 40 = 2500

Multiply the result with **2500** to get the real concentration. The dilution factor depends on the weight of the faeces.

## 11. LIMITATIONS

Samples with Pancreatic Amylase levels greater than the highest calibrator should be diluted and re-assayed.

## 12. QUALITY CONTROL

Immundiagnostik recommends commercial control samples for internal quality control.

Control samples should be analyzed with each run. Results, generated from the analysis of control samples, should be evaluated for acceptability using appropriate statistical methods. The results for the patient samples may not be valid, if within the same assay one or more values of the quality control sample are outside the acceptable limits.

### *Expected values*

|                             |                 |
|-----------------------------|-----------------|
| Pancreatic Amylase (stool): | > 1000 U/l      |
| Grey area                   | 1000 – 1500 U/l |

### 13. PERFORMANCE CHARACTERISTICS

#### *Sensitivity*

The detection limit was defined as  $B_0 + 2SD$ .

n=20

| Sample | Pancreatic Amylase<br>Mean value<br>[OD] | Standard<br>variation | Detection limit<br>[mU/l] |
|--------|--|-----------------------|---------------------------|
| 1      | 0.065                                    | 0.030                 | 237                       |

### 14. REFERENCES

1. Bishop et al.: 1996; Pancreas, 13 (3), 226
2. Katschinsky et al.: 1997; Pancreas, 15 (2), 1991

### 15. GENERAL NOTES ON THE TEST AND TEST PROCEDURE

- This assay was produced and put on the market according to the IVD guidelines of 98/79/EC.
- Human materials used in kit components were tested and found to be negative for HIV, Hepatitis B and Hepatitis C. However, for safety reasons, all kit components should be treated as potentially infectious.
- Kit reagents contain sodium azide or thimerosal as bactericides. Sodium azide and thimerosal are toxic. Substrates for the enzymatic color reactions are toxic and carcinogenic. Avoid contact with skin or mucous membranes.
- All reagents in the kit package are for *in vitro* diagnostic use only.
- Guidelines for medical laboratories should be observed.
- Reagents should not be used beyond the expiration date shown on the kit label.
- Do not interchange different lot numbers of any kit component within the same assay.

- The assay should always be performed according the enclosed manual.
- Incubation time, incubation temperature and pipetting volumes of the components are defined by the producer. Any variation of the test procedure, which is not coordinated with the producer, may influence the results of the test. Immundiagnostik AG can therefore not be held responsible for any damage resulting from wrong use.
- Warranty claims and complaints in respect of deficiencies must be logged within 14 days after receipt of the product. The product shall be send to Immundiagnostik AG along with a written complaint.

3/6/2008 Amylase\_\_15112004.DOC

**Used symbols:**

Temperature limitation

Catalogue Number

In Vitro Diagnostic Medical Device

Contains sufficient for <n> tests

Manufacturer

Use by

Lot number