Chymotrypsin ELISA Kit

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Table of contents	Page
1. INTENDED USE	15
2. CLINICAL RELEVANCE	15
3. PRINCIPLE OF THE TEST	15
4. MATERIAL SUPPLIED	16
5. MATERIAL REQUIRED BUT NOT SUPPLIED	16
6. PREPARATION AND STORAGE OF REAGENTS	17
7. PRECAUTIONS	17
8. SPECIMEN COLLECTION AND PREPARATION	18
9. ASSAY PROCEDURE	18
PROCEDURAL NOTES	18
TEST PROCEDURE	19
10. RESULTS	20
11. LIMITATIONS	20
12. QUALITY CONTROL	21
EXPECTED VALUES	21
13. PERFORMANCE CHARACTERISTICS	21
PRECISION AND REPRODUCIBILITY	21
RECOVERY	22
SENSITIVITY	22
CROSS-REACTIVITY	22
14. GENERAL NOTES ON THE TEST AND TEST PROCEDURE	23

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Manual

Chymotrypsin ELISA Kit

For the in vitro determination of Chymotrypsin in stool


Valid from 22.10.2008

K 6910

96

Immundiagnostik AG, Stubenwald-Allee 8a, D 64625 Bensheim
Tel.: ++49 6251 70190-0
Fax: ++ 49 6251 849430
e.mail: Info@immundiagnostik.com
www.immundiagnostik.com

1. INTENDED USE

The *Immundiagnostik* Assay is a sandwich Enzyme Immuno Assay intended for the quantitative determination of **Chymotrypsin** in stool. For *in vitro* diagnostic use only.

2. CLINICAL RELEVANCE

Chymotrypsin is a serine protease secreted from the pancreas into the duodenum after food intake. **Chymotrypsin** cleaves food proteins preferentially next to aromatic residues. A small part of the active enzyme form is excreted within the stool. In the case of pancreas insufficiency secondary to a chronic pancreatitis, the secretion of the enzyme is reduced markedly, so that a reduced chymotrypsin activity can be measured in stool specimen.

The advantage of this assay is that the concentration of the enzyme is accurately determined and not the enzyme's activity.

Many studies have shown a longer stability of **chymotrypsin** (12 days at room temperature) compared to pancreatic elastase. Therefore, the **chymotrypsin ELISA** is a reliable and low-priced alternative to the pancreatic elastase assay.

Indications

- Chronic pancreatitis
- Exocrine pancreatic insufficiency

3. PRINCIPLE OF THE TEST

This Enzyme Immuno Assay is a sandwich assay for the determination of Chymotrypsin in stool samples. The wells of the microtitre plate are coated with polyclonal anti-Chymotrypsin antibodies. In a first incubation step, the Chymotrypsin in the samples is bound to the coated antibodies (in excess). To remove all unbound substances, a washing step is carried out.

In a second incubation step, a second anti-Chymotrypsin antibody is bound to the antigen. After another washing step, the solid phase is incubated with a peroxidase labelled antibody and then with the peroxidase substrate, Tetramethylbenzidine (TMB). The colour converts to yellow by the addition of an acidic stop solution. The intensity of the yellow color is directly proportional to the Chymotrypsin concentration of the sample. A dose response curve of the absorbance (at 450 nm) unit vs. concentration is generated. Chymotrypsin, present in the patient samples, is determined directly from this calibration curve.

The combination of two specific antibodies in the Chymotrypsin ELISA drastically reduces the possibility of wrong-negative results and offers a secure diagnostic system to the user.

4. MATERIAL SUPPLIED

Cat. No.	Label	Kit Components	Quantity
K 6910MTP	PLATE	one holder with precoated strips	96
K 6910WP	WASHBUF	ELISA wash buffer concentrate 10x	2 x 100 ml
K 6910A2	AB	Detection antibody, (Mouse anti-Chymotrypsin)	1 x 150 µl
K 6910K	CONJ	Conjugate (Streptavidin, peroxidase labeled)	1 x 150 µl
K 6910VP	ABBUF	Detection antibody dilution buffer	1 x 12 ml
K 6910ST	STD	Calibrators, ready-to-use (1000; 250; 125; 62.5; 15.6; 0 ng/ml)	6 x 250 µl
K 6910ko	CTRL	Control, ready-to-use	1 x 250 µl
K 6910TMB	SUB	TMB substrate (Tetramethylbenzidine), ready to use	1 x 15 ml
K 6910DI	DIL	Dilution buffer	1 x 11 ml
K 6910AC	STOP	ELISA stop solution, ready to use	1 x 15 ml

5. MATERIAL REQUIRED BUT NOT SUPPLIED

- Bidistilled water (aqua bidest.)
- Laboratory balance
- Precision pipettors and disposable tips to deliver 10-1000 µl
- Foil to cover the microtiter plate
- Horizontal microtiter plate shaker
- A multi-channel dispenser or repeating dispenser
- Centrifuge capable of 3000 x g
- Vortex-Mixer
- Standard laboratory glass or plastic vials, cups, etc.
- Microtiter plate reader at 450 or 405 nm (reference wave length 620 or 690 nm)

6. PREPARATION AND STORAGE OF REAGENTS

- To run assay more than once, ensure that reagents are stored at conditions stated on the label. **Prepare only the appropriate amount necessary for each assay.** The kit can be used up to 4 times within the expiry date stated on the label.
- Reagents with a volume less than **100 µl** should be centrifuged before use to avoid loss of volume.
- The **WASHBUF** (wash buffer concentrate) should be diluted with aqua bidest. **1:10** before use (100 ml WASHBUF + 900 ml aqua bidest.), mix well. Crystals could occur due to high salt concentration in the stock solutions. The crystals must be redissolved at 37°C in a water bath before dilution. The **WASHBUF** (wash buffer concentrate) is stable at **2-8°C** until the expiry date stated on the label. Diluted **buffer solution** can be stored in a closed flask at **2-8°C for one month**.
- The **AB** (detection antibody, biotinylated) must be diluted **1:100** in **ABBUF** (detection antibody dilution buffer; e.g. 100 µl AB + 10 ml ABBUF). The undiluted AB is stable at -20°C until the expiry date given on the label. **Diluted antibody solution is not stable and can not be stored.**
- The **CONJ** (conjugate, POD-antibody) must be diluted **1:100** in wash buffer (100 µl CONJ + 10 ml wash buffer). The undiluted CONJ is stable at **2-8 °C** until the expiry date stated on the label. **Diluted conjugate is not stable and can not be stored.**
- All other test reagents are ready to use. Test reagents are stable until the expiry date (see label of test package) when stored at **2-8°C**.

7. PRECAUTIONS

- For *in vitro* diagnostic use only.
- Human materials used in kit components were tested and found to be negative for HIV, Hepatitis B and Hepatitis C. However, for safety reasons, all kit components should be treated as potentially infectious.
- Kit reagents contain sodium azide or thimerosal as bactericides. Sodium azide and thimerosal are toxic. Substrates for the enzymatic color reactions are toxic and carcinogenic. Avoid contact with skin or mucous membranes.

- Stop solution is composed of sulfuric acid, which is a strong acid. Even diluted, it still must be handled with care. It can cause acid burns and should be handled with gloves, eye protection, and appropriate protective clothing. Any spills should be wiped out immediately with copious quantities of water.
- Reagents should not be used beyond the expiration date shown on the kit label.

8. SPECIMEN COLLECTION AND PREPARATION

Faeces

Weigh ca. **100 mg** of the sample, note the exact sample amount, add **5 ml** of the wash buffer and mix well.

Centrifuge the sample suspension for **10 min** at 3000 rpm. Transfer **1 ml of the supernatant** into an Eppendorf tube and centrifuge again at 13000 rpm for 5 min. The resulting supernatant can be stored at -20°C for about one month.

Centrifuge the supernatant at 13000 rpm for **2 min** before use. Use **20 µl** of the final centrifuged supernatant in the assay.

Immundiagnostik recommends the use of the sample tubes from Boehringer / Mannheim for sample preparation.

9. ASSAY PROCEDURE

Procedural notes

- Do not mix different lot numbers of any kit component.
- Quality control guidelines should be followed.
- Incubation time, incubation temperature and pipetting volumes of the components are defined by the producer. Any variation of the test procedure, which is not coordinated with the producer, may influence the results of the test. Immundiagnostik AG can therefore not be held responsible for any damage resulting from wrong use.
- The assay should always be performed according the enclosed manual.

Test procedure

The precoated microtiter strips have to be washed **5 x with 250 µl** wash buffer before use. Bring all reagents to room temperature. Mix reagents well before use. Avoid direct sun light during all incubation steps. Covering the microtitre plate during the different incubation steps is recommended. Carrying out the tests in duplicate is also recommended.

1. Add **100 µl** of **DIL** (dilution buffer) into each well.
2. Pipette **20 µl** of **STD** (standards), **CTRL** (control) or **samples** into each well.
3. Incubate for **1 hour** at room temperature shaking on a horizontal mixer.
4. Decant the contents of the plate and wash the wells **5 x with 250 µl** of wash buffer.
5. Add **100 µl** of diluted **AB** (detection antibody) solution.
6. Incubate for **1 hour** at room temperature shaking on a horizontal mixer.
7. Decant the contents of the plate and wash the wells **5 x with 250 µl** of wash buffer.
8. Add **100 µl** of diluted **CONJ** (conjugate) solution.
9. Incubate for **1 hour** at room temperature shaking on a horizontal mixer.
10. Decant the contents of the plate and wash the wells **5 x with 250 µl** of wash buffer.
11. Add **100 µl** of **SUB** (TMB-substrate) solution.
12. Incubate for **5 - 15 minutes** at room temperature in the dark until sufficient coloring is achieved.
13. Add **50 µl** of **STOP** (stop solution) and mix shortly.
14. Determine **absorption immediately** with an ELISA reader at **450 nm** against 620 nm (or 690 nm) as a reference. If no reference wavelength is available, read only at 450 nm. If the extinction of the highest standard exceeds the measurement range of the photometer, absorption must be measured immediately at 405 nm against 620 nm (or 690 nm) as reference.

10. RESULTS

The following algorithms can be used alternatively to calculate the results. We recommend to use the "4-Parameter-algorithm".

1. 4-parameter-algorithm

It is recommended to use a linear ordinate for optical density and a logarithmic abscissa for concentration. When using a logarithmic abscissa, the zero calibrator must be specified with a value less than 1 (e. g. 0.001).

2. Point-to-point-calculation

We recommend a linear ordinate for optical density and a linear abscissa for concentration.

3. Spline-algorithm

We recommend a linear ordinate for optical density and a logarithmic abscissa for concentration. When using a logarithmic abscissa, the zero calibrator must be specified with a value less than 1 (e. g. 0.001).

The plausibility of the pairs of values should be examined before the automatic evaluation of the results. If this option is not available with the used program, a control of the paired values should be done manually.

Faecal samples

The Chymotrypsin concentration of faecal specimen can be calculated as shown in the following example:

Weight: 80 mg (1ml stool = 1g) = 0.08 ml

Dilution step 1: 5 ml / 0.08 ml = 62.5

Dilution factor: 62.5

Multiply the results with the calculated dilution factor (in this case 62.5) to get the Chymotrypsin concentration of the stool samples.

Note: The dilution factor depends on the weight of the used faecal specimen.

11. LIMITATIONS

Samples with Chymotrypsin levels greater than the highest calibrator should be further diluted and re-assayed.

12. QUALITY CONTROL

Control samples or serum pools should be analyzed with each run. Results, generated from the analysis of the control samples, should be evaluated for acceptability using appropriate statistical methods. The results for the patient samples may not be valid, if within the same assay one or more values of the quality control sample are outside the acceptable limits.

Expected values

Reference range

Stool: > 6.4 µg/ml mean made by Immundiagnostik (n=80)

The reference value should be used as a guideline only. It is recommended that each laboratory establishes an own expected range for its patient population.

13. PERFORMANCE CHARACTERISTICS

Precision and reproducibility

Intra-Assay-Variation

The precision (intra-assay variation) of the Immundiagnostik Chymotrypsin ELISA test was calculated from 20 replicate determinations on each one of two samples.

Intra-Assay CV n= 20

Sample	Chymotrypsin, mean value [ng/ml]	Intra-Assay CV [%]
1	24,20	3,7
2	16,35	6,9

Inter-Assay-Variation

The total precision (inter-assay variation) of the Immundiagnostik Chymotrypsin ELISA test was calculated from data on one sample obtained in 31 different assays by three technicians on two different lots of reagents over a period of three months.

Inter-Assay CV n= 31

Sample	Chymotrypsin, mean value [ng/ml]	Inter-Assay CV [%]
1	101,8	5,3

Recovery

Two samples were spiked with different Chymotrypsin calibrator amounts and measured with the assay.

Recovery n=2

Sample [ng/ml]	Spike [ng/ml]	Chymotrypsin expected [ng/ml]	Chymotrypsin measured [ng/ml]
139,5	23,7	163,2	159,2
139,5	26,2	165,7	166,2
139,5	81,3	220,8	206,0
83,9	20,0	103,9	116,1
83,9	50,0	133,9	128,0
83,9	92,5	176,4	157,9

Sensitivity

The sensitivity limit was set as $B_0 + 2 \text{ SD}$. The Zero-standard was measured 28 times.

n=28

Sample	Chymotrypsin mean value [OD]	Standard variation	Detection limit [ng/ml]
1	0,006	0,006	14,2

Cross-reactivity

No cross-reactivity was observed to the following proteins:

Human pancreatic amylase

Human pancreatic elastase

Human lipase

Human PMN-elastase

14. GENERAL NOTES ON THE TEST AND TEST PROCEDURE

- This assay was produced and put on the market according to the IVD guidelines of 98/79/EC.
- Human materials used in kit components were tested and found to be negative for HIV, Hepatitis B and Hepatitis C. However, for safety reasons, all kit components should be treated as potentially infectious.
- Kit reagents contain sodium azide or thimerosal as bactericides. Sodium azide and thimerosal are toxic. Substrates for the enzymatic color reactions are toxic and carcinogenic. Avoid contact with skin or mucous membranes.
- All reagents in the test package are to be used for *in-vitro* diagnostics only.
- Reagents should not be used beyond the expiration date shown on the kit label.
- Do not mix different lot numbers of any kit component.
- Guidelines for medical laboratories should be observed.
- Incubation time, incubation temperature and pipetting volumes of the components are defined by the producer. Any variation of the test procedure, which is not coordinated with the producer, may influence the results of the test. Immundiagnostik AG can therefore not be held responsible for any damage resulting from wrong use.
- Warranty claims and complaints in respect of deficiencies must be logged within 14 days after receipt of the product. The product should be send to Immundiagnostik AG along with a written complaint.

Used symbols:

Temperature limitation

Catalogue Number

In Vitro Diagnostic Medical Device

Contains sufficient for <n> tests

Manufacturer

Use by

Lot number