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# GST- $\pi$

RIA

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*For the in vitro determination of Glutathion-S-Transferase  $\pi$  in Serum, Plasma and tumor tissue extract*

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---

Table of contents	Page
<b>1. INTENDED USE</b>	<b>12</b>
<b>2. SUMMARY AND EXPLANATION OF THE TEST</b>	<b>12</b>
<b>3. PRINCIPLE OF THE TEST</b>	<b>12</b>
<b>4. MATERIAL SUPPLIED</b>	<b>13</b>
<b>5. MATERIAL REQUIRED BUT NOT SUPPLIED</b>	<b>13</b>
<b>6. PREPARATION AND STORAGE OF REAGENTS</b>	<b>13</b>
<b>7. PRECAUTIONS</b>	<b>14</b>
<b>8. SPECIMEN COLLECTION AND PREPARATION</b>	<b>14</b>
PLASMA ODER SERUM	14
TISSUE MATERIAL :	14
<b>9. ASSAY PROCEDURE</b>	<b>15</b>
PROCEDURAL NOTES	15
TEST PROCEDURE	15
<b>10. RESULTS</b>	<b>16</b>
TYPICAL CALIBRATION CURVE	16
<b>11. LIMITATIONS</b>	<b>17</b>
<b>12. QUALITY CONTROL</b>	<b>17</b>
EXPECTED VALUES	17
<b>13. PERFORMANCE CHARACTERISTICS</b>	<b>17</b>
PRECISION AND REPRODUCIBILITY	17
RECOVERY	18
LINEARITY	18
SENSITIVITY	18
<b>14. REFERENCES</b>	<b>18</b>
<b>15. GENERAL NOTES ON THE TEST AND TEST PROCEDURE</b>	<b>19</b>

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# GST- $\pi$

RIA

*For the in vitro determination of Glutathion-S-Transferase  $\pi$  in  
Serum, Plasma and tumor tissue extract*

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CE

## 1. INTENDED USE

The *Immundiagnostik* Assay is intended for the quantitative determination of GST- $\pi$  in serum and tumor tissue extract. For in vitro diagnostic use only.

## 2. SUMMARY AND EXPLANATION OF THE TEST

Glutathion-S-transferases (GST) are a group of enzymes, which are involved in detoxification processes. The GSTs are divided in 5 subgroups: Alpha, Mu, Pi, Teta and microsomal. The classes differ in their physicochemical, immunological, enzymatic and structural properties. Their function in cellular metabolism is the coupling of electrophile compounds, for example most of the carcinogenes, to Glutathion (GSH). The GST Pi is the main group in erythrocytes and in some leucocytes, moreover they are found in most tissues except liver. The GST Pi are dimeric molecules with a molecular weight of 46000 Da and an isoelectric point of about 4,7. High GST Pi levels were found in all neoplastic tissue (fetal, maligne) and in serum of patients with post-operative recidives. Increased serum levels were also found in haematological diseases (haemolytic anaemia and sub-groups of leukemia). Moreover GST Pi concentrations give informations about the multi-drug-resistance-status of tumours.

## 3. PRINCIPLE OF THE TEST

This radioimmunoassay is capable for the quantification of GST- $\pi$ . GST- $\pi$  in the standards and the samples competes with <sup>125</sup>I-labelled GST- $\pi$  (tracer) for the binding sites of an anti-GST- $\pi$  antibody (1. Antibody). Bound antigen is separated from unbound by double antibody precipitation. The quantity of bound tracer is indirectly proportional to the GST- $\pi$  content.

#### 4. MATERIAL SUPPLIED

Catalogue No	Kit Components	Quantity
K 7961T	$^{125}\text{I}$ GST-pi < 110 kBq (red color)	1 x 10 ml
K 7961WP	RIA wash buffer, ready-to-use	1 x 100 ml
K 7961A1	1 <sup>st</sup> Antibody (rabbit anti h GST-pi)	1 x 10 ml
K 7961A2	2 <sup>nd</sup> Antibody (donkey anti rabbit IgG)	1 x 25 ml
K 8021ST	Calibrators, ready-to-use (NSB, 0, 0.74, 2.2, 6.6, 20, 60, 180, 540 ng/ml)	9 x 250 $\mu\text{l}$

#### 5. MATERIAL REQUIRED BUT NOT SUPPLIED

- Precision Pipettors to deliver 5 and 100  $\mu\text{l}$
- Repeating dispenser for 100 and 250  $\mu\text{l}$
- Vortex mixer
- Gamma Counter

#### 6. PREPARATION AND STORAGE OF REAGENTS

All test reagents are ready to use solutions and could be stored at  $-20^{\circ}\text{C}$  until the expiry date printed on the labels of the test components.

## 7. PRECAUTIONS

- The calibrators and controls could contain human source material which was tested and found to be non-reactive to HBsAg, anti-HIV-1/2, and anti-HCV. Since no method can offer complete assurance that hepatitis B virus, HIV-1/2, HVC or other infectious agents are absent, these reagents should be handled as if potentially infectious.
- The tracer contains  $^{125}\text{I}$  and is a radioactive solution. Although diluted, it still must be handled with care. It should be handled with gloves, eye protection, and appropriate protective clothing. Any spill should be wiped up immediately with copious quantities of water. Do not breath vapor and avoid inhalation.
- Reagents should not be used beyond the expiration date shown on kit label.

## 8. SPECIMEN COLLECTION AND PREPARATION

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Blood is collected in plasma- or serum tubes (don't use heparin-tubes and care of hemolysis !), shake carefully and after 15 - 30 minutes centrifuge at 1800 x g for 10 min. The collected serum or plasma should be stored rapidly at -20 C till use. **Do not use hemolysed blood.**

### *Tissue material :*

The tissue has to be homogenised (f.e. mechanically with a dismembrator) and then resuspended in phosphate buffer. Then follows an ultra centrifugation step (1h at 100 000 x g). The supernatant of these centrifugation step is the cytosolic fraction. Protein in the supernatant should be determined according to Lowry et al. or with BCA Protein Assay [Pierce]. Protein should be diluted to about 10 µg/ml and 2 x 100 µl / sample should be used for the RIA.

Alternatively the tissue material can be resuspended after the homogenisation in phosphatebuffer with 1% Triton X100. After mixing at least 2 hours this tissue lysate has to be centrifuged at 10 000 x g for 10 minutes. After determination of proteincontent the supernatant can be used for the RIA.

## 9. ASSAY PROCEDURE

### *Procedural notes*

- Do not interchange different lot numbers of any kit component within the same assay.
- The quality control guidelines should be observed.
- Incubation time, incubation temperature and pipetting volumes of the different components are defined by the producer. Any variations of the test procedure, that are not coordinated with the producer, may influence the test results. Immundiagnostik can therefore not be held reliable for any damage resulting from this.
- Carry out the assay with the actual manual delivered with the kit.

### *Test procedure*

The tests should be carried out in polypropylene RIA tubes in duplicates.


1. Pipette **50 µl** of NSB, standards and samples in RIA tubes.
2. Add **100 µl** tracer (125-I-GSTpi, red coloured)
3. Add **100 µl** 1. Antibody to each tube (not to NSB) and incubate **over night** (16 – 20 h) at room temperature
4. Add **250 µl** 2. antibody precipitation solution (mix vigorously prior to use)
5. Incubate for **60 min** at **room temperature**.
6. Centrifuge for **10 min.** with **1800 x g**, decant the supernatants.
7. Add **1 ml** washing solution to each sediment and mix vigorously.
8. Centrifugate for **10 min.** with **1800 x g**, decant the supernatants.
9. Count the sediment in the gamma-counter for one minute.

## 10. RESULTS

A calibration curve is constructed from the standards. Commercially available software can be used as well as graph paper. Results of the samples are read from this calibration curve.

THE CALIBRATION CURVE IS NOT LINEAR, therefore a spline-algorithm is recommended.

### *Typical calibration curve*

## 11. LIMITATIONS

Samples with GST-pi levels greater than the highest calibrator, should be diluted and re-assayed.

## 12. QUALITY CONTROL

Immundiagnostik recommends commercial control samples for internal quality control.

Control samples or serum pools should be analyzed with each run of calibrators and patient samples. Results, generated from the analysis of control samples, should be evaluated for acceptability using appropriate statistical methods. In assays in which one or more of the quality control sample values lie outside the acceptable limits, the results for the patient sample may not be valid.

### *Expected values*

Normal range (Plasma or Serum):

5-30 ng/ml.

## 13. PERFORMANCE CHARACTERISTICS

### *Precision and reproducibility*

The precision (intra-assay variation) of the Immundiagnostik GST-pi RIA test was calculated from 10 replicate determinations of the sample.

Intra-Assay CV n= 10

Sample	GST-pi Mean value [ng/ml]	Intra-Assay CV [%]
1	125	10

The total precision (inter-assay variation) of the Immundiagnostik GST-pi RIA test was calculated from data one sample obtained in 10 different assays by three technicians on two different lots of reagents over a period of 3 months.

Inter-Assay CV n= 10

Sample	GST-pi Mean value [ng/ml]	Inter-Assay CV [%]
1	125	12

### *Recovery*

One sample were spiked with different GST- $\pi$  concentrations and measured in the kit. The recovery was between 90-105 per cent.

### *Linearity*

The linearity of the RIA is 1 - 250 ng/ml.

### *Sensitivity*

The sensitivity or minimum detection limit of this assay is defined as the smallest single value that can be distinguished from zero at the 3 x standard variation as confidence limit ( $B_0 + 3 SD$ ) and was calculated with 1 ng/ml.

## 14. REFERENCES

1. Schäfer et al. A new Glutathione-S-Transferase assay for serological detection of gastrointestinal tumors. *Klin.Lab.* 11/1993

## 15. GENERAL NOTES ON THE TEST AND TEST PROCEDURE

- This assay was produced and put on the market according to the IVD guidelines of 98/79/EC.
- The test components which are made of human serum are tested for Australia antigen and HIV and found to be negative. However, since no test method can offer complete assurance that infectious agents are absent, these reagents should be handled as recommended for any potentially infectious human serum or blood specimen. The normal precautions for laboratory working should be observed.
- Reagents of the test package contain sodium azide as a bactericide. Contact with skin or mucous membranes has to be avoided.
- All reagents in the test package are to be used for in-vitro diagnostics only.
- The reagents should not be used after the date of expiry stated on the label.
- Single components with different lot numbers should not be mixed or exchanged.
- The guidelines for medical laboratories should be observed.
- Incubation time, incubation temperature and pipetting volumes of the different components are defined by the producer. Any variations of the test procedure, that are not coordinated with the producer, may influence the results of the test. Immundiagnostik can therefore not be held reliable for any damage resulting from this.