# Immunglobulin G ELISA Kit

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
# Immunoglobulin G ELISA Kit

*For the in vitro determination of Immunoglobulin G in serum,  
plasma and urine*

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Table of contents	Page
<b>1. INTENDED USE</b>	<b>11</b>
<b>3. MATERIAL SUPPLIED</b>	<b>11</b>
<b>4. MATERIAL REQUIRED BUT NOT SUPPLIED</b>	<b>12</b>
<b>5. PREPARATION AND STORAGE OF REAGENTS</b>	<b>12</b>
<b>6. PRECAUTIONS</b>	<b>13</b>
<b>7. SPECIMEN COLLECTION AND PREPARATION</b>	<b>13</b>
<b>8. ASSAY PROCEDURE</b>	<b>13</b>
PROCEDURAL NOTES	13
TEST PROCEDURE	14
<b>9. RESULTS</b>	<b>15</b>
<b>10. QUALITY CONTROL</b>	<b>16</b>
EXPECTED VALUES	16
<b>11. GENERAL NOTES ON THE TEST AND TEST PROCEDURE</b>	<b>16</b>

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
# Immunglobulin G ELISA Kit

*For the in vitro determination of Immunglobulin G in serum,  
plasma and urine*

Valid from 29.02.2008

K 6510

## 1. INTENDED USE

The *Immundiagnostik* Assay is intended for the quantitative determination of **Immunoglobulin G (IgG)** in plasma, serum and urine. For *in vitro* diagnostic use only.

## 2. PRINCIPLE OF THE TEST

In a first incubation step, the Immunoglobulin G in the samples is bound to polyclonal rabbit antibodies (in excess) immobilized to the surface of the microtitre wells. After removal of all unbound substances, a Peroxidase-labeled anti Immunoglobulin G antibody is added. The second washing step is followed by incubation with the substrate, tetramethylbenzidine (TMB). The reaction is terminated by an acidic stop solution converting the color from blue to yellow. The intensity of the yellow color is directly proportional to the concentration of Immunoglobulin G in the sample. A dose response curve of the absorbance unit (optical density, OD) vs. concentration is generated using the results obtained from the calibrators. Immunoglobulin G in the patient samples is determined directly from this curve.

## 3. MATERIAL SUPPLIED

Cat. No	Content	Kit Components	Quantity
K 6510MTP	PLATE	One holder with precoated strips	12 x 8
K 6510WB	WASHBUF	ELISA wash concentrate 10x	1 x 100 ml
K 6510K	CONJ	Conjugate, (rabbit-anti-IgG, Peroxidase-labeled)	1 x 50 µl
K 6510KV	CONJBUF	Conjugate dilution buffer, ready to use	22 ml
K 6510ST	STD	Calibrators, lyophilized (0; 0.32; 0.8; 2; 5 mg/l)	5 x
K 6510KO	CTRL	Control, lyophilized	1 x
K 6510NaCl	NACL	0.9 % NaCl solution	25 ml
K 6510TMB	SUB	TMB substrate (Tetramethylbenzidine)	2 x 15 ml
K 6510AC	STOP	ELISA stop solution, ready to use	1 x 15 ml

#### 4. MATERIAL REQUIRED BUT NOT SUPPLIED

- Bidistilled water (aqua bidest.)
- Laboratory balance
- Precision pipettors calibrated and tips to deliver 10-1000  $\mu$ l
- Covering foil for the microtiter plate
- Horizontal microtiter plate shaker with 37 °C incubator
- A multi-channel dispenser or repeating dispenser
- Centrifuge capable of 3000 x g
- Vortex-Mixer
- Standard laboratory glass or plastic vials, cups, etc.

Microtiter plate reader at 450 nm  
(reference wave length 620 or 690 nm)

#### 5. PREPARATION AND STORAGE OF REAGENTS

- To run assay more than once, ensure that reagents are stored at conditions stated on the label. **Prepare only the appropriate amount necessary for each assay.** The kit can be used up to 4 times within the expiry date stated on the label.
- Reagents with a volume less than **100  $\mu$ l** should be centrifuged before use to avoid loss of volume.
- The **WASHBUF** (wash buffer concentrate) should be diluted with aqua bidest. **1:10** before use (100 ml WASHBUF + 900 ml aqua bidest.), mix well. Crystals could occur due to high salt concentration in the stock solutions. The crystals must be redissolved at 37°C in a water bath before dilution of the buffer solutions. The **buffer concentrate** is stable at **2-8°C** until the expiry date stated on the label. **Diluted buffer solution** can be stored in a closed flask at **2-8°C for one month**.
- The **STD** (Calibrators) and the **CTRL** (control) must be reconstituted with **250  $\mu$ l** aqua dest. Allow the vial content to dissolve for 10 minutes and mix thoroughly by gentle inversion to insure complete reconstitution. Reconstituted calibrators and control can be stored at -20 °C until the expiry date given on the label. Repeated thawing and freezing should be avoided.
- The **CONJ** (conjugate) must be diluted **1: 1000** in **CONJBUF** (conjugate dilution buffer) (10  $\mu$ l CONJ + 10 ml CONJBUF). The undiluted **CONJ** (conjugate) is stable at **2-8 °C** until the expiry date stated on the label. **Diluted conjugate is not stable and can not be stored.**

- All other test reagents are ready for use. The test reagents are stable up to the date of expiry (see label of test package) when stored at 2–8 °C.

## **6. PRECAUTIONS**

- For *in vitro* diagnostic use only.
- Human materials used in kit components were tested and found to be negative for HIV, Hepatitis B and Hepatitis C. However, for safety reasons, all kit components should be treated as potentially infectious.
- Stop solution is composed of sulfuric acid, which is a strong acid. Even diluted, it still must be handled with care. It can cause acid burns and should be handled with gloves, eye protection, and appropriate protective clothing. Any spills should be wiped out immediately with copious quantities of water.
- Reagents should not be used beyond the expiration date shown on the kit label.

## **7. SPECIMEN COLLECTION AND PREPARATION**

### **Plasma or serum**

Samples can be stored for two weeks at 2-8°C. For longer storage, samples should be frozen at –20°C. Dilute all plasma and serum samples 1:10.000 with 0.9% NaCl.

### **Urine**

Adjust the urine to a pH of 6 to 8 with 1 NaOH and store samples at 2-8°C until testing. For longer storage, samples should be frozen at –20°C. Samples with IgG-concentration higher than 5 mg/l must be diluted 1:10 with 0.9% NaCl.

## **8. ASSAY PROCEDURE**

### *Procedural notes*

- Do not mix different lot numbers of any kit component.
- Quality control guidelines should be followed.
- Incubation time, incubation temperature and pipetting volumes of the components are defined by the producer. Any variation of the test procedure, which is not coordinated with the producer, may influence the results of the test. Immundiagnostik AG can therefore not be held responsible for any damage resulting from wrong use.
- The assay should always be performed according the enclosed manual.

### *Test procedure*

Wash the precoated microtiter plate 5 x with 250 µl ELISA wash buffer. Carry out the tests in duplicate.

1. Pipette **200 µl** of NACL (0.9% NaCl solution) into each well
2. Add **10 µl STD** (standard) **CTRL** (control) and **patient samples** (urine, plasma and serum diluted, see above).
3. Incubate for **1 hour** shaking on a horizontal mixer at room temperature.
4. Decant the content of the plate and wash the wells **5 x with 250 µl** ELISA wash buffer.
5. Add **200 µl** prediluted **CONJ** (conjugate).
6. Incubate for **1 hour** shaking on a horizontal mixer at room temperature.
7. Decant the content of the plate and wash the wells **5 x with 250 µl** ELISA wash buffer.
8. Add **200 µl** of **SUB** (TMB substrate solution).
9. Incubate for **10-15 minutes** at room temperature.
10. Add **50 µl STOP** (stop solution) and mix shortly.
11. Determine absorption with an ELISA reader at **450 nm** against 620 nm as reference. If no reference wavelength is available, read only at 450 nm. If the extinction of the highest standard exceeds the measurement range of the photometer, absorption must be measured immediately at 405 nm against 620 nm as reference.

## **9. RESULTS**

The following algorithms can be used alternatively to calculate the results. We recommend to use the "4-Parameter-algorithm".

### **1. 4-parameter-algorithm**

It is recommended to use a linear ordinate for optical density and a logarithmic abscissa for concentration. When using a logarithmic abscissa, the zero calibrator must be specified with a value less than 1 (e. g. 0.01).

### **2. Point-to-point-calculation**

We recommend a linear ordinate for optical density and a linear abscissa for concentration.

### **3. Spline-algorithm**

We recommend a linear ordinate for optical density and a logarithmic abscissa for concentration. When using a logarithmic abscissa, the zero calibrator must be specified with a value less than 1 (e. g. 0.01).

The plausibility of the pairs of values should be examined before the automatic evaluation of the results. If this option is not available with the used program, a control of the paired values should be done manually.

### **Serum/Plasma**

The result must be multiplied by **10.000** to calculate the serum value.

### **Urine (1:10 dilution)**

The result must be multiplied by **10** to obtain the IgG concentration in urine.

## 10. QUALITY CONTROL

Control samples should be analyzed with each run. Results, generated from the analysis of control samples, should be evaluated for acceptability using appropriate statistical methods. The results for the patient samples may not be valid, if within the same assay one or more values of the quality control sample are outside the acceptable limits.

### *Expected values*

Plasma or serum:	8-18 g/l
Urine:	0.5 – 3.2 mg/24 hours

Immundiagnostik recommends commercial control samples for internal quality control.

It is recommended that each laboratory should establish its own normal range. Above mentioned values are only for orientation and may vary from other published data.

## 11. GENERAL NOTES ON THE TEST AND TEST PROCEDURE

- This assay was produced and put on the market according to the IVD guidelines of 98/79/EC.
- Human materials used in kit components were tested and found to be negative for HIV, Hepatitis B and Hepatitis C and Australia antigen. However, for safety reasons, all kit components should be treated as potentially infectious.
- Reagents of the test package contain sodium azide as a bactericide. Contact with skin or mucous membranes must be avoided.
- All reagents in the test package are for *in-vitro* diagnostics only.
- Reagents should not be used beyond the expiration date shown on the kit label.
- The test components contain organic solvents. Contact with skin or mucous membranes must be avoided.
- Do not interchange different lot numbers of any kit component within the same assay.
- Quality control guidelines should be observed.

- Incubation time, incubation temperature and pipetting volumes of the components are defined by the producer. Any variation of the test procedure, which is not coordinated with the producer, may influence the results of the test. Immundiagnostik AG can therefore not be held responsible for any damage resulting from wrong use.

**Used symbols:**

Temperature limitation

Catalogue Number

In Vitro Diagnostic Medical Device

Contains sufficient for <n> tests

Manufacturer

Use by

Lot number