Lysozym ELISA Kit

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Hersteller

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Lysozyme ELISA Kit

For the in vitro determination of Lysozyme in serum and stool


Valid from 15.01.2009

K 6900

96

1. INTENDED USE

The *Immundiagnostik* Enzyme-Linked-Immuno-Sorbent-Assay (ELISA) is intended for the quantitative determination of **Lysozyme** in serum and stool. For *in vitro* diagnostic use only.

2. SUMMARY AND EXPLANATION OF THE TEST

Lysozyme (muramidase) is a protein with a molecular weight of approx. 15 kDa and belongs to the group of alkaline glycosidases. **Lysozyme** is produced by granulocytes, monocytes and macrophages. The main source for faecal **Lysozyme** are the intestinal granulocytes. **Lysozyme** can be detected in all cells of the inflammatory infiltrate during an acute attack of Crohn's disease. To some extent, **Lysozyme** is also secreted actively by mononuclear cells into the bowel lumen.

Indication

- Diagnosis and monitoring of Crohn's Disease
- Early diagnosis of rejection reactions in kidney transplantation cases
- Differential diagnosis and monitoring of leukosis
- Diagnosis and treatment monitoring of urinary tract infections in children
- Differential diagnosis between viral and bacterial meningitis in children
- Identification of sepsis in neonates

3. MATERIAL SUPPLIED

Catalogue No	Content	Kit Components	Quantity
K 6900MTP	PLATE	one holder with precoated strips	12 x 8 wells
K 6900WB	WASHBUF	ELISA wash concentrate 10x	2 x 100 ml
K 6900VP	VPBUF	Conjugate dilution buffer	1 x 15 ml
K 6900K	CONJ	Conjugate, (Rabbit-anti-Lysozyme, Peroxidase-labelled)	1 x 50 μ l
K 6900ST	STD	Calibrators, ready to use (0; 1.1; 3.3; 10; 30 ng/ml)	5 x 1 ml
K 6900KO	CTRL	Control	1 x 1 ml
K 6900TMB	SUB	TMB substrate (Tetramethylbenzidine)	1 x 15 ml
K 6900AC	STOP	ELISA stop solution, ready to use	1 x 15 ml

4 MATERIAL REQUIRED BUT NOT SUPPLIED

- Bidistilled water (aqua bidest.)
- Laboratory balance
- Precision pipettors calibrated to deliver 10-1000 μ l
- Covering foil for the microtiter plate
- Horizontal microtiter plate shaker
- Multi-channel dispenser or repeating dispenser
- Centrifuge capable of 3000 x g
- Vortex-Mixer
- Standard laboratory glass or plastic vials, cups, etc.
- Microtiter plate reader 450 nm

5. PREPARATION AND STORAGE OF REAGENTS

- To run the assay more than one time, make sure that the reagents are carefully stored at 2-8°C. Prepare just the appropriate amount necessary for the assay. The kit can be used up to 4 times within the expiry date stated on the label.
- Reagents with a volume less than 100 µl should be centrifuged before use to avoid losses.
- The **ELISA wash buffer concentrate** (WASHBUF) should be diluted with aqua bidest. **1:10** before use (add 900 ml of a. bidest. to 100 ml concentrate). Crystals could occur due to high salt concentration in the stock solutions. The crystals have to be redissolved at room temperature or at 37°C using a water bath **before dilution of the buffer solutions**. The buffer concentrates are stable at 2-8°C until the expiry date stated on the label. Diluted solutions could be stored at 2-8°C for 1 month.
- The **Conjugate** (CONJ, POD antibody) must be diluted **1:1000** in conjugate dilution buffer (10 µl CONJ and 10 ml PVBUF). The antibody is stable at 2-8 °C until expiry date given on the label. **Diluted antibody solution is not stable and can not be stored.**
- All other test reagents are ready for use. The test reagents are stable up to the date of expiry (see label of test package) when stored at 2-8°C.

6. SAMPLE PREPARATION

Extraction step of the stool sample

The **diluted wash buffer** is used as sample extraction buffer.

1. We recommend the use of a stool sample preparation kit for dosing of 100 mg of stool sample (e.g. Sample preparation kit from Roche Diagnostics, Mannheim, Germany; cat # 745804). Add 5 ml buffer to the stool sample in the preparation kit.

Constant buffer volume: 5 ml

Constant dilution factor: 1 : 50

2. Alternatively, stool samples can be manually weighted within the range of 80 – 120 mg. Please note the exact sample amount!
 - a. Add **5 ml** buffer to the stool sample independent of the sample amount.

Constant buffer volume: 5 ml

Variable dilution factor

The **dilution factor varies** depending on the sample amount which must be considered in the subsequent calculations. Please see the list of the correction factors below:

Weight [mg]	Dilution factor
80	62.5
82	60.9
84	59.5
86	58.1
88	56.8
90	55.6
92	54.3
94	53.2
96	52.1
98	51.0
100	50

Weight [mg]	Dilution factor
102	49.0
104	48.1
106	47.2
108	46.3
110	45.5
112	44.6
114	43.9
116	43.1
118	42.4
120	41.6

- b. **The buffer volume** for the individual samples **varies** depending on the sample amount (see table).

Variable buffer volume

Constant dilution factor: 1 : 50

Therefore, the same dilution factor can be used for all samples in the subsequent evaluation of the results .

Weight [mg]	Buffer Volume [ml]
80	4.0
82	4.1
84	4.2
86	4.3
88	4.4
90	4.5
92	4.6
94	4.7
96	4.8
98	4.9
100	5.0

Weight [mg]	Buffer Volume [ml]
102	5.1
104	5.2
106	5.3
108	5.4
110	5.5
112	5.6
114	5.7
116	5.8
118	5.9
120	6.0

Afterwards, mix the weighed stool sample with the buffer and vortex for 30 sec. Transfer ca. 1 ml stool suspension (dilution step I) in an Eppendorf-tube and centrifuge for 5 min at 13000 rpm.

Dilution of samples

The supernatant of the extraction (dilution step I) must be diluted 1:2 with wash buffer (e.g. 500 µl supernatant of the extraction + 500 µl wash buffer = dilution step II)

For analysis, pipet **100 µl** of the supernatant of dilution step II per well.

The supernatant is not stable and can not be stored. We recommend to weight fresh sample amount for a new assay, if the analysis should be repeated.

Serum

Serum should be centrifuged within one hour after collection. Store samples at -20 °C if not assayed on the same day. Lipemic or hemolytic samples may give false results. Samples should be mixed well before assaying. We recommend that STD, CTRL and SAMPLE are assayed in duplicate..

Samples are diluted between **1:500** and **1:1000** with wash buffer.

Use this dilution factor to calculate the Lysozyme concentration.

7. ASSAY PROCEDURE

Principle of the test

This Enzyme-Linked-Immuno-Sorbent-Assay (ELISA) is suitable for the quantitative determination of Lysozyme in human serum and stool. In a first incubation step, the Lysozyme in the patient samples is captured by the precoated antibody. The quantification of the bound Lysozyme is carried out by adding an rabbit anti Lysozym- Peroxidase-labelled conjugate. Finally the solid phase is incubated with the substrate, Tetramethylbenzidine (TMB). An acidic stop solution is used to stop the reaction. The color converts from blue to yellow. The intensity of the yellow color is directly proportional to the concentration of Lysozyme in the sample. A dose response curve of the absorbance unit (optical density, OD) vs. concentration is generated, using the results obtained from the calibrators. Lysozyme, present in the patient samples, is determined directly from this curve.

Test procedure

<p>Prior to use in the assay allow all reagents and samples to come to room temperature and mix by gentle swirling and inversion.</p>
<p>Mark the positions of STD (Standards)/SAMPLE/CTRL (Control) on a protocol sheet</p>
<p>Take microtiter strips out of the kit. Store unused strips covered at 2-8° C. Strips are stable until the expiry date stated on the label</p>
<p>Wash the wells 5x with 250 µl of diluted WASHBUF (Wash buffer), remove remaining WASHBUF by hitting the plate against paper towel after the last wash</p>
<p>Add 100 µl of STD/SAMPLE/CTRL (Standard/Sample/Control) in duplicate into respective well</p>
<p>Cover the plate tightly and incubate for 1 hour at room temperature (18-26°C) on a horizontal mixer</p>
<p>Aspirate and wash the wells 5x with 250 µl of diluted WASHBUF (Wash buffer), remove remaining WASHBUF by hitting the plate against paper towel after the last wash</p>
<p>Add 100 µl CONJ (Conjugate) into each well</p>
<p>Cover the plate tightly and incubate for 1 hour at room temperature (18-26°C) on a horizontal mixer</p>

Aspirate and wash the wells 5x with 250 µl of diluted WASHBUF (Wash buffer), remove remaining WASHBUF by hitting the plate against paper towel after the last wash

Add 100 µl of SUB (Substrate) into each well

Incubate for 10 - 20 minutes at room temperature (18-26°C) in the dark

Add 50 µl of STOP (Stop solution) into each well, shake well

Determine the absorption with an ELISA reader at 450 nm against 620 nm (or 690 nm) as a reference. If no reference wavelength is available, read only at 450 nm. If the extinction of the highest standard exceeds the range of the photometer, absorption must be measured immediately at 405 nm against 620 nm as a reference

8. RESULTS

The following algorithms can be used alternatively to calculate the results. However, we recommend to use the option 1: 4-parameter-algorithm.

1. 4-parameter-algorithm

It is recommended to use a linear ordinate for the optical density and a logarithmic abscissa for the concentration. When using a logarithmic abscissa, the zero calibrator has to be specified with a value smaller than 1 (e. g. 0.01).

2. Point-to-point-calculation

We recommend for the optical density a linear ordinate and for the concentration a linear abscissa.

3. Spline-algorithm

We recommend for the optical density a linear ordinate and for the concentration a logarithmic abscissa. When using a logarithmic abscissa, the zero calibrator has to be specified with a value smaller than 1 (e. g. 0.01).

The plausibility of the pairs of values should be examined before the automatic evaluation of the results. If this option is not available with the used program, a manual control of the paired values should be made.

Stool

To obtain the Lysozyme concentration in stool samples, multiply the estimated value by the dilution factor according to the sample preparation:

Sample preparation 1 and 2b: dilution factor constant : 1 : 50

Multiply the obtained result by **100** (dilution step I x dilution step II) to get the concentration.

Sample preparation 2a: The dilution factor is variable

Multiply the obtained result by the corresponding dilution factor taken from the table (dilution step I) and by 2 (dilution step II).

Serum

The result should be multiplied by **500 or 1000** to calculate the serum value.

9. LIMITATIONS

Samples with Lysozyme value greater than the highest calibrator should be diluted and re-assayed.

10. QUALITY CONTROL

Control samples should be analyzed with each run. Results, generated from the analysis of control samples, should be evaluated for acceptability using appropriate statistical methods. If one or more of the quality sample values are outside the acceptable limits, the results of the patient samples may be invalid.

Expected values

Lysozyme concentration (faeces): < 600 ng/ml
(obviously healthy persons)

It is recommended for each laboratory to establish its own normal range.

11. PERFORMANCE CHARACTERISTICS

Precision and reproducibility

Intra-Assay (n=20)		
Probe	Lysozyme [ng/ml]	Vk [%]
1	6.8	8
2	1.4	13

Inter-Assay (n=20)		
Probe	Lysozyme [ng/ml]	Vk [%]
1	5.6	17

Recovery

Two samples were spiked with Lysozyme calibrator and measured using this assay.

Recovery n=2

Sample [ng/ml]	Spike [ng/ml]	Lysozyme expected [ng/ml]	Lysozyme measured [ng/ml]
1.5	3	4.5	5.5
1.5	6	7.5	7.1
0.75	5	5.8	5.7
0.75	1.7	2.4	2.7
0.75	0.7	1.4	1.5

Sensitivity

n=20

Sample	Lysozyme mean value [OD]	Standard variation	Detection limit [ng/ml]
1	0.184	0.029	0.5

Cross reactivity

No cross reactivity to other proteins in serum and stool was observed.

Sample dilution

Linearity n= 2

Two patient serum samples were diluted with ELISA wash buffer. The results are shown below:

Sample	Dilution	Expected [ng/ml]	Measured [ng/ml]
A	undiluted	15.3	15.3
	1:2	7.7	7.5
	1:4	3.8	3.9
	1:8	1.9	<detection limit
B	undiluted	9.3	9.3
	1:2	4.65	6.0
	1:4	2.3	4.1
	1:8	1.16	1.6

12. PRECAUTIONS

- For *in vitro* diagnostic use only.
- The quality control guidelines should be observed.
- Human material used in the kit components was tested and found to be negative for HIV, Hepatitis B and Hepatitis C. However, for safety reasons, all kit components should be treated as potentially infectious.
- Reagents of the kit package contain sodium azide and thimerosal as bactericides. Sodium azide and thimerosal are toxic. The substrates for the enzymatic color reactions are described to be also toxic and carcinogenic. Contact with skin or mucous membranes has to be avoided.
- Stop solution consists of sulfuric acid, a strong acid. Although diluted, it still must be handled with care. It can cause burns and should be handled with gloves, eye protection, and appropriate protective clothing. Any spill should be wiped out immediately with copious quantities of water.

13. TECHNICAL HINTS

- Do not interchange different lot numbers of any kit component within the same assay.
- Reagents should not be used beyond the expiration date shown on the kit label.
- Substrate solution should remain colourless until added to the plate.
- To ensure accurate results, proper adhesion of plate sealers during incubation steps is necessary.
- Avoid foaming when mixing reagents.
- The assay should always be performed according the enclosed manual.
- Incubation time, incubation temperature and the volumes of the different components are defined by the producer. Any variation of the test procedure, which is not coordinated with the producer, may influence the test results. Immundiagnostik can therefore not be held responsible for any damage resulting from this.

14. GENERAL NOTES ON THE TEST AND TEST PROCEDURE

- This assay was produced and put on the market according to the IVD guidelines of 98/79/EC.
- The kit components made of human serum are tested for Hepatitis B, Hepatitis C and HIV and found to be negative. However, since no test method can offer complete assurance that infectious agents are absent, these reagents should be handled as recommended for any potentially infectious human serum or blood specimen. The normal precautions for laboratory work should be observed.
- All reagents in the kit package are for *in vitro* diagnostic use only.
- The guidelines for medical laboratories should be observed.
- Incubation time, incubation temperature and volumes of the different components are defined by the producer. Any variation of the test procedure, which is not coordinated with the producer, may influence the results of the test. Immundiagnostik can therefore not be held responsible for any damage resulting from this.
- Warranty claims and complaints in respect of deficiencies must be lodged within 14 days of receipt of the product. The product shall be send to Immundiagnostik together with the complaint in writing.

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Temperature limitation

Catalogue Number

In Vitro Diagnostic Medical Device

Contains sufficient for <n> tests

Manufacturer

Use by

Lot number