Myostatin EIA Kit

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Zachwieja JJ et al. (1999) J Gravit Physiol 6:11-15

Yarasheski KE (2002) J Nutr Health Aging 6: 343-348

Roth SM and Walsh S (2004) Curr Opin Clin Nutr Metab Care 7: 259-263

Shi X et al. (2007) PPAR Res 92501

Zhao B et al. (2005) Biochem Biophys Res Comm 337(1):248–255

Hittel et al. (2010) Medicine & Science in Sports & Exercise

DOI: 10.1249/MSS.0b013e3181e0b9a8

Verwendete Symbole:

Temperaturbegrenzung

Bestellnummer

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Hersteller

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
Myostatin EIA Kit

For the determination of Myostatin in human serum and plasma

Valid from 08.07.2010

K 1012

1. INTENDED USE

This enzyme immunoassay is intended for the quantitative determination of myostatin in human serum and plasma. It is for research use only.

2. SUMMARY

Myostatin belongs to the transforming growth differentiation factor- β (TGF- β) super family. The molecule is a negative regulator of muscle growth, but details about the actions of myostatin are uncertain (Roth and Walsh, 2004).

Myostatin was first identified in 1997 by McPherron et al. They found out that null-mutant knockout mice were significantly larger than wild-type animals and exhibited a large and widespread increase in skeletal muscle mass due to an increase of muscle fiber number (hyperplasia) and thickness (hypertrophy). Other groups identified mutations in the myostatin gene in naturally bred "double-muscles" cattle breeds.

Similar to the findings in animal models, increased myostatin immunoreactivity or expression has been observed in HIV-infected men with muscle wasting (Gonzales-Cadavid et al. 1998), after prolonged bed rest in young men (Zachwieja et al. 1999) and in older men and women with muscle wasting (Yarasheski KE et al. 2002).

Shi et al. (2007) and others have found that myostatin deficiency inhibits adipogenesis *in vivo*, even when mice are fed a high-fat diet. Transgenic overexpression of myostatin pro-peptide, which inhibits myostatin signaling, also inhibits body fat gain with a high-fat diet (Zhao et al. 2005). Similar alterations in myostatin signaling are associated with changes in body fat among humans.

Hittel et al. (2010) report that Myostatin-levels are regulated by aerobic exercise. Moreover, Myostatin is in the causal pathway of acquired insulin resistance with physical inactivity.

Indications

- Regulation of muscle growth
- Muscle atrophy
- Muscle wasting
- Acquired insulin resistance

4. MATERIAL SUPPLIED

Catalogue No	Content	Kit Components	Quantity
K 1012MTP	PLATE	One holder with precoated strips	12 x 8 wells
K 1012WP	WASHBUF	ELSA wash concentrate 10x	1 x 100 ml
K 1012VP	ABBUF	Antibody dilution buffer	1 x 22 ml
K 1012PV	DIL	Dilution buffer	1 x 22 ml
K 1012AK	AB	Anti-human Myostatin antibody (lyophilized)	2 x 10 ml
K 1012ST	STD	Standard concentrate (lyophilized)	2 x 380 µl (lyo.)
K 1012K	CONJ	Conjugate, peroxidase labeled, ready to use	1 x 15 ml
K 1012TMB	SUB	TMB substrate (Tetramethylbenzidine), ready to use	1 x 15 ml
K 1012AC	STOP	EIA stop solution, ready to use	1 x 15 ml

4. MATERIAL REQUIRED BUT NOT SUPPLIED

- Bidistilled water (aqua bidest.)
- Precision pipettors calibrated to deliver 10-1000 µl
- Horizontal microtiter plate shaker
- Multi-channel dispenser or repeating dispenser
- Vortex-Mixer
- Standard laboratory glass vials, cups, etc. (disposable products)
- Microtiter plate reader 450 nm
(reference wave length 620 or 690 nm)

5. PREPARATION AND STORAGE OF REAGENTS

- The **WASHBUF** (wash buffer concentrate) should be diluted with aqua bidest. **1:10** before use (100 ml WASHBUF + 900 ml aqua bidest.), mix well. Crystals could occur due to high salt concentration in the stock solutions. The crystals must be redissolved at room temperature or at 37°C in a water bath before dilution of the buffer solutions. The **buffer concentrate** is stable at **2-8°C** until the expiry date stated on the label. **Diluted buffer solution** can be stored in a closed flask at **2-8°C for two weeks**.
- Reconstitute **AB** (antibody) with 10 ml of **ABBUF** (antibody dilution buffer) and mix gently for 10 minutes at room temperature (do not vortex). After use, freeze the remaining reconstituted antibody and store at -20 °C up to maximal 2 weeks.
- Reconstitute **STDKONZ** (standard concentrate) with 380 µl of bidistilled water and mix gently for 10 minutes at room temperature (do not vortex); then use immediately. Reconstituted standard cannot be stored.
- Using the reconstituted standard prepare dilution series (1+2, respectively 1:3) with **DIL** (dilution buffer) in glass vials.

e. g.:

STD7: 150 µl STDKONZ results in	195 ng/ml
STD6: 75 µl STD 7 + 150 µl DIL results in	65 ng/ml
STD5: 75 µl STD 6 + 150 µl DIL results in	21,6 ng/ml
STD4: 75 µl STD 5 + 150 µl DIL results in	7,2 ng/ml
STD3: 75 µl STD 4 + 150 µl DIL results in	2,4 ng/ml
STD2: 75 µl STD 3 + 150 µl DIL results in	0,8 ng/ml
STD1: 150 µl DIL results in	0 ng/ml

(Please note: take out of STD2 a volume of 75 µl and discard it in order to have the same volume of 150 µl in all glass vials.)

- Unreconstituted **AB** (antibody) und unreconstituted **STDKONZ** (standard concentrate) must be stored at least at **-18°C** until use.
- All other test reagents are ready for use. The test reagents are stable up to the date of expiry (see label of test package) when stored at **2-8 °C**.

6. SAMPLE COLLECTION AND PREPARATION

Either EDTA plasma or serum can be used for analysis.

7. ASSAY PROCEDURE

Principle of the test

The test utilizes a competitive immunoassay technique. Standards with defined myostatin concentrations and patient samples which are assayed for myostatin are incubated with polyclonal antibodies against human recombinant myostatin in glass vials. (Do not use plastic vials, use glass vials only!). During the incubation, free myostatin can be bound by the myostatin antibodies. The pre-incubates are then transferred into the microtiter wells coated with recombinant myostatin.

The unbound antibodies in the pre-incubates can now bind to the immobilized antigen and are detected using a peroxidase-conjugated secondary antibody and TMB as a substrate.

Test procedure

1. Prior to use in the assay allow all reagents and samples to come to room temperature (18 - 26 °C) and mix gentle, avoid foam formation
2. Mark the positions of STD (Standards)/SAMPLE on a protocol sheet
3. Reconstitute AB (antibody) with 10 ml of ABBUF
4. Reconstitute STDKONZ (standard concentrate) with 380 µl bidistilled water
5. Label glass vials
6. Place glass vials containing the prepared standard series to the other glass vials
7. Pipette 30 µl of each sample in a new glass vial and dilute with 120 µl DIL (results in a dilution factor of 5)
8. Pipette 150 µl of diluted AB (inclusive standard series) into each glass vial, vortex well
9. Cover all glass vials and incubate for 2 hours at room temperature
10. Take as many microtiter strips as needed from the kit. Wash each well 2 times with 250 µl of diluted wash buffer. Use of an 8-channel pipette is recommended. After the final washing step, the inverted microtiter plate should be firmly tapped on absorbent paper to remove excess solution
11. Transfer 100 µl from each glass vial into a respective well in duplicate

12. Incubate microtiter plate while shaking (ca. 350 rpm) for 2 hours at room temperature (18-26°C)
13. Aspirate and wash the wells 5x with 250 µl of diluted wash buffer. The use of 8-channel pipette is recommended. Remove remaining wash buffer by hitting the plate against paper towel after the last wash
14. Add 100 µl of CONJ (conjugate) into each well
15. Cover strips tightly and incubate for 1 hour at room temperature (18-26°C) under shaking
16. Aspirate and wash the wells 5x with 250 µl of diluted wash buffer. The use of 8-channel pipette is recommended. Remove remaining wash buffer by hitting the plate against paper towel after the last wash
17. Add 100 µl of SUB (Substrate) into each well
18. Incubate for 10 - 20 minutes at room temperature (18-26°C) in the dark*
19. Add 50 µl of STOP (Stop solution) into each well
20. Determine the absorption with an EIA reader at 450 nm against 620 nm (or 690 nm) as a reference. If no reference wave length is available, read only at 450 nm. If the extinction of the highest standard exceeds the range of the photometer, absorption must be measured immediately at 405 nm against 620 nm (690 nm) as a reference

*The intensity of the color change is temperature sensitive. We recommend to observe the procedure of the color change and to stop the reaction upon good differentiation.

8. RESULTS

The following algorithms can be used alternatively to calculate the results. We recommend the use of the "4-Parameter-algorithm".

1. 4-parameter-algorithm

It is recommended a linear ordinate for optical density and a logarithmic abscissa for concentration. When using a logarithmic abscissa, the zero calibrator must be specified with a value less than 1 (e. g. 0.01).

2. Point-to-point-calculation

We recommend a linear ordinate for optical density and a linear abscissa for concentration.

The plausibility of the pairs of values should be examined before the automatic evaluation of the results. If this option is not available with the used program, a control of the paired values should be done manually.

9. LIMITATIONS

If the absorption of any sample is outside the measurement range, the sample should be diluted with DIL (dilution buffer) and re-assayed.

Whole blood is not suitable. Untreated lipemic samples may produce incorrect results.

10. QUALITY CONTROL


Immundiagnostik AG recommends the use of commercial control samples for internal quality control if available.

Control samples should be analyzed with each run. Results, generated from the analysis of control samples should be evaluated for acceptability using appropriate statistic methods. The results for the patient samples may not be valid, if within the same assay one or more values of the quality control sample are outside the acceptable limits.

11. PERFORMANCE CHARACTERISTICS

Linearity

One patient's sample with high myostatin concentration was diluted with DIL (dilution buffer) and analyzed. The obtained values are shown on the graph below.

linear signal response was obtained at dilutions of 1:256, 1:128, 1:64 and 1:32.

12. PRECAUTIONS

- All kit reagents are for research use only.
- The quality control guidelines should be observed.
- Human material used in the kit components was tested and found to be negative for HIV, Hepatitis B and Hepatitis C. However, for safety reasons, all kit components should be treated as potentially infectious.
- Reagents of the kit package contain sodium azide and thimerosal as bactericides. Sodium azide and thimerosal are toxic. The substrates for the enzymatic color reactions are described to be also toxic and carcinogenic. Contact with skin or mucous membranes has to be avoided.
- Stop solution consists of sulfuric acid, a strong acid. Although diluted, it still must be handled with care. It can cause burns and should be handled with gloves, eye protection, and appropriate protective clothing. Any spill should be wiped out immediately with copious quantities of water.

13. TECHNICAL HINTS

- Do not interchange different lot numbers of any kit component within the same assay.
- Reagents should not be used beyond the expiration date shown on the kit label.
- Substrate solution should remain colorless until use.
- To ensure accurate results, proper adhesion of plate sealers during incubation steps is necessary.
- Avoid foaming when mixing reagents.
- The assay should always be performed according the enclosed manual.
- Incubation time, incubation temperature and pipetting volumes of the components are defined by the producer. Any variation of the test procedure, which is not coordinated with the producer, may influence the results of the test. Immundiagnostik AG can therefore not be held responsible for any damage resulting from wrong use.

14. GENERAL NOTES ON THE TEST AND TEST PROCEDURE

- All reagents in the kit package are for research use only.
- The guidelines for medical laboratories should be observed.
- Incubation time, incubation temperature and pipetting volumes of the components are defined by the producer. Any variation of the test procedure, which is not coordinated with the producer, may influence the results of the test. Immundiagnostik AG can therefore not be held responsible for any damage resulting from wrong use.
- Warranty claims and complaints in respect of deficiencies must be lodged within 14 days of receipt of the product. The product shall be send to Immundiagnostik together with the complaint in writing.

15. REFERENCES

- Mc Pherron et al. (1997) Nature 387:83-90
Turner CH (2000) Bone 27 : 339-340
Yarasheski KE et al. (2002) J Nutrition, Health & Aging 6 : 343-348
Gonzales-Cadavid NF et al. (1998) Proc Natl Acad Sci USA 95 : 14938-14943
Hamrick MW et al. (2006) International J Obesity.30(5):868–870
Zachwieja JJ et al. (1999) J Gravit Physiol 6:11-15
Yarasheski KE (2002) J Nutr Health Aging 6: 343-348
Roth SM and Walsh S (2004) Curr Opin Clin Nutr Metab Care 7: 259-263
Shi X et al. (2007) PPAR Res 92501
Zhao B et al. (2005) Biochem Biophys Res Comm 337(1):248–255
Hittel et al. (2010) Medicine & Science in Sports & Exercise DOI:
10.1249/MSS.0b013e3181e0b9a8

08/07/2010 02082005_Myostatin.DOC

Used symbols:

Temperature limitation

Catalogue Number

For research use only

Contains sufficient for <n> tests

Manufacturer

Use by

Lot number