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Relaxin RIA Kit

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*For the in vitro determination of Relaxin in plasma, serum, urine,
cell culture supernatant, seminal plasma and tissue extract*

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Table of content	Page
1. INTENDED USE	12
2. SUMMARY AND EXPLANATION OF THE TEST	12
3. PRINCIPLE OF THE TEST	12
4. MATERIAL SUPPLIED	13
5. MATERIAL REQUIRED BUT NOT SUPPLIED	13
6. PREPARATION AND STORAGE OF REAGENTS	14
7. PRECAUTIONS	15
8. SPECIMEN COLLECTION AND PREPARATION	15
9. ASSAY PROCEDURE	16
PROCEDURAL NOTES	16
TEST PROCEDURE	17
10. RESULTS	18
TYPICAL CALIBRATION CURVE	18
11. LIMITATIONS	19
12. QUALITY CONTROL	19
CROSS REACTIVITY	19
13. PERFORMANCE CHARACTERISTICS	19
EXPECTED VALUES	19
14. REFERENCES	19
15. GENERAL NOTES ON THE TEST AND TEST PROCEDURE	20

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Manual

Relaxin RIA Kit

For the in vitro determination of Relaxin in plasma, serum, urine, cell culture supernatant, seminal plasma and tissue extract

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CE

1. INTENDED USE

The *Immundiagnostik* Relaxin Assay is intended for the quantitative determination of **Relaxin** in serum, plasma, seminal plasma, tissue extract and cell culture supernatant.

2. SUMMARY AND EXPLANATION OF THE TEST

Relaxin belongs to the insulin family and was originally described as a hormone produced in the ovaries, which was mainly detectable during pregnancy and birth. However, latest studies have shown, that synthesis of relaxin could also occur in the placenta, endometrium, mammary glands, or prostata. The binding of **relaxin** to specific receptors has been observed in the brain, uterus and the heart, which speaks in favour for the multi-functionality of this hormone. In the first animal models (rat) **relaxin** functions as a strong chronotropic and inotropic agent in normotensive as well as in spontane hypertensive rats. The mechanism, that leads to the increase in heart frequency and promote contractions of the heart muscles, is still unclear. The effects were however stronger in comparison to angiotensin II and it is assumed that this effect is mediated via specific receptors in the heart.

3. PRINCIPLE OF THE TEST

In this radioimmunoassay for the determination of relaxin unlabeled antigen in standards and samples compete with the ¹²⁵I-labeled relaxin (tracer) for the binding site of the rabbit antibody against human relaxin.

Bound antigen is separated from unbound antigen by means of a specific double antibody precipitation. After a washing standards and samples are counted for 1 minute in a γ -counter. The radioactivity of the bound fraction is indirectly proportional to the amount of relaxin.

4. MATERIAL SUPPLIED

Catalogue No	Kit Components	Quantity
K 9211AP	Assay buffer, ready-to-use	1 x 10 ml
K 9211T	¹²⁵ I labeled Relaxin (<110 kBq)	1 x 10 ml
K 9211A1	1 st Antibody (Rabbit anti hRelaxin)	1 x 10 ml
K 9211A2	2 nd Antibody (Donkey anti Rabbit IgG)	1 x 10 ml
K 9211VP	Dilution buffer, ready-to-use	1 x 50 ml
K 9211ST	Calibrator (Concentrate), lyophilisiert (4000 pg/ml)	1 x 1 vial

5. MATERIAL REQUIRED BUT NOT SUPPLIED

- Deionized water
- Precision pipettes calibrated to deliver 10-100 µl
- A multi-channel dispenser or repeating dispenser
- Centrifuge capable of 1800 x g
- Vortex-Mixer
- 0.9% NaCl Solution
- Gamma-Counter

6. PREPARATION AND STORAGE OF REAGENTS

- To run the assay more than once, make sure that the reagents are carefully stored as mentioned. Prepare just the appropriate amount necessary for the assay.
- The **calibrator concentrate** has to be resuspended in **1 ml aqua dest.**, resulting in a concentration of 4000 pg/ml. Diluted resuspended standard solution is stable at -20 °C, even if frozen and thawed several times. Further diluted standards should not be stored.

S1 = 4000 pg/ml

S2 = 250 µl **S1** + 250 µl dilution buffer = 2000 pg/ml

S3 = 250 µl **S2** + 250 µl dilution buffer = 500 pg/ml

S4 = 250 µl **S3** + 250 µl dilution buffer = 250 pg/ml

S5 = 250 µl **S4** + 250 µl dilution buffer = 125 pg/ml

S6 = 250 µl **S5** + 250 µl dilution buffer = 62 pg/ml

S7 = 250 µl **S6** + 250 µl dilution buffer = 31 pg/ml

S8 = 250 µl **S7** + 250 µl dilution buffer = 16 pg/ml

Dilution buffer is used as standard 0 pg/ml. For the NSB S8 is used.

- All other test reagents are ready for use. The test reagents are stable up to the date of expiry (see label of test package) when stored at 4°C.

7. PRECAUTIONS

- For in vitro diagnostic use only.
- The calibrators and controls could contain human source material which was tested and found to be non-reactive to HBsAg, anti-HIV-1/2, and anti-HCV. Since no method can offer complete assurance that hepatitis B virus, HIV-1/2, HVC or other infectious agents are absent, these reagents should be handled as if potentially infectious.
- The tracer contains ^{125}I and is a radioactive solution. Although diluted, it still must be handled with care. It should be handled with gloves, eye protection, and appropriate protective clothing. Any spill should be wiped up immediately with copious quantities of water. Do not breath vapor and avoid inhalation.
- Reagents should not be used beyond the expiration date shown on kit label.

8. SPECIMEN COLLECTION AND PREPARATION

Serum, plasma, urine or seminal plasma should be stored at -20°C until use.

Serum/Plasma

Serum and plasma samples have to be diluted at least 1:10 with sample dilution buffer if the samples have a concentration >4000 pg/ml.

Urine

Urine samples should be diluted at least 1:4 with sample dilution buffer

Seminal plasma

Seminal plasma has to be diluted at least 1:10 with sample dilution buffer

Tissue extract

- 200 mg or more deep frozen tissue samples, taken from liquid nitrogen were pulverized in the shaking holder of a micro-dismembrator (30 sec/1500 rpm).
- Homogenize the powder in 1 ml phosphate buffer((0,14 M NaCl; 2,6 mM KCl; 8 mM Na₂HPO₄; 1,4 mM KH₂PO₄), 1 % Triton-X 100; pH 7,4). After ultra-centrifugation (1 h/100.000 x g) the protein concentration was determined in the supernatant by commercially available Pierce-BCA or Peterson-Lowry Protein Assay.
- The supernatant is used in the assay.

9. ASSAY PROCEDURE*Procedural notes*

- Do not interchange different lot numbers of any kit component within the same assay.
- The quality control guidelines should be observed.
- Incubation time, incubation temperature and pipette volumes of the different components are defined by the producer. Any variations of the test procedure, that are not coordinated with the producer, may influence the test results. Immundiagnostik can therefore not be held reliable for any damage resulting from this.
- Carry out the assay with the actual manual delivered with the kit.

Test procedure

Carry out the assay in duplicate in polypropylen RIA tubes (12 x 75 mm)


1. Pipette **100 µl** of standard solutions, NSB, control solutions and patient samples.
2. Add **100 µl** 1st antibody to each tube (except NSB).
3. Incubate **20-24 h** at room temperature.
4. Add **100 µl** of tracer (¹²⁵I-labeled relaxin, red).
5. Incubate **3 h** at room temperature.
6. Add **100 µl** 2nd antibody to each tube (mix vigorously prior to use).
7. Add **100 µl** assay buffer to each tube.
8. Incubate **1 h** at room temperature.
9. Centrifugate for **10 min** at 1800 x g and decant the supernatant.
10. Add **1 ml** of 0,9 % NaCl to the precipitants.
11. Centrifugate for **10 min** at 1800 x g and decant the supernatant.
12. Count the precipitants in a gamma-counter for 1 min.

10. RESULTS

A calibration curve is constructed from the standards. Commercially available software can be used as well as graph paper. Results of the samples are read from this calibration curve.

THE CALIBRATION CURVE IS NOT LINEAR, therefore a spline-algorithm is recommended.

Typical calibration curve

Concentration [pg/ml]	2000	1000	500	250	125	62	31	16
CPM	1414	2699	4777	5671	7001	7447	8147	9022

The data is for demonstration only and cannot be used in place of data generations at the time of the assay.

For the actual concentration of the samples, the measured values have to be multiplied with the used dilution factor.

11. LIMITATIONS

Samples with Relaxin levels greater than the highest calibrator, should be further diluted and re-assayed.

12. QUALITY CONTROL

Cross reactivity

No cross reactivity with insulin

13. PERFORMANCE CHARACTERISTICS

Expected values

Plasma/Serum: < 120 pg/ml

14. REFERENCES

1. Dschietzig T, et al (2001) *FASEBJ* 15:2187-2195
2. Dschitzig et al (2000) 13. congress of the european society of cardiology Vo21, Netherlands
3. Echner et al (2000) 3rd international conference on Relaxin & Related peptides, Broome , Australien
4. Lippert T H et al.(1996) *Clin. Exp. Obst. Gyn.* 23:65-69
5. Seeger H et al (2000) *J Clin Pharm Therapeutics* 25:221-226
6. Ida Vogel Relaxin and Peterm Delivery (2000) PH.D. Thesis
7. Piedras-Renteria, E.S., et al. (1997) *Am J Physiol* 272:H1791-1799

15. GENERAL NOTES ON THE TEST AND TEST PROCEDURE

- The test components which are made of human serum are tested for Australia antigen and HIV and found to be negative. However, since no test method can offer complete assurance that infectious agents are absent, these reagents should be handled as recommended for any potentially infectious human serum or blood specimen. The normal precautions for laboratory working should be observed.
- Reagents of the test package contain sodium azide as a bactericide. Contact with skin or mucous membranes is to be avoided.
- All reagents in the test package are to be used only for in-vitro diagnostics.
- The reagents should not be used after the date of expiry (see label on the test package).
- Single components with different lot numbers should not be mixed or exchanged.
- For quality control, the guidelines for medical laboratories should be observed.
- The characteristic test data, such as incubation time, incubation temperature and pipetting volumina of the different components are defined by the producer. Any variations of the test procedure, that are not coordinated with the producer, may influence the results of the test. Our company can therefore not be held reliable for any damage resulting from this.