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# Resistin

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*For the determination of Resistin in serum and plasma*

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Table of content	Page
<b>1. INTENDED USE</b>	<b>14</b>
<b>2. CLINICAL IMPORTANCE OF THE TEST</b>	<b>14</b>
<b>3. PRINCIPLE OF THE TEST</b>	<b>14</b>
<b>4. MATERIAL SUPPLIED</b>	<b>15</b>
<b>5. MATERIAL REQUIRED BUT NOT SUPPLIED</b>	<b>15</b>
<b>6. PRECAUTIONS</b>	<b>16</b>
<b>7. PREPARATION AND STORAGE OF REAGENTS</b>	<b>16</b>
<b>8. SPECIMEN COLLECTION AND PREPARATION</b>	<b>16</b>
<b>9. ASSAY PROCEDURE</b>	<b>17</b>
PROCEDURAL NOTES	17
PREINCUBATION	17
PIPETTING SCHEME	18
<b>10. RESULTS</b>	<b>19</b>
TYPICAL CALIBRATION CURVE	19
<b>11. LIMITATIONS</b>	<b>20</b>
<b>12. QUALITY CONTROL</b>	<b>20</b>
<b>13. PERFORMANCE CHARACTERISTICS</b>	<b>20</b>
SENSITIVITY	20
<b>14. GENERAL NOTES ON THE TEST AND TEST PROCEDURE</b>	<b>21</b>

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Lagerung/Storage: 2-8 °C

## **1. INTENDED USE**

The Immundiagnostik assay is a sandwich ELISA intended for the quantitative determination of resistin in serum and plasma. For *in vitro* diagnostic use only.

## **2. CLINICAL IMPORTANCE OF THE TEST**

## **3. PRINCIPLE OF THE TEST**

The resistin test kit is a competitive enzyme immunoassay. The resistin from the sample competes with the tracer (immobilised resistin) for the binding site of the specific antibodies.

In the first step, the sample or the calibrator is preincubated with the antibodies in a reaction vessel. In the second incubation step, the antigen-antibody complex competes with the immobilized tracer for free binding sites of the anti-resistin antibody.

After a washing step, which removes all non-specific bound material, peroxidase-labeled conjugate is added to the wells. After removal of unbound conjugate by washing, tetramethylbenzidine (TMB) is added to the wells as substrate. Resistin is quantitated by an enzyme catalysed colour change detectable on a standard ELISA reader. The colour of the sample changes from blue to yellow after the addition of a stop reagent. The amount of colour developed is inversely proportional to the amount of resistin present in the sample. A calibration curve is drawn up in order to obtain the concentration.

#### 4. MATERIAL SUPPLIED

Catalogue no.	Kit Components	Quantity
K8029MTP	Pre-coated microtiter plate	12 x 8
K8029WP	ELISA washing buffer 10x	100 ml
K8029PV	Sample dilution buffer ready to use	1 x 15 ml
K8029AK	Antibody (goat anti-resistin) ready to use	1 x 10 ml
K8029K	Conjugate, (peroxidase labeled), ready- to-use	1x11 ml
K8029KO/1/2/3	Controls (3 ranges)	3 x 50 µl
K8029ST	Calibrators, lyophilized (NSB; 0.08, 0.4, 2, 10, 50 ng/ml)	6 vials
K8029TMB	TMB substrate (Tetramethylbenzidine), ready-to-use	1 x 15 ml
K8029AC	ELISA stop solution, ready-to-use	1 x 7 ml

#### 5. MATERIAL REQUIRED BUT NOT SUPPLIED

- Distilled or deionized water
- Precision pipettes with disposable tips calibrated to deliver 50 -1000 µl
- Multichannel or Multipipette
- Centrifuge 3000 x g
- Freezer (-20°C)
- Vortex mixer
- Conventional glass or plastic tubes (disposable articles)
- ELISA reader equipped with 450 nm filter

## 6. PRECAUTIONS

- Calibrator and controls are based on human serum. They have underwent testing for HIV and HBV and were found to be negative. All components from this test kit should, however, be treated as potentially infectious.
- Stop solution contains diluted sulphuric acid. Sulphuric acid is a strong acid and should be handled with care even when diluted. As sulphuric acid causes burns upon contact with skin, goggles and gloves should be worn. Flush immediately with ample amounts of water after contact.
- Do not use reagents beyond expiration date.

## 7. PREPARATION AND STORAGE OF REAGENTS

- **Washing buffer:** Dilute washing buffer concentrate to **1000 ml** (add 900 ml of distilled water to 100 ml of concentrate) and mix well. Crystals in the buffer concentrate will dissolve at room temperature or in a water bath (37°C). Buffer is stable for 1 month at 2-8 °C.
- **Calibrators:** Dissolve each vial in **500 µl** distilled water, mix well. The reconstituted calibrators are to be stirred for at least **10 min.** at room temperature to ensure a complete dissolution. The **sample dilution buffer** is used as the zero calibrator (**0 ng/ml**).
- **Controls** are dissolved in **500 µl** sample dilution buffer each and their test results multiplied by **10**.
- All the other reagents (antibody, conjugate, sample dilution buffer and the substrate) can be stored at 2-8°C until their respective expiry date.

## 8. SPECIMEN COLLECTION AND PREPARATION

Serum or plasma are stored at -20°C. Repeated freezing and thawing must be avoided.

## 9. ASSAY PROCEDURE

### *Procedural notes*

- Do not interchange different lot numbers of any kit component within the same assay.
- The quality control guidelines should be followed.
- Incubation time, incubation temperature and pipetting volumes of the different components are defined by the manufacturer. Any variations of the test procedure, that are not coordinated with the producer, may influence the test results. Immundiagnostik can therefore not be held responsible for any damage.
- Perform the assay with the **latest** manual delivered with the kit.

### *Preincubation*

1. Samples are diluted **1:10** with the sample dilution buffer.
2. Pipet **150 µl** calibrator, control or prediluted sample into a reaction vessel.
3. Add **150 µl** antibody. Mix well.
4. Incubate for **20-22 hours** at 2-8°C.

## *Pipetting scheme*

The precoated microtiterplate has to be washed 5 times with 350 µl washing buffer each. After the last washing step, dry by knocking the plates on absorbent tissue paper. All determinations are to be made in duplicate.


1. Pipet **100 µl** of the preincubate (sample+antibody or calibrator+antibody) into the respective microtiter wells.
2. Incubate and shake plate for **2 hours** at room temperature.
3. Discard contents of the wells, knock the plates on absorbent tissue paper and wash **5 times with 350 µl** washing buffer each. Knock again on absorbent tissue paper.
4. Add **100 µl** conjugate to the wells.
5. Incubate and shake plate for **1 hour** at room temperature.
6. Discard contents of the wells, knock the plates on absorbent tissue paper and wash **5 times with 350 µl** washing buffer each. Knock again on absorbent tissue paper.
7. Add **100 µl** of substrate solution (TMB).
8. Incubate for **15-20 min.** in the dark at room temperature (18-24°C).
9. Pipet **50 µl** stop solution into all wells.
10. Determine absorption **immediately** with an ELISA reader at **450 nm** against 690 or 620 nm as reference. Should the extinction of the zero calibrator exceed the detection range of the photometer, the measurement is to be repeated immediately at 405 nm (against 690 or 620 nm as reference).

## 10. RESULTS

A calibration curve is constructed from the calibrator results. Results of the samples are read from this calibration curve.

The calibration curve is not linear. We recommend a point-to-point or a 4-parameter algorithm for analysis. The extinction of the NSB can be deducted from the values of the measured extinctions.

### Typical calibration curve

concentration [pmol/l]	0	0.08	0.4	2.0	10.0	50.0
B/B <sub>0</sub> mean value	100	82	77	48	32	12

The data is supplied for demonstration purposes only and cannot be used in evaluating data generated from the assay.

## 11. LIMITATIONS

Samples with resistin levels greater than the highest calibrator should be diluted and re-assayed.

## 12. QUALITY CONTROL

Immundiagnostik recommends to use control samples from pooled serum samples as internal quality control.

Control samples or serum pools should be analyzed with each run of calibrators and patient samples. Results, generated from the analysis of control samples, should be evaluated for acceptability using appropriate statistical methods. In assays in which one or more of the quality control sample values are located outside the acceptable limits, the results for the patient sample may not be valid.

We recommend that every laboratory should establish its own reference range.

## 13. PERFORMANCE CHARACTERISTICS

### *Sensitivity*

The detection limit was determined to be 0.068 ng/ml ( $B_0 + 2 SD$ )

sample	resistin mean value [OD]	standard deviation	detection limit [ng/ml]
1	1.952	0.049	0.068

## **14. GENERAL NOTES ON THE TEST AND TEST PROCEDURE**

- This assay was produced and put on the market according to the IVD guidelines of 98/79/EC.
- The test components which are made of human serum are tested for HIV and HBV and found to be negative. However, since no test method can offer complete assurance that infectious agents are absent, these reagents should be handled as recommended for any potentially infectious human serum or blood specimen. The normal precautions for laboratory working should be observed.
- Reagents of the test package contain sodium azide as a bactericide. Contact with skin or mucous membranes has to be avoided.
- All reagents in the test package are to be used for *in-vitro* diagnostics only.
- The reagents should not be used after the date of expiry (see label on the test package).
- Single components with different lot numbers should not be mixed or exchanged.
- The guidelines for medical laboratories should be observed.
- Incubation time, incubation temperature and pipetting volumes of the different components are defined by the producer. Any variations of the test procedure, that are not coordinated with the producer, may influence the results of the test. Immundiagnostik can therefore not be held responsible for any damage resulting from this.