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# Zonulin ELISA Kit

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<b>1. INTENDED USE</b>	<b>14</b>
<b>2. SUMMARY AND EXPLANATION OF THE TEST</b>	<b>14</b>
<b>3. PRINCIPLE OF THE TEST</b>	<b>14</b>
<b>4. MATERIAL SUPPLIED</b>	<b>15</b>
<b>5. MATERIAL REQUIRED BUT NOT SUPPLIED</b>	<b>15</b>
<b>6. PREPARATION AND STORAGE OF REAGENTS</b>	<b>16</b>
<b>7. PRECAUTIONS</b>	<b>17</b>
<b>8. SPECIMEN COLLECTION AND PREPARATION</b>	<b>17</b>
<b>9. ASSAY PROCEDURE</b>	<b>17</b>
PROCEDURAL NOTES	17
TEST PROCEDURE	18
<b>10. RESULTS</b>	<b>19</b>
TYPICAL CALIBRATION CURVE	19
<b>11. LIMITATIONS</b>	<b>20</b>
<b>12. QUALITY CONTROL</b>	<b>20</b>
EXPECTED VALUES	20
<b>12. REFERENCES</b>	<b>20</b>
<b>14. GENERAL NOTES ON THE TEST AND TEST PROCEDURE</b>	<b>21</b>

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7. Watts T, Berti I, Sapone A et al. Role of the intestinal tight junction modulator zonulin in the pathogenesis of type 1 diabetes in BB diabetic-prone rats. *Proc Nat Acad Sci USA* 2005; 102(8): 2916-21.
8. Sapone A, de Magistris L, Pietzak M, Clemente MG, Tripathi A, Cucca F, Lampis R, Kryszak D, Carteni M, Generoso M, Iafusco D, Prisco F, Laghi F, Riegler G, Carratu R, Counts D, Fasano A. Zonulin upregulation is associated with increased gut permeability in subjects with type 1 diabetes and their relatives. *Diabetes* 2006; 55: 1443-9.
9. De Magistris MT. Zonula occludens toxin as a new promising adjuvant for mucosal vaccines. *Vaccine* 2006;24 Suppl 2: S2-60-1.

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Hersteller

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# Zonulin ELISA Kit


*For the in vitro determination of Zonulin in serum and EDTA-plasma*

*For research use only*

Valid from 06.11.2008

K 5600

## 1. INTENDED USE

The *Immundiagnostik* Assay is intended for the quantitative determination of zonulin in serum and EDTA-plasma. For research use only.

## 2. SUMMARY AND EXPLANATION OF THE TEST

Zonulin is a novel human protein analogue to the *Vibrio cholerae* derived Zonula occludens toxin, which participates in tight junctions between cells of the wall of the digestive tract. Zonulin binds to a specific receptor on the surface of intestinal epithelia and triggers a cascade of biochemical events which induces tight junction disassembly and a subsequent permeability increase of the intestinal epithelia, allowing some substances to pass through and activate immune reactions.

Dr. Fasano and his co-workers found out that the zonulin-zonulin-receptor-system is more activated in celiac disease and type 1 diabetes mellitus patients. Patients with active celiac disease showed higher levels of zonulin and anti-zonulin antibodies compared to non-celiac patients and patients in remission, who were eating a gluten-free diet.

Concerning the autoimmune type 1 diabetes, in experiments with rats could be demonstrated, that elevated zonulin levels as well as increased intestinal permeability precede a type 1 diabetes disease. Conversely, type 1 diabetes could be prevented by inhibition of zonulin in animal experiments.

In addition, it was reported that many people who suffer from celiac disease also suffer from other autoimmune disorders. It is suggested that increased levels of zonulin are a contributing factor to the development of celiac disease and other autoimmune disorders such as insulin dependent diabetes, multiple sclerosis, and rheumatoid arthritis

## 3. PRINCIPLE OF THE TEST

This enzyme immuno assay can be used for the quantitative determination of Zonulin in serum and EDTA-plasma.

The test principle is based on a competition between the free antigen in the samples or standards and the antigen coated on the wells of the microplate. Standards, samples and the primary anti-zonulin antibody are transferred directly into the precoated microplate wells. The antigen in the samples competes with the antigen immobilized on the wells of the microplate for the binding sites of the specific anti-zonulin antibody. A peroxidase-conjugated antibody is used for detection, and tetramethylbenzidine (TMB) as a peroxidase substrate. The enzymatic reaction is terminated by acidic stop solution. The quantification is based on the optical density at 450 nm.

#### 4. MATERIAL SUPPLIED

Cat. No	Label	Kit Components	Quantity
K 5600MTP	MTP	Microtiter plate, precoated	12 x 8 wells
K 5600WP	WASHBUF	BM-wash buffer concentrate 10 x	100 ml
K 5600EP	EXBUF	Extraction buffer	25 ml
K 5600KV	CONJBUF	Conjugate dilution buffer	15 ml
K 5600AK	AB	Zonulin-antibody	200 µl
K 5600VP	ABBUF	Zonulin-antibody dilution buffer	15 ml
K 5600K	2. AB	Conjugate	200 µl
K 5600ST	STD	Standards (lyophilized)	2 x 7 vials
K 5600KO	CTRL	Control (lyophilized)	2 x 1 vial
K 5600TMB	SUB	TMB substrate (Tetramethylbenzidine)	15 ml
K 5600AC	STOP	ELISA stop solution, ready to use	15 ml

#### 5. MATERIAL REQUIRED BUT NOT SUPPLIED

- Bidistilled water (aqua bidest.)
- Precision pipettors calibrated and tips to deliver 5-1000 µl
- Foil to cover the microtiter plate
- Horizontal microtiter plate shaker
- A multi-channel dispenser or repeating dispenser
- Centrifuge capable of 3000 x g
- Eppendorf-cups, 1,5 ml
- 15 ml tubes (e. g. Falcon)
- Vortex-Mixer
- Standard laboratory glass or plastic vials, cups, etc.
- Microtiter plate reader at 450 or 405 nm (reference wave length 620 or 690 nm)

## 6. PREPARATION AND STORAGE OF REAGENTS

- To run assay more than once ensure that reagents are stored at conditions stated on the label. **Prepare only the appropriate amount necessary for each assay.** The kit can be used up to 4 times within the expiry date stated on the label.
- Reagents with a volume less than **100 µl** should be centrifuged before use to avoid loss of volume.
- The **WASHBUF** (BM wash buffer concentrate) should be diluted with aqua bidest. **1:10** before use (100 ml WASHBUF + 900 ml aqua bidest.), mix well. Crystals could occur due to high salt concentration in the stock solutions. The crystals must be redissolved at 37°C in a water bath before dilution of the buffer solutions. The **WASHBUF** is stable at **2-8°C** until the expiry date stated on the label. **Diluted buffer solution** can be stored in a closed flask at **2-8°C for one month**.
- The lyophilized **STD** (standards) and **CTRL** (control) must be reconstituted with **500 µl** aqua bidest. Allow the vial content to dissolve for 10 minutes and mix thoroughly by gentle inversion to insure complete reconstitution. **Reconstituted standards and control can be stored at -20°C for one month. Avoid repeated freeze-thaw cycles.**
- The **AB** (Zonulin-antibody) must be diluted **1:100** in **ABBUF** (Zonulin-Zonulin-antibody dilution buffer) (50 µl AB + 4.95 ml ABBUF). The **AB** (Zonulin-antibody) is stable at **2 -8 °C** until expiry date given on the label. **Diluted AB is not stable and can not be stored.**
- The **2.AB** (Conjugate) must be diluted **1:100** in **CONJBUF** (Conjugate dilution buffer) (50 µl 2.AB + 4.95 ml CONJBUF). The **2.AB** (Conjugate) is stable at **2 -8 °C** until expiry date given on the label. **Diluted 2.AB is not stable and can not be stored.**
- All other test reagents are ready to use. Test reagents are stable until the expiry date (see label of test package) when stored at **2-8°C**.

## 7. PRECAUTIONS

- Human materials used in kit components were tested and found to be negative for HIV, Hepatitis B and Hepatitis C. However, for safety reasons, all kit components should be treated as potentially infectious.
- Kit reagents contain sodium azide or thimerosal as bactericides. Sodium azide and thimerosal are toxic. Substrates for the enzymatic color reactions are toxic and carcinogenic. Avoid contact with skin or mucous membranes.
- Stop solution is composed of sulfuric acid, which is a strong acid. Even diluted, it still must be handled with care. It can cause acid burns and should be handled with gloves, eye protection, and appropriate protective clothing. Any spills should be wiped out immediately with copious quantities of water.
- Reagents should not be used beyond the expiration date shown on the kit label.

## 8. SPECIMEN COLLECTION AND PREPARATION

### EDTA-Plasma and Serum

Store samples until use at -20 °C.

Pipette **250 µl** freshly collected EDTA-plasma or serum samples in Eppendorf-cups, add **250 µl EXBUF** (extraction buffer), vortex well and centrifuge for 20 min at 13000 rpm. Transfer the supernatant into an Eppendorf tube and centrifuge for 10 minutes at 13000 rpm. For analysis, pipet **50 µl** of the obtained supernatant per well.

**Dilution factor 1:2**

## 9. ASSAY PROCEDURE

### *Procedural notes*

- Do not interchange different lot numbers of any kit component within the same assay.
- Substrate solution should remain colorless until use.
- To ensure accurate results, proper adhesion of plate sealers during incubation steps is necessary.
- Avoid foaming when mixing reagents.
- The assay should always be performed according the enclosed manual.

### *Test procedure*

Wash the precoated microtiter plate **5 x with 250 µl** of diluted wash buffer. After the final washing step, the inverted microtiter plate should be firmly tapped on absorbent paper to remove excess solution. Run the test in duplicate.


1. Add **50 µl** of **STD** (standards), **CTRL** (control) or **samples** into each well in duplicate
2. Add **50 µl** of **AB** (anti-zonulin antibody) into each well
3. Cover the strips and incubate for **1,5 hours** shaking on a horizontal mixer in the dark at room temperature (18-26°C)
4. Aspirate the contents of each well. Wash the microtiter plate **5 x with 250 µl** of diluted wash buffer. After the final washing step, the inverted microtiter plate should be firmly tapped on absorbent paper to remove excess solution
5. Add **100 µl** of **2.AB** (conjugate) into each well
6. Cover the strips and incubate for **1,5 hours** shaking on a horizontal mixer in the dark at room temperature (18-26°C)
7. Aspirate the contents of each well. Wash the microtiter plate **5 x with 250 µl** of diluted wash buffer. After the final washing step, the inverted microtiter plate should be firmly tapped on absorbent paper to remove excess solution
8. Add **100 µl** of **SUB** (substrate solution) into each well
9. Incubate for **5-10 minutes** at room temperature (18-26°C)
10. Add **50 µl** of **STOP** (stop solution) into each well and mix shortly in the ELISA reader using the shake option
11. Determine absorption **immediately** with an ELISA reader at **450 nm** against 620 nm (alternatively at 690 nm) as a reference. If no reference wavelength is available, read only at 450 nm. If the extinction of the highest standard exceeds the range of the photometer, absorption must be measured immediately at 405 nm against 620 nm as a reference

## 10. RESULTS

A calibration curve is constructed from the calibrators. Commercially available software can be used as well as graph paper. Results of the samples are read from this calibration curve.

THE CALIBRATION CURVE IS NOT LINEAR, therefore a spline- or 4PL algorithm is recommended.

### *Typical calibration curve*

The calibration curve is for demonstration only. The data cannot be used for the evaluation of test results.

## 11. LIMITATIONS

If the absorption of any sample is outside the measurement range, the sample should be diluted with EXBUF (extraction buffer) and re-assayed.

## 12. QUALITY CONTROL

**Control samples** should be analyzed with each run. Results, generated from the analysis of control samples, should be evaluated for acceptability using appropriate statistical methods. The results for the patient samples may not be valid, if within the same assay one or more values of the quality control sample are outside the acceptable limits.

### *Expected values*

A **dilution factor of 2** is used to obtain the Zonulin concentration of the samples.

### **Normal range**

We recommend each laboratory to establish its own baseline values.

## 12. REFERENCES

1. Wang W, Uzzau S, Goldblum SE, Fasano A. Human zonulin, a potential modulator of intestinal tight junctions. *J Cell Sci* 2000;113 Pt 24: 4435-40.
2. Fasano A. Intestinal zonulin: open sesame! *Gut* 2001; 49: 159-62.
3. Fasano A, Not T, Wang W et al. Zonulin, a newly discovered modulator of intestinal permeability, and its expression in celiac disease. *Lancet* 2000; 355(9214): 1518-19.
4. Thomas KE, Sapone A, Fasano A, Vogel SN. Gliadin stimulation of murine macrophage inflammatory gene expression and intestinal permeability are MyD88-dependent: role of the innate immune response in Celiac disease. *J Immunol* 2006; 176: 2512-21.
5. Freemark M, Lynne LL. Screening for Celiac Disease in Children with Type 1 Diabetes. *Diabetes Care* 2003; 26: 1932-39.

6. Lazzarotto F, Basso D, Plebani M et al. Celiac Disease and Type 1 Diabetes. *Diabetes Care* 2003; 26: 248-49.
7. Watts T, Berti I, Sapone A et al. Role of the intestinal tight junction modulator zonulin in the pathogenesis of type 1 diabetes in BB diabeticprone rats. *Proc Nat Acad Sci USA* 2005; 102(8): 2916-21.
8. Sapone A, de Magistris L, Pietzak M, Clemente MG, Tripathi A, Cucca F, Lampis R, Kryszak D, Carteni M, Generoso M, Iafusco D, Prisco F, Laghi F, Riegler G, Carratu R, Counts D, Fasano A. Zonulin upregulation is associated with increased gut permeability in subjects with type 1 diabetes and their relatives. *Diabetes* 2006; 55: 1443-9.
9. De Magistris MT. Zonula occludens toxin as a new promising adjuvant for mucosal vaccines. *Vaccine* 2006;24 Suppl 2: S2-60-1.

## 14. GENERAL NOTES ON THE TEST AND TEST PROCEDURE

- All reagents in the kit package are for research use only.
- Guidelines for medical laboratories should be observed.
- Incubation time, incubation temperature and pipetting volumes of the components are defined by the producer. Any variation of the test procedure, which is not coordinated with the producer, may influence the results of the test. Immundiagnostik AG can therefore not be held responsible for any damage resulting from wrong use.
- Warranty claims and complaints in respect of deficiencies must be logged within 14 days after receipt of the product. The product shall be send to Immundiagnostik AG along with a written complaint.

**Used symbols:**

Temperature limitation

Catalogue Number

In Vitro Diagnostic Medical Device

Contains sufficient for &lt;n&gt; tests

Manufacturer

Use by

Lot number

For research use only