
# $\alpha_2$ -Makroglobulin ELISA Kit

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# $\alpha_2$ -Macroglobulin ELISA Kit

*For the in vitro determination of  $\alpha_2$ -Macroglobulin in urine, serum and plasma*

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---

Table of contents	Page
<b>1. INTENDED USE</b>	<b>16</b>
<b>2. SUMMARY AND EXPLANATION OF THE TEST</b>	<b>16</b>
<b>3. PRINCIPLE OF THE TEST</b>	<b>17</b>
<b>4. MATERIAL SUPPLIED</b>	<b>17</b>
<b>5. MATERIAL REQUIRED BUT NOT SUPPLIED</b>	<b>18</b>
<b>6. PREPARATION AND STORAGE OF REAGENTS</b>	<b>18</b>
<b>7. PRECAUTIONS</b>	<b>19</b>
<b>8. SPECIMEN COLLECTION AND PREPARATION</b>	<b>19</b>
<b>9. ASSAY PROCEDURE</b>	<b>19</b>
PROCEDURAL NOTES	19
TEST PROCEDURE	20
<b>10. RESULTS</b>	<b>21</b>
TYPICAL CALIBRATION CURVE	21
<b>11. LIMITATIONS</b>	<b>22</b>
<b>12. QUALITY CONTROL</b>	<b>22</b>
EXPECTED VALUES	22
<b>13. PERFORMANCE CHARACTERISTICS</b>	<b>23</b>
PRECISION AND REPRODUCIBILITY	23
RECOVERY	24
SENSITIVITY	24
SAMPLE DILUTION	25
<b>14. REFERENCES</b>	<b>25</b>
<b>15. GENERAL NOTES ON THE TEST AND TEST PROCEDURE</b>	<b>26</b>

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
Manual

# $\alpha_2$ -Macroglobulin ELISA Kit

*For the in vitro determination of  $\alpha_2$ -Macroglobulin in urine,  
serum and plasma*

Valid from 12.12.2007

## 1. INTENDED USE

The *Immundiagnostik* Assay is intended for the quantitative determination of  $\alpha$  2-Macroglobulin in urine, serum and plasma. For *in vitro* diagnostic use only.

## 2. SUMMARY AND EXPLANATION OF THE TEST

**Alpha 2-Macroglobulin** is one of the biggest plasma proteins with a molecular weight of 650-900kDa, dependent on the degree of glycosylation. It consists of 4 subunits. **Alpha 2-Macroglobulin** functions as protease inhibitor of all known classes of endopeptidase (serine-, cysteine-, aspartat- and metalloproteases).

The measurement of urinary proteins in the urine allows the diagnosis of proteinuria. In case of a microalbuminuria (selective proteinuria) increased levels (> 150 mg protein/24h) of albumin and transferrin (marker proteins of the tubulus) could be measured. In case of severe diseases of the kidney with pathological changes of the basal membran in the glomerulus (Chronic Glomerulonephritis, Tubulopathy, Nephrotic Syndrom) the protein level will rise up to > 3 g/24h and also the proteintyp will change to high molecular proteins like  $\beta$ -Lipoproteins and  $\alpha$  2-Macroglobulin (non-selective proteinuria).

Additionally previous investigations have shown that the determination of the two acute-phase proteins in the urins, C-reactive protein (CRP) and  $\alpha$  2-Macroglobulin, allows a noninvasive diagnosis of acute renal transplant dysfunction and rejection.

### Indication

- Marker for glomerulonephropathy
- Monitoring of renal transplant dysfunction and rejection (together with CRP)

### 3. PRINCIPLE OF THE TEST

In a first incubation step, the  $\alpha$  2-Macroglobulin in the samples is bound to polyclonal rabbit antibodies (in excess), which are immobilized to the surface of the microtitre wells. To remove all unbound substances, a washing step is carried out. In a second incubation step, a Peroxidase-labeled anti  $\alpha$  2-Macroglobulin (POD-Antibody) antibody is added. After another washing step, to remove all unbound substances, the solid phase is incubated with the substrate, tetramethylbenzidine (TMB). An acidic stop solution is then added to stop the reaction. The color converts from blue to yellow. The intensity of the yellow color is directly proportional to the concentration of  $\alpha$  2-Macroglobulin in the sample. A dose response curve of the absorbance unit (optical density, OD) vs. Concentration is generated, using results obtained from the calibrators.  $\alpha$  2-Macroglobulin, **present in the patient samples, is determined directly from this curve.**

### 4. MATERIAL SUPPLIED

Cat. No	Content	Kit Components	Quantity
K 6610MTP	PLATE	One holder with pre-coated strips	12 x 8 wells
K 6610WB	WASHBUF	ELISA wash buffer concentrate (10x)	1 x 100 ml
K 6610K	CONJ	Conjugate, (rabbit anti $\alpha$ -2-Macroglobulin peroxidase-labeled)	1 x 20 $\mu$ l
K 6610ST	STD	Calibrators, lyoph.	6 vials
K 6610KO	CTRL	Control, lyoph.	1 vial
K6610VP	SAMPLEBUF	Sample dilution buffer, ready-to-use	1 x 100 ml
K 6610NaCl	NACL	0.9 % NaCl-solution, ready-to-use	25 ml
K 6610TMB	SUB	TMB substrate (Tetramethylbenzidine), ready to use	2 x 15 ml
K 6610AC	STOP	ELISA stop solution, ready-to-use	1 x 15 ml

## 5. MATERIAL REQUIRED BUT NOT SUPPLIED

- Bidistilled water (aqua bidest.)
- Precision pipettors calibrated to deliver 10-200  $\mu$ l
- A multi-channel dispenser or repeating dispenser
- Centrifuge capable of 3000 x g
- Vortex-Mixer
- Microplate reader 450 nm

## 6. PREPARATION AND STORAGE OF REAGENTS

- To run the assay more than one time, please make sure that the reagents are carefully stored as mentioned. Prepare just the appropriate amount necessary for the assay.
- The **ELISA wash buffer concentrate** (WASHBUF) should be diluted with Aqua dest. **1:10** before use (add 450 ml A. dest. to 50 ml concentrate). Crystals could occur due to high salt concentration. The crystals have to be resuspended **before dilution of the buffer solutions** using a water bath (37°C). The buffer concentrates are stable at 2-8°C until the expiry date stated on the label. Diluted solutions could be stored at 2-8°C for 1 month.
- The reagent tubes with the lower amounts should be centrifuged before use, to guarantee that the required amount is available for use.
- The calibrators (STD) and control (CTRL) must be reconstituted with **250 $\mu$ l** of aqua bidest. Allow the vial to stand for 10 minutes and then mix thoroughly by gentle inversion to ensure complete reconstitution. Reconstituted calibrators and control are stable at -20°C until the expiry date stated on the label.
- The **Conjugate** (CONJ; POD antibody) has to be diluted **1:2000** in ELISA wash buffer (10  $\mu$ l POD antibody and 20 ml ELISA wash buffer). The undiluted antibody is stable at 2 -8 °C until expiry date given on the label. **Diluted antibody solution is not stable and could not be stored.**

## 7. PRECAUTIONS

- For *in vitro* diagnostic use only.
- The calibrators and controls contain human source material which was tested and found to be non-reactive to HBsAg, anti-HIV-1/2, and anti-HCV. Since no method can offer complete assurance that hepatitis B virus, HIV-1/2, HVC or other infectious agents are absent, these reagents should be handled as if potentially infectious.
- Stop Solution consists of diluted Sulfuric Acid. This is a strong acid. Although diluted, it still must be handled with care. It can cause burns and should be handled with gloves, eye protection, and appropriate protective clothing. Any spill should be wiped out immediately with copious quantities of water. Do not breathe vapor and avoid inhalation.
- Reagents should not be used beyond the expiration date shown on kit label.

## 8. SPECIMEN COLLECTION AND PREPARATION

### Plasma and serum

Plasma and sera must be diluted **1:2000** with sample dilution buffer (e.g. 10  $\mu$ l sample + 990  $\mu$ l sample dilution buffer. Take 40  $\mu$ l of this dilution and add 760  $\mu$ l dilution buffer to get the final dilution.)

Plasma and sera are stable at 2-8°C for about 14 days. For long time storage we recommend -20°C:

### Urine

Urine samples can be measured directly.

## 9. ASSAY PROCEDURE

### *Procedural notes*

- Do not interchange different lot numbers of any kit component within the same assay.
- The quality control guidelines should be observed.
- Incubation time, incubation temperature and pipetting volumes of the different components are defined by the producer. Any variations of the test procedure, that are not coordinated with the producer, may influence the test results. Immundiagnostik can therefore not be held reliable for any damage resulting from this.
- Carry out the assay with the actual manual delivered with the kit.

### *Test procedure*

Wash the precoated microtiter plate 5 x with 250  $\mu$ l ELISA wash buffer (WASHBUF). Carry out the tests in duplicate.


1. Add **200  $\mu$ l** 0.9% NACL (NaCl) solution into each well.
2. Add **10  $\mu$ l** of each STD (standard), CTRL (controls) and SAMPLE (patient samples).
3. Incubate for **1 hour**, shaking on a horizontal mixer, at room temperature.
4. Decant the content of the plate and wash the wells **5 x with 250 $\mu$ l** ELISA WASHBUF (wash buffer).
5. Add **200  $\mu$ l** prediluted CONJ (Conjugate).
6. Incubate for **1 hour**, shaking on a horizontal mixer, at room temperature.
7. Decant the content of the plate and wash the wells **5 x with 250 $\mu$ l** ELISA WASHBUF.
8. Add **200  $\mu$ l** SUB (TMB substrate) solution.
9. Incubate for **5-10 minutes** at room temperature.
10. Add **50  $\mu$ l** STOP (stop solution) and mix shortly.
11. Determine absorption with an ELISA reader at **450 nm** against 620 nm as reference. If no reference wavelength is available, read only at 450 nm. If the extinction of the highest standard exceeds the measurement range of the photometer, absorption must be measured immediately at 405 nm against 620 nm as reference.

## 10. RESULTS

A calibration curve is constructed from the standards. Commercially available software can be used as well as graph paper. Results of the samples are read from this calibration curve.

THE CALIBRATION CURVE IS NOT LINEAR, therefore a spline- or 4PL algorithm is recommended.

### *Typical calibration curve*

Concentration [mg/l]	5	1.66	0.56	0.19	0.07	0
OD mean value	1.711	1.021	0.501	0.260	0.160	0.148

The data is for demonstration only and cannot be used in place of data generations at the time of the assay.

## Plasma and Serum

The measured plasma and serum concentration have to be multiplied with 2000 to get the right concentration.

## 11. LIMITATIONS

Samples with  $\alpha_2$ -Macroglobulin levels greater than the highest calibrator, should be diluted and re-assayed.

## 12. QUALITY CONTROL

Immundiagnostik recommends commercial control samples for internal quality control.

Control samples should be analyzed with each run. Results, generated from the analysis of control samples, should be evaluated for acceptability using appropriate statistical methods. The results for the patient samples may not be valid, if within the same assay one or more values of the quality control sample are outside the acceptable limits.

### *Expected values*

Plasma and serum: 1.3 – 3.0 g/l

Urine: < 0,18 mg/l

We recommend each laboratory to establish its own norm concentration range.

## 13. PERFORMANCE CHARACTERISTICS

### *Precision and reproducibility*

The intra-assay variation of the Immundiagnostik  $\alpha_2$ -Makroglobulin ELISA test was calculated from 32 determinations on each one of two samples.

Intra-Assay CV n= 32

Sample	$\alpha_2$ -Makroglobulin mean value [mg/L]	Intra-Assay CV [%]
1	1125.53	6.20
2	676.92	7.96

The inter-assay variation of the Immundiagnostik  $\alpha_2$ -Makroglobulin ELISA test was calculated from data on 2 samples obtained in 10 different assays by three technicians on two different lots of reagents over a period of three months.

Inter-Assay CV n= 10

Sample	$\alpha_2$ -Makroglobulin mean value [mg/L]	Inter-Assay CV [%]
1	1937.12	9.23
2	1298.07	8.85

### Recovery

Two samples were spiked with three different  $\alpha_2$ -Makroglobulin-concentrations and measured with the assay.

Recovery n=3

Sample [ng/ml]	Spike [ng/ml]	$\alpha_2$ -Makroglobulin expected [ng/ml]	$\alpha_2$ -Makroglobulin measured [ng/ml]
344	521	865	736
344	415	759	730
344	293	637	585
176	521	697	571
176	415	591	540
176	293	469	420

### Sensitivity

The sensitivity was set as  $B_0 + 2SD$ . The zero-standard was measured 22 times.

Sample	$\alpha_2$ -Makroglobulin mean value [OD]	Standard variation	Detection limit [mg/L]
1	0.019	0.011	0.058

*Sample dilution*

Linearity n= 2

Two patient serum samples were diluted with sample dilution buffer. The results are shown below:

Sample	Dilution	Expected [mg/L]	Measured [mg/L]
A	1:2000	2528.3	2528.3
	1:4000	1264.2	1346.1
	1:8000	632.1	675.9
	1:16000	316.0	343.4
B	1:2000	2036.8	2036.8
	1:4000	1018.4	1046.6
	1:8000	509.2	509.2
	1:16000	254.6	285.4

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## 15. GENERAL NOTES ON THE TEST AND TEST PROCEDURE

- This assay was produced and put on the market according to the IVD guidelines of 98/79/EC.
- The test components which are made of human serum are tested for Australia antigen and HIV and found to be negative. However, since no test method can offer complete assurance that infectious agents are absent; these reagents should be handled as recommended for any potentially infectious human serum or blood specimen. The normal precautions for laboratory working should be observed.
- Reagents of the test package contain sodium azide as a bactericide. Contact with skin or mucous membranes has to be avoided.
- All reagents in the test package are to be used for in-vitro diagnostics only.
- The reagents should not be used after the date of expiry (see label on the test package).
- Single components with different lot numbers should not be mixed or exchanged.
- The guidelines for medical laboratories should be observed.
- Incubation time, incubation temperature and pipetting volumes of the different components are defined by the producer. Any variations of the test procedure, that are not coordinated with the producer, may influence the results of the test. Immundiagnostik can therefore not be held reliable for any damage resulting from

**Used symbols:**

Temperature limitation

Catalogue Number

In Vitro Diagnostic Medical Device

Contains sufficient for <n> tests

Manufacturer

Use by

Lot number