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# $\beta_2$ -Mikroglobulin ELISA Kit

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
# $\beta_2$ -Microglobulin ELISA Kit

*For the in vitro determination of  $\beta_2$ -Microglobulin in serum, plasma and urine*

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Table of content	Page
	2
<b>1. INTENDED USE</b>	<b>14</b>
<b>2. SUMMARY AND EXPLANATION OF THE TEST</b>	<b>14</b>
<b>3. PRINCIPLE OF THE TEST</b>	<b>14</b>
<b>4. MATERIAL SUPPLIED</b>	<b>15</b>
<b>5. MATERIAL REQUIRED BUT NOT SUPPLIED</b>	<b>15</b>
<b>6. PREPARATION AND STORAGE OF REAGENTS</b>	<b>16</b>
<b>7. PRECAUTIONS</b>	<b>17</b>
<b>8. SPECIMEN COLLECTION AND PREPARATION</b>	<b>17</b>
<b>9. ASSAY PROCEDURE</b>	<b>17</b>
PROCEDURAL NOTES	17
TEST PROCEDURE	18
<b>10. RESULTS</b>	<b>19</b>
TYPICAL CALIBRATION CURVE	19
<b>11. LIMITATIONS</b>	<b>20</b>
<b>12. QUALITY CONTROL</b>	<b>20</b>
NORMAL RANGE	20
<b>13. PERFORMANCE CHARACTERISTICS</b>	<b>20</b>
DETECTION LIMIT	21
LINEARITY	21
RECOVERY	21
PRECISION AND REPRODUCIBILITY	20
<b>14. GENERAL NOTES ON THE TEST AND TEST PROCEDURE</b>	<b>21</b>

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Manual


# $\beta_2$ -Microglobulin ELISA Kit

*For the in vitro determination of  $\beta_2$ -Microglobulin in serum,  
plasma and urine*

Gültig ab / Valid from 06.04.2010

K 6210

Immundiagnostik AG, Stubenwald-Allee 8a, D 64625 Bensheim  
Tel.: ++49 6251 70190-0  
Fax: ++ 49 6251 849430  
e.mail: [Info@immundiagnostik.com](mailto:Info@immundiagnostik.com)  
[www.Immundiagnostik.com](http://www.Immundiagnostik.com)

## 1. INTENDED USE

The *Immundiagnostik* Assay is intended for the quantitative determination of **beta-2- microglobulin** in serum, plasma and urine. For *in vitro* diagnostic use only.

## 2. SUMMARY AND EXPLANATION OF THE TEST

$\beta$ -2-Microglobulin is a light chain protein (11.8 kD) of the HLA-class-I antigens and is found on the cell membrane of all nucleated cells. This protein is metabolised extensively in the kidney. The serum concentration is influenced by the rates of synthesis and metabolism and is usually stable in healthy persons. Changes in the serum concentrations are indicative of disorders in glomerular and tubular functions.

Indications

- Early detection of a renal transplant rejection
- Assessment of the glomerular filtration rate (GFR)

## 3. PRINCIPLE OF THE TEST

This Enzyme-Linked Immunosorbent Assay (ELISA) allows the quantitative determination of  $\beta$ -2-microglobulin from plasma, serum and urine.

$\beta$ -2-microglobulin of the samples is bound to an excess of polyclonal rabbit anti  $\beta$ -2-microglobulin antibodies immobilised on the microtitre plate surface. After a washing step, to remove all unbound substances, the quantification of the bound  $\beta$ -2-microglobulin is carried out by adding an enzyme (horseradish peroxidase, POD) labeled antibody (POD-antibody), which also binds to the  $\beta$ -2-microglobulin. The amount of bound enzyme is directly proportional to the  $\beta$ -2-microglobulin content. The substrate tetramethylbenzidine (TMB) is converted by the enzyme to a chromogenic compound, which is determined photometrically at 450 nm (if the extinction is out of range, measure at 410 nm).

**4. MATERIAL SUPPLIED**

<b>Cat. No</b>	<b>Content</b>	<b>Kit Components</b>	<b>Quantity</b>
K 6210MTP	PLATE	One holder with precoated strips	96
K 6210WP	WASHBUF	ELISA wash buffer concentrate 10x	100 ml
K 6210PV	SAMPLEBUF	Sample dilution buffer, ready-to-use	1 x 100 ml
K 6210K	CONJ	Conjugate (rabbit anti $\beta$ -2-M, peroxidase-labelled), ready-to-use	2 x 15 ml
K 6210ST	STD	Calibrators, lyophilized (0; 0.6; 1.2; 2.5; 5; 10 mg/l)	6 x 1 vial
K 6210KO1	CTRL	Control, lyophilized	1 x 1 vial
K 6210KO2	CTRL	Control, lyophilized	1 x 1 vial
K 6210NaCl	NACL	0.9 % NaCl - solution, ready-to-use	25 ml
K 6210TMB	SUB	TMB substrate (Tetramethylbenzidine), ready-to-use	2 x 15 ml
K 6710AC	STOP	ELISA stop solution, ready to use	1 x 15 ml

**5. MATERIAL REQUIRED BUT NOT SUPPLIED**

- Bidistilled water (aqua bidest.)
- Laboratory balance
- Precision pipettors and disposable tips to deliver 10-1000  $\mu$ l
- Foil to cover the microtiter plate
- Horizontal microtiter plate shaker
- A multi-channel dispenser or repeating dispenser
- Centrifuge capable of 3000 x g
- Vortex-Mixer
- Standard laboratory glass or plastic vials, cups, etc.
- Microtiter plate reader at 450 or 405 nm (reference wave length 620 or 690 nm)

## 6. PREPARATION AND STORAGE OF REAGENTS

- To run assay more than once, ensure that reagents are stored at conditions stated on the label. **Prepare only the appropriate amount necessary for each assay.** The kit can be used up to 4 times within the expiry date stated on the label.
- Reagents with a volume less than **100 µl** should be centrifuged before use to avoid loss of volume.
- The **WASHBUF** (wash buffer concentrate) should be diluted with aqua bidest. **1:10** before use (100 ml WASHBUF + 900 ml aqua bidest.), mix well. Crystals could occur due to high salt concentration in the stock solutions. The crystals must be redissolved at 37°C in a water bath before dilution. The WASHBUF (wash buffer concentrate) is stable at **2-8°C** until the expiry date stated on the label. Diluted **buffer solution** can be stored in a closed flask at **2-8°C for one month.**
- The lyophilized **STD** (standards) and **CTRL** (controls) must be reconstituted with **250 µl** aqua bidest. Allow the vial content to dissolve for 10 minutes and mix thoroughly by gentle inversion to insure complete reconstitution. The reconstituted standards and controls must be diluted **1:50** with SAMPLEBUF (sample dilution buffer) (e.g. 20 µl sample + 980 µl SAMPLEBUF) before use. Diluted standards are not stable. Reconstituted standards are stable for 4 weeks at 2 – 8 ° C and for long term storage at –20 °C until the expiry date stated on the label.
- All other test reagents are ready to use. Test reagents are stable until the expiry date (see label of test package) when stored at **2-8°C.**

## 7. PRECAUTIONS

- For *in vitro* diagnostic use only.
- Human materials used in kit components were tested and found to be negative for HIV, Hepatitis B and Hepatitis C. However, for safety reasons, all kit components should be treated as potentially infectious.
- Kit reagents contain sodium azide or thimerosal as bactericides. Sodium azide and thimerosal are toxic. Substrates for the enzymatic color reactions are toxic and carcinogenic. Avoid contact with skin or mucous membranes.
- Stop solution is composed of sulfuric acid, which is a strong acid. Even diluted, it still must be handled with care. It can cause acid burns and should be handled with gloves, eye protection, and appropriate protective clothing. Any spills should be wiped out immediately with copious quantities of water.
- Reagents should not be used beyond the expiration date shown on the kit label.

## 8. SPECIMEN COLLECTION AND PREPARATION

All samples must be diluted **1:50** with SAMPLEBUF (sample dilution buffer) before use, for example 10 µl sample + 490 µl SAMPLEBUF.

Samples with a-2-microglobulin content higher than 10 mg/l should be further diluted 1:10 with SAMPLEBUF (sample dilution buffer).

### **Plasma or serum**

Samples can be stored for two weeks at 2-8 °C. For longer storage, samples should be frozen at -20 °C.

### **Urine**

Urine should be adjusted to a pH of 6 to 8 with 1 N NaOH. Adjusted samples can be stored at 2-8 °C for 14 days. For longer storage, non-treated samples should be frozen at -20 °C.

## 9. ASSAY PROCEDURE

### *Procedural notes*

- Do not mix different lot numbers of any kit component.
- Quality control guidelines should be followed.

- Incubation time, incubation temperature and pipetting volumes of the components are defined by the producer. Any variation of the test procedure, which is not coordinated with the producer, may influence the results of the test. Immundiagnostik AG can therefore not be held responsible for any damage resulting from wrong use.
- The assay should always be performed according the enclosed manual.

### *Test procedure*

Wash the precoated microtiter plate 5 x with 250 µl ELISA wash buffer. Carry out the tests in duplicate.


1. Pipette **200 µl** of **NACL** (0.9% NaCl solution) into each well of the microtitre plate.
2. Add **10 µl** of the diluted **STD** (standards), **CTRL** (controls) and **patient samples**.
3. Incubate for **1 hour**, shaking on a horizontal mixer, at room temperature.
4. Decant the contents of the plate and wash the cavities **5 x with 250 µl** of diluted washing buffer solution.
5. Add **200 µl** of **CONJ** (conjugate, POD-antibody) into each well.
6. Incubate for **15 minutes**, shaking on a horizontal mixer, at room temperature.
7. Decant the content of the plate and wash the cavities **5 x with 250 µl** of washing buffer solution.
8. Add **200 µl** of **SUB** (TMB-solution) into each well.
9. Incubate for **5-15 minutes** at room temperature, shaking slightly, until color differences are sufficient.
10. Add **50 µl** of **STOP** (stop solution) and mix shortly.
11. Determine **absorption immediately** with an ELISA reader at **450 nm** against 620 nm (or 690 nm) as reference. If no reference wavelength is available, read only at 450 nm. If the extinction of the highest standard exceeds the measurement range of the photometer, absorption must be measured immediately at 405 nm against 620 nm (or 690 nm) as reference.

## 10. RESULTS

A calibration curve is constructed from the calibrators. Commercially available software can be used as well as graph paper. Results of the samples are read from this calibration curve.

THE CALIBRATION CURVE IS NOT LINEAR, therefore a spline- or 4PL algorithm is recommended.

### *Typical calibration curve*

<b>Concentration [mg/l]</b>	10	5	2.5	1.2	0.6	0
<b>OD mean value</b>	1.834	1.429	0.990	0.642	0.412	0.168

The data is for demonstration only and cannot be used for the evaluation of test results.

## Evaluation

The results are **not** multiplied by any factor at 1:50 sample dilution.

If samples are diluted more than 1:50 (at concentrations higher than 10 mg/ml), the corresponding dilution factor must be considered.

## 11. LIMITATIONS

Samples with levels greater than the highest standard should be further diluted and re-assayed.

## 12. QUALITY CONTROL

**Immundiagnostik AG recommends the use of commercial control samples for internal quality control if available.**

Control samples should be analyzed with each run. Results, generated from the analysis of control samples, should be evaluated for acceptability using appropriate statistical methods. The results for the patient samples may not be valid, if within the same assay one or more values of the quality control sample are outside the acceptable limits.

### *Reference range*

Plasma or serum: < 2.5 mg/l

Urine: < 0.4 mg/l

## 13. PERFORMANCE CHARACTERISTICS

### *Precision and reproducibility*

Through repeated measurements (n=12) of plasma and urine samples, containing  $\beta$ -2-microglobulin, the following results were obtained:

Intra-assay n = 12:

Sample	$\beta$ -2-M mean [mg/l]	Intra-Assay CV [%]
Urine	5.7	9
Plasma	1.1	11

Inter-assay n = 12:

Sample	$\beta$ -2-M mean [mg/l]	Inter-Assay CV [%]
Urine	5.9	15
Plasma	1.0	12

### Recovery

Recovery studies were made with urine and plasma samples. The mean value of the plasma and urine samples was 98% and 93%, respectively.

### Detection limit

The detection limit was defined as  $B_0 + 2 \text{ SD}$  and set to 0,1 mg/l.

### Linearity

The linearity of this assay was estimated by dilution with sample buffer of  $\beta$ -2-microglobulin-containing material. The linearity extends from 0,2 to 10 mg/L

## 14. GENERAL NOTES ON THE TEST AND TEST PROCEDURE

- For *in vitro* diagnostic use only.
- Quality control guidelines should be followed.
- Human materials used in kit components were tested and found to be negative for HIV, Hepatitis B and Hepatitis C. However, for safety reasons, all kit components should be treated as potentially infectious.
- Kit reagents contain sodium azide or thimerosal as bactericides. Sodium azide and thimerosal are toxic. Substrates for the enzymatic color reactions are toxic and carcinogenic. Avoid contact with skin or mucous membranes.
- Stop solution is composed of sulfuric acid, which is a strong acid. Even diluted, it still must be handled with care. It can cause acid burns and should be handled with gloves, eye protection, and appropriate protective clothing. Any spills should be wiped out immediately with copious quantities of water.
- Do not mix different lot numbers of any kit component.

- Reagents should not be used beyond the expiration date shown on the kit label.
- The assay should always be performed according the enclosed manual.
- Incubation time, incubation temperature and pipetting volumes of the components are defined by the producer. Any variation of the test procedure, which is not coordinated with the producer, may influence the results of the test. Immundiagnostik AG can therefore not be held responsible for any damage resulting from wrong use.
- Warranty claims and complaints in respect of deficiencies must be logged within 14 days after receipt of the product. The product should be send to Immundiagnostik AG along with a written complaint.

**Used symbols:**

Temperature limitation

Catalogue Number

In Vitro Diagnostic Medical Device

Contains sufficient for &lt;n&gt; tests

Manufacturer

Use by

Lot number