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β-Defensin 2 ELISA Kit

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TABLE OF CONTENT	PAGE
1. INTENDED USE	14
2. SUMMARY AND EXPLANATION OF THE TEST	14
3. TEST PRINCIPLE	14
4. MATERIAL SUPPLIED	15
5. MATERIAL REQUIRED BUT NOT SUPPLIED	15
6. PREPARATION AND STORAGE OF REAGENTS	16
7. PRECAUTIONS	16
8. SPECIMEN COLLECTION AND PREPARATION	17
STOOL EXTRACTION	17
DILUTION OF SAMPLES	19
9. ASSAY PROCEDURE	19
PROCEDURAL NOTES	19
TEST PROCEDURE	20
10. RESULTS	21
11. LIMITATIONS	22
12. QUALITY CONTROL	22
EXPECTED VALUES	22
13. TEST CHARACTERISTICS	23
SENSITIVITY	23
INTER AND INTRA VARIATION COEFFICIENT	23
14. PUBLICATIONS WITH THE USE OF THE ID β -DEFENSIN ELISA KIT	23
15. GENERAL NOTES ON THE TEST AND TEST PROCEDURE	23

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Manual


β-Defensin 2 ELISA Kit

For the in vitro determination of β-Defensin 2 in stool

Valid from 25.01.2010

K 6500

Immundiagnostik AG, Stubenwald-Allee 8a, D 64625 Bensheim
Tel.: ++49 6251 70190-0
Fax: ++ 49 6251 849430
e.mail: Info@immundiagnostik.com
www.immundiagnostik.com

1. INTENDED USE

This assay is intended for the quantitative determination of β -Defensin 2 in stool. For *in vitro* diagnostic use only.

2. SUMMARY AND EXPLANATION OF THE TEST

The β -defensins are an integral part of the congenital immune system and contribute with their antimicrobial effect to the barrier function of intestinal epithelial cells.

Defensins exert a variable degree of antimicrobial activity against bacteria, fungi, and some enveloped viruses. Vertebrate defensins are classified as alpha- or beta-defensins, based on their pattern of disulfide bridges. Nine human defensins of epithelial origin have been found, three of them being β -defensins (HBD-1, -2 and -3). The expression of β -defensins is induced by the pro-inflammatory cytokines and also through microorganisms (e.g. *E. coli*, *H. pylori* or *P. aeruginosa*).

A β -defensin-2-deficiency can, for example, be observed in the intestinal mucous of patients with Crohn's disease. The defense system of the mucous membrane is therefore restricted and allows an increased invasion of bacteria, which could possibly lead to a typical infection in Crohn's disease patients.

Whether the β -defensin-2 deficiency could even play a role in the development of Crohn's disease is currently being researched. As is the possibility that it is the probiotic bacterium, which produces β -defensin.

Indications:

- Reduced β -defensin levels with Crohn's disease (HBD-2)
- Increased β -defensin levels with Colitis Ulcerosa (HBD-2)

3. TEST PRINCIPLE

The β -defensin 2 in standard and samples is bound to an available excess of polyclonal antibodies against β -defensin 2, which are immobilized on the surface of the microtiter plate. After a washing step, to remove all interfering substances, the quantification of bound β -defensin 2 is carried out by adding a polyclonal anti β -defensin 2 antibody, which is biotinylated. This antibody is detected with horseradish peroxidase labeled streptavidin. The amount of converted substrate by Peroxidase is directly proportional to the amount of bound β -defensin 2 and can be determined photometrically at 450 nm.

4. MATERIAL SUPPLIED

Catal. No.	Content	Kit Components	Quantity
K 6500MTP	PLATE	One holder with precoated strips	12 x 8 wells
K 6500WP	WASHBUF	ELISA wash buffer concentrate 10x	2 x 100 ml
K 6500A2	AB	Detection antibody (Goat anti β -Defensin 2, biotinylated)	200 μ l
K 6500K	CONJ	Conjugate, (Streptavidin, HRP-conjugated)	200 μ l
K 6500ST	STD	Calibrator, lyophilized	4 x 5 vials x 500 μ l
K 6500KO	CTRL	Control, lyophilized	4 x 1 vial x 500 μ l
K 6500TMB	SUB	TMB substrate (Tetramethylbenzidine), ready to use	15 ml
K 6500AC	STOP	ELISA stop solution, ready to use	15 ml
K 6500VP	ABBUF	Detection antibody dilution buffer	15 ml

5. MATERIAL REQUIRED BUT NOT SUPPLIED

- Bidistilled water (aqua bidest.)
- Laboratory balance
- Precision pipettors and disposable tips to deliver 10-1000 μ l
- Foil to cover the microtiter plate
- Horizontal microtiter plate shaker
- A multi-channel dispenser or repeating dispenser
- Centrifuge capable of 3000 x g
- Vortex-Mixer
- Standard laboratory glass or plastic vials, cups, etc.
- Microtiter plate reader at 450 or 405 nm (reference wave length 620 or 690 nm)

6. PREPARATION AND STORAGE OF REAGENTS

- To run assay more than once, ensure that reagents are stored at conditions stated on the label. **Prepare only the appropriate amount necessary for each assay.** The kit can be used up to 4 times within the expiry date stated on the label.
- Reagents with a volume less than **100 µl** should be centrifuged before use to avoid loss of volume.
- The **WASHBUF** (wash buffer concentrate) should be diluted with aqua bidest. **1:10** before use (100 ml concentrate + 900 ml aqua bidest.), mix well. Crystals could occur due to high salt concentration in the stock solutions. The crystals must be redissolved at room temperature or at 37°C before dilution. The **WASHBUF** (wash buffer concentrate) is stable at **2-8°C** until the expiry date stated on the label. Diluted **buffer solution** can be stored in a closed flask at **2-8°C for one month.**
- The lyophilized **STD** (standards) and **CTRL** (control) must be reconstituted with **500 µl** aqua bidest. Allow the vial content to dissolve for 10 minutes and mix thoroughly by gentle inversion to insure complete reconstitution. Reconstituted standards and control are **not stable.**
- The **AB** (detection antibody, biotinylated) must be diluted **1:100** in **ABBUF** (detection antibody dilution buffer; e.g. 100 µl AB + 10 ml ABBUF). The undiluted AB is stable at -20°C until the expiry date given on the label. **Diluted antibody solution is not stable and can not be stored.**
- The **CONJ** (conjugate, POD-antibody) must be diluted **1:100** in wash buffer (100 µl CONJ + 10 ml wash buffer). The undiluted CONJ is stable at **2-8 °C** until the expiry date stated on the label. **Diluted conjugate is not stable and can not be stored.** All other test reagents are ready to use. Test reagents are stable until the expiry date (see label of test package) when stored at **2-8°C.**
- All other test reagents are ready to use. Test reagents are stable until the expiry date (see label of test package) when stored at **2-8°C.**

7. PRECAUTIONS

- For *in vitro* diagnostic use only.
- Human materials used in kit components were tested and found to be negative for HIV, Hepatitis B and Hepatitis C. However, for safety reasons, all kit components should be treated as potentially infectious.

- Kit reagents contain sodium azide or thimerosal as bactericides. Sodium azide and thimerosal are toxic. Substrates for the enzymatic color reactions are toxic and carcinogenic. Avoid contact with skin or mucous membranes.
- Stop solution is composed of sulfuric acid, which is a strong acid. Even diluted, it still must be handled with care. It can cause acid burns and should be handled with gloves, eye protection, and appropriate protective clothing. Any spills should be wiped out immediately with copious quantities of water.
- Reagents should not be used beyond the expiration date shown on the kit label.

8. SPECIMEN COLLECTION AND PREPARATION

Stool Extraction

Please use the **washing buffer** as a sample extraction buffer.

1. We recommend the use of a stool sample preparation kit for dosing of 100 mg of stool sample (e.g. Sample preparation kit from Roche Diagnostics, Mannheim, Germany; cat # 10 745804332). Add 5 ml buffer to the stool sample in the preparation kit.

Constant buffer volume: 5 ml

Constant dilution factor: 1 : 50

2. Alternatively, stool samples can be manually weighted within the range of 80 – 120 mg. Please note the exact sample amount!
 - a. Add **5 ml** buffer to the stool sample independent of the sample amount.

Constant buffer volume: 5 ml

Variable dilution factor

The dilution factor varies depending on the sample amount which must be considered in the subsequent calculations. Please see the list of the correction factors below:

Weight [mg]	Dilution factor
80	62.5
82	60.9
84	59.5
86	58.1
88	56.8
90	55.6
92	54.3
94	53.2
96	52.1
98	51.0
100	50

Weight [mg]	Dilution factor
102	49.0
104	48.1
106	47.2
108	46.3
110	45.5
112	44.6
114	43.9
116	43.1
118	42.4
120	41.6

- b. The **buffer volume** for the individual samples **varies** depending on the sample amount (see table).

Variable buffer volume

Constant dilution factor: 1 : 50

Therefore, the same dilution factor can be used for all samples in the subsequent evaluation of the results .

Weight [mg]	Buffer Volume [ml]
80	4.0
82	4.1
84	4.2
86	4.3
88	4.4
90	4.5
92	4.6
94	4.7
96	4.8
98	4.9
100	5.0

Weight [mg]	Buffer Volume [ml]
102	5.1
104	5.2
106	5.3
108	5.4
110	5.5
112	5.6
114	5.7
116	5.8
118	5.9
120	6.0

Afterwards, mix the weighed stool sample with the buffer and vortex for 30 sec. Transfer ca. 1 ml stool suspension in an Eppendorf-cup and centrifuge for 5 min at 13000 rpm.

The supernatant is not stable and can not be stored. We recommend to weight fresh sample amount for a new assay, if the analysis should be repeated.

Dilution of samples

Stool samples

Dilute the supernatant of the extraction (**dilution step I**) 1:2 with **wash buffer**, for example:

300 µl supernatant (dilution I) + 300 µl wash buffer = 1:2 (**dilution step II**)

For analysis, pipet **100 µl** of the supernatant of **dilution step II** per well.

The supernatant is not stable and can not be stored. We recommend to weight fresh sample amount for a new assay, if the analysis should be repeated.

9. ASSAY PROCEDURE

Procedural notes

- Do not mix different lot numbers of any kit component.
- Quality control guidelines should be followed.
- Incubation time, incubation temperature and pipetting volumes of the components are defined by the producer. Any variation of the test procedure, which is not coordinated with the producer, may influence the results of the test. Immundiagnostik AG can therefore not be held responsible for any damage resulting from wrong use.
- The assay should always be performed according the enclosed manual.

Test procedure

We recommend to carry out the tests in duplicate.

1. Wash the precoated microtiter plate **5 x with 250 µl** ELISA wash buffer.
2. Pipette **100 µl** of **STD** (standards), **CTRL** (control) or **samples** into each well.
3. Incubate for **1 hour** at room temperature, shaking on a horizontal mixer.
4. Decant the contents of the plate and wash the cavities **5 x with 250 µl** of washing buffer solution.
5. Add **100 µl** diluted **AB** (detection antibody solution).
6. Incubate for **1 hour** at room temperature, shaking on a horizontal mixer.
7. Decant the contents of the plate and wash the cavities **5 x with 250 µl** of washing buffer solution.
8. Add **100 µl** diluted **CONJ** (conjugate solution).
9. Incubate for **1 hour** at room temperature, shaking on a horizontal mixer.
10. Decant the contents of the plate and wash the cavities **5x with 250 µl** of washing buffer solution.
11. Add **100 µl** of **SUB** (TMB-substrate) solution
12. Incubate approximately for **15-20 minutes** at room temperature, shaking slightly, until sufficient coloring is achieved.
13. Add **50 µl** of **STOP** (stop solution) and mix shortly.
14. Determine **absorption immediately** with an ELISA reader at **450 nm** against 620 nm (or 690 nm) as reference. If no reference wavelength is available, read only at 450 nm. If the extinction of the highest standard exceeds the measurement range of the photometer, absorption must be measured immediately at 405 nm against 620 nm (or 690 nm) as reference.

10. RESULTS

The following algorithms can be used alternatively to calculate the results. We recommend to use the "4-Parameter-algorithm".

1. 4-parameter-algorithm

It is recommended a linear ordinate for optical density and a logarithmic abscissa for concentration. When using a logarithmic abscissa, the zero calibrator must be specified with a value less than 1 (e. g. 0.01).

2. Point-to-point-calculation

We recommend a linear ordinate for optical density and a linear abscissa for concentration.

3. Spline-algorithm

We recommend a linear ordinate for optical density and a logarithmic abscissa for concentration. When using a logarithmic abscissa, the zero calibrator must be specified with a value less than 1 (e. g. 0.01).

The plausibility of the pairs of values should be examined before the automatic evaluation of the results. If this option is not available with the used program, a control of the paired values should be done manually.

Faeces

To obtain the β -defensin 2 concentration in fecal samples multiply the estimated value by the dilution factor according to the sample preparation:

Sample preparation 1 and 2b: **Constant dilution factor : 50**

Multiply the obtained result by **100** to get the real concentration.

Sample preparation 2a: **Variable dilution factor**

The corresponding dilution factor is taken from the table and should be multiplied by 2 (dilution step II) to get the final concentration.

11. LIMITATIONS

Samples with β -defensin-2 levels greater than the highest calibrator, should be diluted and re-assayed.

12. QUALITY CONTROL

Immundiagnostik recommends commercial control samples for internal quality control.

Control samples or serum pools should be analyzed with each run of calibrators and patient samples. Results, generated from the analysis of control samples, should be evaluated for acceptability using appropriate statistical methods. In assays in which one or more of the quality control sample values lie outside the acceptable limits, the results for the patient sample may not be valid.

Expected values

Reference range

(1 g stool is equivalent to 1ml)

Stool (n = 101): **35 ng/ml**

Based on Immundiagnostik studies of evidently healthy persons (n = 101) a mean value of 35 ng/ml stool was estimated. This value is consistent with the results published by Langhorst et al. (2007) obtained using the β -Defensin ELISA Kit of Immundiagnostik.

Stool (n = 23 healthy controls): **31.0 \pm 15.4 ng/g Stool***

*Langhorst J, Wieder A, Rueffer A, Michalsen A, Musial F, Dobos GJ. (2007) Activated innate immune system in irritable bowel syndrome? *Gut Online First*, published on April 11, 2007 as 10.1136/gut.2007.125005

We recommend each laboratory to establish its own reference range.

13. TEST CHARACTERISTICS

*Sensitivity***

Detection limit 0.077 ng/ml.

*Inter and Intra variation coefficient***

Inter variation coefficient - 5.11 %

Intra variation coefficient - 5.26 %

**Soto E, Espinoza J, Nien JK, Kusanovic JP, Erez O, Richani K, Santolaya-Forgas J, Romero R. (2007) Human beta-defensin-2: A natural antimicrobial peptide present in amniotic fluid participates in the host response to microbial invasion of the amniotic cavity. *J Matern Fetal Neonatal Med.* Jan; 20(1):15-22.

14. PUBLICATIONS WITH THE USE OF THE K 6500 ELISA KIT

1. Langhorst J, Wieder A, Rueffer A, Michalsen A, Musial F, Dobos GJ (2007) Activated innate immune system in irritable bowel syndrome? *Gut* Apr 11, 2007 [Epub ahead of print]
2. Soto E, Espinoza J, Nien JK, Kusanovic JP, Erez O, Richani K, Santolaya-Forgas J, Romero R. (2007) Human beta-defensin-2: A natural antimicrobial peptide present in amniotic fluid participates in the host response to microbial invasion of the amniotic cavity. *J Matern Fetal Neonatal Med.* Jan; 20(1):15-22.
3. Langhorst J, Elsenbruch S, Rueffer A, Michalsen A, Dobos GJ (2006) Higher Human Beta-Defensin 2 levels in feces in active Ulcerative Colitis and Ileal but not Colonic Crohn's disease compared to IBS. Suppl 2 to *J Gastroenterol* (4):Abstract A205, S1340, *107th Annual Meeting of AGA*, May 20-25, 2006, Los Angeles
4. Kapel N, Benahmed N, Morali A, Svahn J, Canioni D, Goulet O, Ruemmele FM (2009) Fecal beta-Defensin-2 in Children With Inflammatory Bowel Diseases. *J Pediatric Gastroenterol and Nutr* 48:117–120.

15. GENERAL NOTES ON THE TEST AND TEST PROCEDURE

- This assay was produced and put on the market according to the IVD guidelines of 98/79/EC.
- The test components which are made of human serum are tested for HVB and HIV and found to be negative. However, since no test method can

offer complete assurance that infectious agents are absent, these reagents should be handled as recommended for any potentially infectious human serum or blood specimen. The normal precautions for laboratory working should be observed.

- Reagents of the test package contain sodium azide as a bactericide. Contact with skin or mucous membranes has to be avoided.
- All reagents in the test package are to be used for in-vitro diagnostics only.
- The reagents should not be used after the date of expiry (see label on the test package).
- Single components with different lot numbers should not be mixed or exchanged.
- The guidelines for medical laboratories should be observed.
- Incubation time, incubation temperature and pipetting volumes of the different components have been defined by the producer. Any alterations of the test procedure, that are not coordinated with the producer, may influence the results of the test. Immundiagnostik can therefore not be held responsible for any damage.

Used symbols:

Temperature limitation

Catalogue Number

In Vitro Diagnostic Medical Device

Contains sufficient for <n> tests

Manufacturer

Use by

Lot number