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# Thyreoglobulin ELISA KIT

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
# Thyroglobulin ELISA KIT

*For the in vitro Determination of Thyroglobulin in Serum and Plasma*

Valid from 13.02.2008

K 7510

## 1. INTENDED USE

The described Enzyme-Linked-Immuno-Sorbent-Assay (ELISA) is intended for the quantitative determination of **human thyroglobulin (hTG)** in plasma and serum. It is for *in vitro* diagnostic use only.

## 2. INTRODUCTION

**Human Thyroglobulin (hTG)** is a large glycoprotein (MW 660.000 Da) that is stored in the follicular colloid of the thyroid gland. It functions as a prohormone in the intrathyroid synthesis of T<sub>3</sub> and T<sub>4</sub>. Traditionally circulating levels of **hTG** have been measured using double antibody immunoassays. **HTG** is evaluated in thyroid follicular and papillary carcinoma, thyroid adenoma, subacute thyroiditis, Hashimoto's thyroiditis and Grave's disease. **HTG** levels are not increased in patients with medullary thyroid carcinoma. Low **hTG** concentrations are an indication that thyrotoxicosis factitia may be present. Measurement of **hTG** is most useful in detecting recurrence of differentiated thyroid carcinoma following surgical resection or radioactive iodine ablation. **HTG** determination is used as an adjunct to iodine 131 scanning but not as a replacement for it. Assessment of serum **hTG** also aids in the management of infants with congenital hypothyroidism.

### Indications:

- Postoperative monitoring of differentiated thyroid carcinoma
- Congenital hypothyroidism
- Hyperthyreosis factitia
- Tumor recurrence in patients with thyroid carcinoma

### 3 MATERIAL SUPPLIED

Catalogue No	Content	Kit Components	Quantity
K 7510MTP	PLATE	One holder with precoated strips	12 x 8 wells
K 7510WP	WASHBUF	ELISA wash buffer concentrate 10x	1 x 100 ml
K 7510K	CONJ	Conjugate, rabbit-anti-hTG antibodies, peroxidase-labeled, ready to use	1 x 15 ml
K 7510KO	CTRL	Control, ready-to-use	2 vials
K 7510ST	STD	Standards, ready to use (0; 4; 8; 16; 31; 62; 125; 250 µg/l)	2 x 8 vials
K 7510TMB	SUB	Tetramethylbenzidin (TMB) substrate solution, ready to use	1 x 15 ml
K 7510AC	STOP	ELISA stop solution, ready to use	1 x 15 ml
K 7510AP	ASYBUF	Assay buffer, ready to use	1 x 10 ml

### 4 MATERIAL REQUIRED BUT NOT SUPPLIED

- Bidistilled water (aqua bidest.)
- Laboratory balance
- Precision pipettors and disposable tips to deliver 10-1000 µl
- Foil to cover the microtiter plate
- Horizontal microtiter plate shaker
- A multi-channel dispenser or repeating dispenser
- Centrifuge capable of 3000 x g
- Vortex-Mixer
- Standard laboratory glass or plastic vials, cups, etc.
- Microtiter plate reader at 450 nm  
(reference wave length 620 or 690 nm)

## 5. PREPARATION AND STORAGE OF REAGENTS

- To run the assay more than one time, make sure that the reagents are stored at the conditions stated on the label. **Prepare just the appropriate amount necessary for the assay.** The kit can be used up to 4 times within the expiry date stated on the label.
- Reagents with a volume less than **100 µl** should be centrifuged before use to avoid loss of volume.
- The **WASHBUF** (ELISA wash buffer concentrate) should be diluted with aqua bidest. **1:10** before use (100 ml WASHBUF + 900 ml aqua bidest.), mix well. Crystals could occur due to high salt concentration in the stock solutions. The crystals have to be redissolved at 37°C using a water bath before dilution. The **WASHBUF** is stable at **2-8°C** until the expiry date stated on the label. Diluted **buffer solution** can be stored in a closed flask at **2-8°C for one month**.
- The **STD** (standards) and the **CTRL** (control) can be stored at 2-8°C for 4 weeks. Long term storage until the expiry date given on the label at – 20°C only.
- All other test reagents are ready to use. The test reagents are stable until the expiry date (see label of test package) when stored at **2-8°C**.

## 6. SAMPLE PREPARATION

### Serum and plasma

Fresh collected blood should be centrifuged within one hour. Store samples at -20 °C if not assayed on the same day. Lipemic or hemolytic samples may give erroneous results. Samples should be mixed well before assaying. We recommend duplicate analyses for each sample.

The **serum** and **plasma** samples should be diluted **1:10** with **ASYBUF** (assay buffer) if their hTG concentration is higher than 250 µg/L.

For example:

**50 µl** sample + **450 µl ASYBUF** (assay buffer)

## 7. ASSAY PROCEDURE

### *Principle of the test*

The assay utilizes the "sandwich" technique with two selected polyclonal antibodies that bind to human Thyroglobulin.

Standards, controls and patient samples which are assayed for human thyroglobulin are added into the wells of a microplate coated with a high affine polyclonal anti-human thyroglobulin antibody. During the first incubation step, thyroglobulin is bound by the immobilized antibody. Then a peroxidase-conjugated polyclonal rabbit anti human thyroglobulin antibody is added into each microtiter well and a "sandwich" of capture antibody - human thyroglobulin – peroxidase-conjugate is formed. Tetra-methylbenzidine (TMB) is used as peroxidase substrate. Finally, an acidic stop solution is added to terminate the reaction. The color changes from blue to yellow. The intensity of the yellow color is directly proportional to the concentration of thyroglobulin. A dose response curve of the absorbance unit (optical density, OD at 450 nm) vs. concentration is generated, using the values obtained from the standard. Thyroglobulin present in the patient samples is determined directly from this curve.

### *Test procedure*

Prior to use in the assay allow all reagents and samples to come to room temperature (18-26 °C) and mix well. **Perform analysis always in duplicate**

Mark the positions of **STD** (standards), **CTRL** (control) or **SAMPLE** (samples) on a protocol sheet

Take microtiter strips out of the kit. Store unused strips covered at 2-8° C. Strips are stable until the expiry date stated on the label

Wash each well **5 times by dispensing 250 µl of diluted wash buffer** into each well. After the final washing step remove residual buffer by tapping the plate on absorbent paper

Add **100 µl** of **STD** (standards), **CTRL** (control) or **SAMPLE** (sample) in duplicate into respective well

Cover the plate tightly and incubate for **3 hours** at **37°C** shaking on a horizontal mixer or **4 hours** at **37°C** without shaking. Alternatively, incubate **overnight at room temperature** on a horizontal mixer.

Discard the contents of each well. Wash **5 times by dispensing 250 µl of diluted wash buffer** into each well. After the final washing step remove residual buffer by tapping the plate on absorbent paper

Add **100 µl of CONJ** (conjugate) into each well

Cover the plate tightly and **incubate for 1 hour at room temperature (18-26°C)** on a horizontal mixer

Discard the contents of each well. Wash **5 times by dispensing 250 µl of diluted wash buffer** into each well. After the final washing step remove residual buffer by tapping the plate on absorbent paper

Add **100 µl of SUB** (substrate) into each well

Incubate for **10 - 20 minutes at room temperature (18-26°C)** in the dark\*

Add **50 µl of STOP** (stop solution) into each well, mix thoroughly

Determine absorption immediately with an ELISA reader at **450 nm** against 620 nm (or 690 nm) as a reference. If no reference wavelength is available, read only at 450 nm. If the extinction of the highest standard exceeds the range of the photometer, absorption must be measured immediately at 405 nm against 620 nm as a reference

\*The intensity of the color change is temperature sensitive. We recommend to observe the color change and to stop the reaction upon good differentiation.

### **Note**

The assay can be run automatically, e. g. using the ETILAB instrument of Sorin/Biomedica.

## 8. RESULTS

The following algorithms can be used alternatively to calculate the results. We recommend to use the "4-Parameter-algorithm".

1. 4-Parameter-algorithm

It is recommended to use a linear ordinate for the optical density and a logarithmic abscissa for the concentration. When using a logarithmic abscissa, the zero calibrator has to be specified with a value smaller than 1 (e. g. 0.01).

2. Point-to-point-calculation

We recommend a linear ordinate for the optical density and a linear abscissa for the concentration.

3. Spline-algorithm

We recommend for the optical density a linear ordinate and for the concentration a logarithmic abscissa. When using a logarithmic abscissa, the zero calibrator has to be specified with a value smaller than 1 (e. g. 0.01).

THE CALIBRATION CURVE IS NOT LINEAR, therefore a point-to-point-calculation - algorithm is recommended.

The plausibility of the pairs of values should be examined before the automatic evaluation of the results. If this option is not available with the used program, a control of the paired values should be done manually.

### Serum and plasma

For the calculation of the thyroglobulin concentration in diluted **plasma and serum** samples the result has to be multiplied by **10**.

## 9. QUALITY CONTROL

Control samples should be analyzed with each run. Results, generated from the analysis of control samples, should be evaluated for acceptability using appropriate statistical methods. The results for the patient samples may not be valid, if within the same assay one or more values of the quality control sample are outside the acceptable limits.

### *Expected values*

#### **Normal ranges**

Serum/plasma: < 50µg/l

We recommend each laboratory to establish its own norm concentration range.

## 10. PERFORMANCE CHARACTERISTICS

### *Sensitivity*

The detection limit of this hTG ELISA was determined to  $B_0 + 2SD$ . The limit is 1 µg/l.

### *Cross reactivity*

The described hTG assay is highly specific. No cross reactivity with other proteins in plasma and serum was observed.

### *Linearity*

The linearity of the assay was determined by serial dilutions of material with known hTG concentration with the assay buffer. The linearity is in the range of 2.5 - 250 µg/l.

### *Precision data*

Intra- (CV%) and inter-assay (CV%) variations were determined using hTG containing samples within different concentration ranges. The intra-assay CV was 10.2 %, whereas 14.8 % was obtained for the inter-assay CV using three different kit batches.

## 11. PRECAUTIONS

- For *in vitro* diagnostic use only.
- The quality control guidelines should be followed.
- Human material used in the kit components was tested and found to be negative for HIV, Hepatitis B and Hepatitis C. However, for safety reasons, all kit components should be treated as potentially infectious.
- Reagents of the kit package contain sodium azide or thimerosal as bactericides. Sodium azide and thimerosal are toxic. The substrates for the enzymatic color reactions are toxic and carcinogenic. Avoid contact with skin or mucous membranes.
- Stop solution is composed of sulfuric acid, which is a strong acid. Even diluted, it still must be handled with care. It can cause acid burns and should be handled with gloves, eye protection, and appropriate protective clothing. Any spill should be wiped out immediately with copious quantities of water.

## 12. TECHNICAL HINTS

- Do not mix different lot numbers of any kit component.
- Reagents should not be used beyond the expiration date shown on the kit label.
- Substrate solution should remain colourless until use.
- To ensure accurate results, proper adhesion of plate sealers during incubation steps is necessary.
- Avoid foaming when mixing reagents.
- The assay should always be performed according the enclosed manual.

### 13. GENERAL NOTES ON THE TEST AND TEST PROCEDURE

- This assay was produced and distributed according to the IVD guidelines of 98/79/EC.
- All reagents in the kit package are for *in vitro* diagnostic use only.
- The guidelines for medical laboratories should be followed.
- Incubation time, incubation temperature and pipetting volumes of the components are defined by the producer. Any variation of the test procedure, not in accordance with manufacturer's recommendations, may influence the results of the test. Therefore, Immundiagnostik AG accepts no responsibility for experimental data resulting from any modifications or incorrect use of the product.
- Warranty claims and complaints in respect of deficiencies must be lodged within 14 days after receipt of the product. The product should be send to Immundiagnostik AG together with a written complaint.

13.02.2008 13.02.2008\_hTG.DOC

#### Used symbols:

Temperature limitation

Catalogue Number

In Vitro Diagnostic Medical Device

Contains sufficient for <n> tests

Manufacturer

Use by

Lot number