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Super sensitive

Präalbumin/Transthyretin ELISA Kit

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Table of contents	Page
1. INTENDED USE	14
2. SUMMARY AND EXPLANATION OF THE TEST	14
3. PRINCIPLE OF THE TEST	14
4. MATERIAL SUPPLIED	15
5. MATERIAL REQUIRED BUT NOT SUPPLIED	15
6. PREPARATION AND STORAGE OF REAGENTS	16
7. PRECAUTIONS	17
8. SPECIMEN COLLECTION AND PREPARATION	17
9. ASSAY PROCEDURE	18
PROCEDURAL NOTES	18
TEST PROCEDURE	19
10. RESULTS	20
TYPICAL CALIBRATION CURVE	20
11. LIMITATIONS	21
12. QUALITY CONTROL	21
EXPECTED VALUES	21
13. PERFORMANCE CHARACTERISTICS	21
14. REFERENCES	22
15. GENERAL NOTES ON THE TEST AND TEST PROCEDURE	22

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Super sensitive


Prealbumin/Transthyretin ELISA Kit

*For the in vitro determination of prealbumin/transthyretin in
serum, urine and stool*

Gültig ab / Valid from 17.12.2007

K 6331

1. INTENDED USE

The Prealbumin Enzyme-Linked Immunosorbent Assay (ELISA) kit is intended for the quantitative determination of **prealbumin (transthyretin)** in serum, urine and stool. It is for *in vitro* diagnostic use only.

2. SUMMARY AND EXPLANATION OF THE TEST

soon available

3. PRINCIPLE OF THE TEST

This Enzyme-Linked Immunosorbent Assay (ELISA) is a two step assay for the ultrasensitive determination of human **prealbumin (transthyretin)** in serum, urine and stool.

Standards, controls and samples containing human **prealbumin (transthyretin)** are added to the wells of a microplate coated with polyclonal anti-human prealbumin antibodies. The antibodies immobilized on the walls of the microtiter wells capture prealbumin in the patient samples during the first incubation step. After washing away the unbound components, a detection peroxidase-conjugated anti-prealbumin antibody is added to each well. During a second incubation, the detection antibody is bound to the captured prealbumin. A "sandwich" of capture antibody - human prealbumin - peroxidase conjugate is formed. Tetramethylbenzidine (TMB) is used as a peroxidase substrate. Finally, an acidic stop solution is added to terminate the reaction. The intensity of the yellow color is directly proportional to the prealbumin concentration of the sample. A dose response curve of the absorbance unit (optical density, OD at 450 nm) vs. concentration is generated using the values obtained from the standards. Prealbumin present in the patient samples is determined directly from this curve.

4. MATERIAL SUPPLIED

Cat. No	Label	Kit Components	Quantity
K 6331MTP	PLATE	One holder with precoated strips	12 x 8 wells
K 6331WP	WASHBUF	ELISA wash buffer concentrate (10x)	2 x 100 ml
K 6331K	CONJ	Conjugate, (rabbit-anti-prealbumin peroxidase-labeled)	1 x 50 μ l
K 6331VP	CONJBUF	Conjugate dilution buffer	1 x 15 ml
K 6331ST	STD	Calibrators, lyophilized (0; 1.6; 6.25; 25; 100 ng/ml)	4 x 5 vials
K 6331KO	CTRL	Control, lyophilized	4 x 1 vial
K 6331TMB	SUB	TMB substrate (Tetramethylbenzidine), ready-to-use	1 x 15 ml
K 6331AC	STOP	ELISA stop solution, ready to use	1 x 15 ml

5. MATERIAL REQUIRED BUT NOT SUPPLIED

- Bidistilled (aqua bidest.) and sterile water
- Laboratory balance
- Precision pipettors calibrated and tips to deliver 5-1000 μ l
- Foil to cover the microtiter plate
- Horizontal microtiter plate shaker
- A multi-channel dispenser or repeating dispenser
- Centrifuge capable of 3000 x g
- Vortex-Mixer
- Standard laboratory glass or plastic vials, cups, etc.
- Microtiter plate reader at 450 or 405 nm
(reference wave length 620 or 690 nm)

6. PREPARATION AND STORAGE OF REAGENTS

- To run assay more than once, ensure that reagents are stored at conditions stated on the label. **Prepare only the appropriate amount necessary for each assay.** The kit can be used up to 4 times within the expiry date stated on the label.
- Reagents with a volume less than **100 µl** should be centrifuged before use to avoid loss of volume.
- The **WASHBUF** (wash buffer concentrate) should be diluted with aqua bidest. **1:10** before use (50 ml WASHBUF + 450 ml aqua bidest.), mix well. Crystals could occur due to high salt concentration in the stock solutions. The crystals must be redissolved at 37°C in a water bath before dilution of the buffer solutions. The **WASHBUF** (wash buffer concentrate) is stable at **2-8°C** until the expiry date stated on the label. Diluted **buffer solution** can be stored in a closed flask at **2-8°C for one month**.
- The lyophilized **STD** (standards) and **CTRL** (control) must be reconstituted with **500 µl** aqua bidest. Allow the vial content to dissolve for 10 minutes and mix thoroughly by gentle inversion to insure complete reconstitution. The lyophilized **STD** (standards) and **CTRL** (control) are stable at **2-8°C** until the expiry date stated on the label. Reconstituted standards and control are not stable.
- The **CONJ** (conjugate) must be diluted **1: 500** in conjugate dilution buffer (CONJBUF) (20 µl CONJ + 10 ml CONJBUF). The undiluted **CONJ** (conjugate) is stable at **2-8 °C** until expiry date stated on the label. **Diluted conjugate is not stable and can not be stored.**
- All other test reagents are ready to use. Test reagents are stable until the expiry date (see label of test package) when stored at **2-8°C**.

7. PRECAUTIONS

- For *in vitro* diagnostic use only.
- Quality control guidelines should be observed.
- Human materials used in kit components were tested and found to be negative for HIV, Hepatitis B and Hepatitis C. However, for safety reasons, all kit components should be treated as potentially infectious.
- Kit reagents contain sodium azide or thimerosal as bactericides. Sodium azide and thimerosal are toxic. Substrates for the enzymatic color reactions are toxic and carcinogenic. Avoid contact with skin or mucous membranes.
- Stop solution is composed of sulfuric acid, which is a strong acid. Even diluted, it still must be handled with care. It can cause acid burns and should be handled with gloves, eye protection, and appropriate protective clothing. Any spills should be wiped out immediately with copious quantities of water.

8. SPECIMEN COLLECTION AND PREPARATION

Faeces

Add about **100 mg** of the sample (note the sample weight for the calculation) to **5 ml** of wash buffer and mix.

Centrifuge the sample suspension for 10 min at 3000 rpm.

Transfer 1 ml of the supernatant into an Eppendorf tube and centrifuge once more at 13.000 rpm for 5 min.

Dilute supernatant **1:40** in wash buffer (25 µl supernatant + 975 µl wash buffer). Use **100 µl** of the solution for the experiment.

We recommend to weight fresh stool samples for each run. Supernatant is not stable and can't be stored.

Stool samples can be stored at -20°C for 4 weeks. Avoid repeated freezing and thawing.

Immundiagnostik recommends the use of tubes from Roche Diagnostics/Mannheim (No. 745804) for sample preparation.

Note: We recommend the determination of the optimal sample dilution (1:500 – 1:5.000) in a preliminary experiment.

Serum samples

Fresh collected blood should be centrifuged within one hour. Store samples at -20 °C if not assayed on the same day. Lipemic or hemolytic samples may give erroneous results. Samples should be mixed well before assaying. We recommend duplicate analyses for each sample.

Normal samples are diluted **1:20,000**.

1:20,000 dilution

The dilution can be performed in three steps. For example:

1. 50 µl serum + 950 µl wash buffer, mix well;
2. of the 1. dilution
100 µl + 900 µl wash buffer, mix well;
3. of the 2. dilution
10 µl + 990 µl wash buffer, mix well (dilution **1:20,000**)

For analysis, pipette in each well **100 µl** of the **final dilution step**.

Note: We recommend the determination of the optimal sample dilution (1:10,000 – 1:30,000) in a preliminary experiment.

Use the corresponding dilution factor to calculate the concentration.

9. ASSAY PROCEDURE

Procedural notes

- Do not interchange different lot numbers of any kit component within the same assay.
- Reagents should not be used beyond the expiration date shown on the kit label.
- Guidelines for medical laboratories should be observed.
- Incubation time, incubation temperature and pipetting volumes of the components are defined by the producer. Any variation of the test procedure, which is not coordinated with the producer, may influence the results of the test. Immundiagnostik AG can therefore not be held responsible for any damage resulting from wrong use.
- Substrate solution should remain colorless until use.
- To ensure accurate results, proper adhesion of plate sealers during incubation steps is necessary.
- Avoid foaming when mixing reagents.
- The assay should always be performed according the enclosed manual.

Test procedure

Wash the precoated PLATE (microtiter plate) **5 x with 250 µl diluted wash buffer**. After the final washing step, the inverted PLATE should be firmly tapped on absorbent paper to remove excess solution.

Carry out the tests in duplicate.


1. Add **100 µl STD** (Standard), **CTRL** (Control) and pre-diluted sample into respective well.
2. Incubate for **1 hour** shaking on a horizontal mixer at room temperature.
3. Decant the content of the PLATE and wash the wells **5 x with 250µl** diluted wash buffer.
4. Add **100 µl** pre-diluted **CONJ** (Peroxidase-labeled antibody).
5. Incubate for **1 hour** shaking on a horizontal mixer at room temperature.
6. Decant the content of the PLATE and wash the wells **5 x with 250µl** diluted wash buffer.
7. Add **100 µl SUB** (TMB substrate).
8. Incubate for **5-15 minutes** at room temperature.
9. Add **50 µl STOP** (stop solution) and mix shortly.
10. Determine **absorption immediately** with an ELISA reader at **450 nm** against 620 nm (or 690 nm) as reference. If no reference wavelength is available, read only at 450 nm. If the extinction of the highest standard exceeds the measurement range of the photometer, absorption must be measured immediately at 405 nm against 620 nm (or 690 nm) as reference.

10. RESULTS

A calibration curve is constructed from the standards. Commercially available software can be used as well as graph paper. Results of the samples are read from this calibration curve.

THE CALIBRATION CURVE IS NOT LINEAR, therefore a spline-algorithm is recommended.

Typical calibration curve

Concentration [ng/ml]	100	25	6.25	1.6	0
OD mean values	2.088	1.145	0.318	0.134	0.081

The data is for demonstration only and cannot be used for the evaluation of test results.

Faeces

To determine the prealbumin concentration of faecal samples, calculate as described in the following example:

Weight: 80 mg (1ml stool = 1g) = 0,08 ml

Dilution step 1: 5ml / 0,08ml = 62,5

Dilution step 2: 40

Dilution factor: 2500

Multiply the result by **2500** to get the real concentration. The dilution factor depends on the weight of the faeces.

Serum samples

For the calculation of the prealbumin concentration in **serum** the result must be multiplied by **20,000**.

11. LIMITATIONS

Samples with prealbumin levels greater than the highest standard value should be further diluted and re-assayed.

12. QUALITY CONTROL

Immundiagnostik recommends commercial control samples for internal quality control.

Control samples should be analyzed with each run. Results, generated from the analysis of control samples, should be evaluated for acceptability using appropriate statistical methods. The results for the patient samples may not be valid, if within the same assay one or more values of the quality control sample are outside the acceptable limits.

Expected values

Each laboratory should evaluate its own baseline value.¹

13. PERFORMANCE CHARACTERISTICS

soon available

14. REFERENCES

1. Schmidt P N, et al. (1995) Scand J Lab Invest 55:35-45

15. GENERAL NOTES ON THE TEST AND TEST PROCEDURE

- This assay was produced and put on the market according to the IVD guidelines of 98/79/EC.
- The test components which are made of human serum are tested for HVB and HIV and found to be negative. However, since no test method can offer complete assurance that infectious agents are absent; these reagents should be handled as recommended for any potentially infectious human serum or blood specimen. The normal precautions for laboratory working should be observed.
- Reagents of the test package contain sodium azide as a bactericide. Contact with skin or mucous membranes has to be avoided.
- All reagents in the test package are to be used for in-vitro diagnostics only.
- The reagents should not be used after the date of expiry (see label on the test package).
- Single components with different lot numbers should not be mixed or exchanged.
- The guidelines for medical laboratories should be observed.
- Incubation time, incubation temperature and pipetting volumes of the different components have been defined by the producer. Any alterations of the test procedure, that are not coordinated with the producer, may influence the results of the test. Immundiagnostik can therefore not be held responsible for any damage.

Used symbols:

Temperature limitation

Catalogue Number

In Vitro Diagnostic Medical Device

Contains sufficient for <n> tests

Manufacturer

Use by

Lot number