# slgA ELISA Kit

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
***For the in vitro determination of Secretory IgA in saliva and stool***

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| <b>TABLE OF CONTENT</b>                                 | <b>PAGE</b> |
|---|-------------|
| <b>1. INTENDED USE</b>                                  | <b>19</b>   |
| <b>2. INTRODUCTION</b>                                  | <b>19</b>   |
| <b>3. PRINCIPLE OF THE TEST</b>                         | <b>19</b>   |
| <b>4. MATERIAL SUPPLIED</b>                             | <b>20</b>   |
| <b>5. MATERIAL REQUIRED BUT NOT SUPPLIED</b>            | <b>20</b>   |
| <b>6. PREPARATION AND STORAGE OF REAGENTS</b>           | <b>21</b>   |
| <b>7. PRECAUTIONS</b>                                   | <b>21</b>   |
| <b>8. SPECIMEN COLLECTION AND PREPARATION</b>           | <b>22</b>   |
| EXTRACTION OF THE STOOL SAMPLE                          | 22          |
| DILUTION OF SAMPLES                                     | 24          |
| <b>9. ASSAY PROCEDURE</b>                               | <b>24</b>   |
| PROCEDURAL NOTES  | 24          |
| TEST PROCEDURE  | 25          |
| <b>10. RESULTS</b>                                      | <b>26</b>   |
| <b>11. LIMITATIONS</b>                                  | <b>27</b>   |
| <b>12. QUALITY CONTROL</b>                              | <b>27</b>   |
| EXPECTED VALUES   | 27          |
| <b>13. PERFORMANCE CHARACTERISTICS</b>                  | <b>28</b>   |
| PRECISION AND REPRODUCIBILITY                           | 28          |
| RECOVERY  | 29          |
| SENSITIVITY   | 29          |
| CROSS REACTIVITY  | 29          |
| SAMPLE DILUTION   | 30          |
| <b>14. REFERENCES</b>                                   | <b>30</b>   |
| <b>15. GENERAL NOTES ON THE TEST AND TEST PROCEDURE</b> | <b>32</b>   |

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Wiederfindung n=2

| Probe [ng/ml] | Spike [ng/ml] | slgA Erwartet [ng/ml] | slgA Gemessen [ng/ml] |
|---------------|---------------|-----------------------|-----------------------|
| 103.7         | 150           | 253.7                 | 279.7                 |
| 103.7         | 75            | 178.7                 | 194.7                 |
| 103.7         | 50            | 153.7                 | 158.7                 |
| 103.7         | 25            | 128.7                 | 141.2                 |
| 100.3         | 150           | 250.3                 | 272.4                 |
| 100.3         | 75            | 175.3                 | 212.9                 |
| 100.3         | 50            | 150.3                 | 165.4                 |
| 100.3         | 25            | 125.3                 | 126.5                 |

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Brandtzaeg P, Johansen FE (2005) Mucosal B cells: phenotypic characteristics, transcriptional regulation, and homing properties. *Immunol Rev.* Aug; **206**: 32-63. Review

Burnett D, Crocker J, Stockley RA. (1987) Cells containing IgA subclasses in bronchi of subjects with and without chronic obstructive lung disease. *J Clin Pathol.* Oct; **40(10)**: 1217-20

- Corthésy B and Spertini F (1999) Secretory immunoglobulin A: from mucosal protection to vaccine development. *Biol. Chem. Nov*; **380**: 1251–1262
- Hein M, Petersen A C, Helmig Rb, Uldbjerg N, Reinholdt J (2005) Immunoglobulin levels and phagocytes in the cervical mucus plug at term of pregnancy. *Acta Obstet Gynecol Scand.* **84**: 734–742
- Hocini H, Iscaki S, Bouvet JP, Kazatchkine MD, Bélec L (2000) An ELISA method to measure total and specific human secretory IgA subclasses based on selective degradation by IgA1-protease. *J Immunol Methods.* Feb 21; **235(1-2)**: 53-60
- Kitz R, Ahrens P, Zielen S (2000) Immunoglobulin levels in bronchoalveolar lavage fluid of children with chronic chest disease. *Pediatric Pulmonology* **29**:443–451
- Könönen E, Jousimies-Somer H, Bryk A, Kilpi T, Kilian M (2002) Establishment of streptococci in the upper respiratory tract: longitudinal changes in the mouth and nasopharynx up to 2 years of age. *J. Med. Microbiol.* **51**: 23-730
- Majkowska-Skrobek G, Augustyniak D, Jankowski A (2003) Assessment of IgA subclasses synthesis in children with selective and partial IgA deficiency. *Centr Eur J Immunol* **28(3)**: 110–118
- Mestecky J, Russell MW, Elson CO (1999) Intestinal IgA: novel views on its function in the defence of the largest mucosal surface intestinal. *Gut* **44**: 2-5
- Rose MA, Schubert R, Schmitt-Grohe S, Reichenbach J, Zielen S (2006) Immunoglobulins and inflammatory cytokines in nasal secretions in humoral immunodeficiencies. *Laryngoscope* **116**: 239–244
- Rüssmann H, Lissner R, Schmidt H, Karch H (1999) IgA/IgM and secretory immunity. *Sepsis* **3**: 219–224

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Manual


# slgA ELISA Kit

*For the in vitro determination of Secretory IgA in saliva and stool*

Valid from 30.03.2010

K 8870

## 1. INTENDED USE

The *Immundiagnostik* Assay is intended for the quantitative determination of **secretory IgA (sIgA)** in saliva and stool. For *in vitro* diagnostic use only.

## 2. INTRODUCTION

Secretory IgA (sIgA) consists of two IgA monomers joined by the J-chain and an additional secretory component. It is secreted in plasma cells located in the lamina propria of mucosal membranes. Synthesis of sIgA is independent from the synthesis of serum IgA. This means that lack of serum IgA does not necessarily correlate with a lack of sIgA. Secretory IgA is the major immunoglobulin in saliva, tears, colostrum, nasal mucous, mother's milk, tracheobronchial and gastrointestinal secretes. It plays a major role in preventing adherence of microorganisms to mucosal sites, in activation of the alternative complement pathway and in activating inflammatory reactions. Newborns are provided with sIgA by mother's milk and are passively immunized against gastrointestinal infections.

### Indications

- Proof of an imbalanced immunological barrier on the intestinal mucosa
- Autoimmune disease

## 3. PRINCIPLE OF THE TEST

This Enzyme-Linked-Immuno-Sorbent-Assay (ELISA) is intended for the quantitative determination of **secretory IgA** in stool and saliva. In a first incubation step, the sIgA in the samples is bound to polyclonal antibodies (rabbit anti human IgA), which are immobilized to the surface of the microtiter wells. To remove all unbound substances, a washing step is carried out. In a second incubation step, a Peroxidase-labeled conjugate (mouse anti-sIgA) is added which recognizes specifically the bound secretory IgA. After another washing step, to remove all unbound substances, the solid phase is incubated with the substrate, Tetramethylbenzidine (TMB). An acidic stop solution is then added to stop the reaction. The color converts from blue to yellow. The intensity of the yellow color is directly proportional to the concentration of secretory IgA. A dose response curve of the absorbance unit (optical density, OD) vs. concentration is generated, using the results obtained from the calibrators. Secretory IgA in the patient samples is determined directly from this curve.

#### 4. MATERIAL SUPPLIED

| Cat. No   | Content | Kit Components  | Quantity     |
|-----------|---------|---|--------------|
| K 8870MTP | PLATE   | One holder with precoated strips                                | 12 x 8 wells |
| K 8870WP  | WASHBUF | ELISA wash buffer concentrate 10x                               | 1 x 100 ml   |
| K 8870K   | CONJ    | Conjugate (mouse anti-slgA,<br>Peroxidase-labeled)              | 1 x 200 µl   |
| K 8870ST  | STD     | Calibrators, lyophilized<br><br>(0; 22.2; 66.6; 200; 600 ng/ml) | 2 x 5 vials  |
| K 8870KO1 | CTRL    | Control, lyophilized  | 2 x 1 vial   |
| K 8870KO2 | CTRL    | Control, lyophilized  | 2 x 1 vial   |
| K 8870TMB | SUB     | TMB substrate (Tetramethylbenzidine),<br>ready-to-use           | 1 x 15 ml    |
| K 8870AC  | STOP    | ELISA stop solution, ready-to-use                               | 1 x 15 ml    |

#### 5. MATERIAL REQUIRED BUT NOT SUPPLIED

- Bidistilled water (aqua bidest.)
- Laboratory balance
- Precision pipettors and disposable tips to deliver 10-1000 µl
- Foil to cover the microtiter plate
- Horizontal microtiter plate shaker
- A multi-channel dispenser or repeating dispenser
- Centrifuge capable of 3000 x g
- Vortex-Mixer
- Standard laboratory glass or plastic vials, cups, etc.
- Microtiter plate reader at 450 nm  
(reference wave length 620 or 690 nm)

## 6. PREPARATION AND STORAGE OF REAGENTS

- To run assay more than once, ensure that reagents are stored at conditions stated on the label. **Prepare only the appropriate amount necessary for each assay.** The kit can be used up to 4 times within the expiry date stated on the label.
- The **ELISA WASHBUF** (wash buffer concentrate) should be diluted with aqua dest. **1:10** before use (100 ml WASHBUF + 900 ml aqua bidest.). Crystals could occur due to high salt concentration. The crystals must be resuspended **before dilution of the buffer solutions** using a water bath (37°C). The buffer concentrate is stable at 2-8°C until the expiry date stated on the label. Diluted solutions can be stored at 2-8°C for 1 month.
- The **STD** (standards) and **CTRL** (control) must be reconstituted with **500 µl aqua bidest.** Allow the vial content to dissolve for **10 minutes** and mix thoroughly by gentle inversion to insure complete reconstitution. **Reconstituted standards and control are stable at -20 °C** until the expiry date stated on the label and can be subjected to a **maximum of two freeze-thaw cycles.**
- The **CONJ** (conjugate; POD-labeled antibody) must be diluted **1:101** in wash buffer (100 µl CONJ + 10 ml wash buffer). The undiluted conjugate is stable at 2-8 °C until the expiry date stated on the label. Diluted conjugate is not stable and can not be stored.
- All other test reagents are ready for use. The test reagents are stable up to the date of expiry (see label of test package) when stored at **2 – 8 °C.**

## 7. PRECAUTIONS

- For *in vitro* diagnostic use only.
- Human materials used in kit components were tested and found to be negative for HIV, Hepatitis B and Hepatitis C. However, for safety reasons, all kit components should be treated as potentially infectious.
- Stop solution is composed of sulfuric acid, which is a strong acid. Even diluted, it still must be handled with care. It can cause acid burns and should be handled with gloves, eye protection, and appropriate protective clothing. Any spills should be wiped out immediately with copious quantities of water.
- Reagents should not be used beyond the expiration date shown on the kit label.

## 8. SPECIMEN COLLECTION AND PREPARATION

### Saliva

To avoid variation in slgA content, take saliva samples always at the same time of the day.

No food or liquid should be consumed 30 min before sample collection. Centrifuge the samples at 3000 rpm for 10 min. Sample supernatant can be stored at -20°C.

For analysis, the supernatant is diluted **1:2000** in ELISA wash buffer, e.g.

10 µl supernatant + 990 µl wash buffer; dilute the obtained solution again:  
50 µl diluted supernatant + 950 µl wash buffer.

Use **100 µl** of the final dilution per well.

### Faeces

#### *Extraction of the stool sample*

The **wash buffer** is used as a sample extraction buffer. We recommend the following sample preparation:

1. We recommend the use of a stool sample preparation kit for dosing 100 mg of stool sample (e.g. Sample preparation kit from Roche Diagnostics, Mannheim, Germany; cat # 745804). The stool sample has to be suspended in 5 ml buffer.

**Constant buffer volume: 5 ml**

**Constant dilution factor: 1:50**

2. Alternatively, stool samples can be manually weighted within a range of 80 – 120 mg. Please note the exact sample amount!
  - a. Add **5 ml** buffer to the stool sample independent from the sample amount.

**Constant buffer volume: 5 ml**

**The dilution factor varies** depending on the sample amount which must be considered in the subsequent calculations as shown below:

| Weight [mg] | Dilution factor |
|-------------|-----------------|
| 80          | 62.5            |
| 82          | 60.9            |
| 84          | 59.5            |
| 86          | 58.1            |
| 88          | 56.8            |
| 90          | 55.6            |
| 92          | 54.3            |
| 94          | 53.2            |
| 96          | 52.1            |
| 98          | 51.0            |
| <b>100</b>  | <b>50</b>       |

| Weight [mg] | Dilution factor |
|-------------|-----------------|
| 102         | 49.0            |
| 104         | 48.1            |
| 106         | 47.2            |
| 108         | 46.3            |
| 110         | 45.5            |
| 112         | 44.6            |
| 114         | 43.9            |
| 116         | 43.1            |
| 118         | 42.4            |
| 120         | 41.6            |

- b. **The buffer volume** for the individual samples varies depending on the sample amount (see table). The dilution factor remains constant.

#### **Variable buffer volume**

##### **Constant dilution factor: 1:50**

Therefore, the same dilution factor can be used for all samples in the subsequent evaluation of the results.

| Weight [mg] | Buffer Volume [ml] |
|-------------|--------------------|
| 80          | 4.0                |
| 82          | 4.1                |
| 84          | 4.2                |
| 86          | 4.3                |
| 88          | 4.4                |
| 90          | 4.5                |
| 92          | 4.6                |
| 94          | 4.7                |
| 96          | 4.8                |
| 98          | 4.9                |
| <b>100</b>  | <b>5.0</b>         |

| Weight [mg] | Buffer Volume [ml] |
|-------------|--------------------|
| 102         | 5.1                |
| 104         | 5.2                |
| 106         | 5.3                |
| 108         | 5.4                |
| 110         | 5.5                |
| 112         | 5.6                |
| 114         | 5.7                |
| 116         | 5.8                |
| 118         | 5.9                |
| 120         | 6.0                |

Afterwards, mix stool sample and buffer; vortex for at least 30 sec. depending on the stool consistency.

Transfer approx. 1 ml stool suspension to an Eppendorf-tube and centrifuge for 5 minutes at 13000 rpm ( $\cong$  13000 g).

**The supernatant is not stable and can not be stored.** We recommend to weight fresh sample amount for a new assay, if the analysis should be repeated.

### *Dilution of samples*

#### **Stool samples**

After centrifugation, the supernatant is diluted **1:250 in wash buffer**.

For example:

**40 $\mu$ l** supernatant + **960  $\mu$ l** wash buffer, mix well (**dilution I**) (**1:25**)

**100  $\mu$ l** of this dilution I + **900 $\mu$ l** wash buffer, mix well (**dilution II**) (**1:10**)

For analysis, pipette **100  $\mu$ l** of **dilution step II** solution per well.

## **9. ASSAY PROCEDURE**

### *Procedural notes*

- Do not mix different lot numbers of any kit component.
- Quality control guidelines should be followed.
- Incubation time, incubation temperature and pipetting volumes of the components are defined by the producer. Any variation of the test procedure, which is not coordinated with the producer, may influence the results of the test. Immundiagnostik AG can therefore not be held responsible for any damage resulting from wrong use.
- The assay should always be performed according the enclosed manual.

## Test procedure

Wash the precoated PLATE (microtiter plate) 5 x with 250 µl ELISA wash buffer. Carry out the tests in duplicate.

1. Add **100 µl** STD (standards), CTRL (controls) and patient samples (faeces and saliva diluted, see above).
2. Incubate for **1 hour**, shaking on a horizontal mixer, at room temperature.
3. Aspirate and wash the wells **5 x with 250 µl** ELISA wash buffer.
4. Add **100 µl** CONJ (conjugate; POD antibody).
5. Incubate for **1 hour**, shaking on a horizontal mixer, at room temperature.
6. Decant the content of the plate and wash the wells **5 x with 250 µl** wash buffer.
7. Add **100 µl** SUB (TMB substrate).
8. Incubate for **10-20 minutes** at room temperature.
9. Add **50 µl** STOP (ELISA stop solution) and mix shortly.
10. Determine absorption with an ELISA reader at **450 nm** against 620 nm as reference. If no reference wavelength is available, read only at 450 nm. If the extinction of the highest standard exceeds the measurement range of the photometer, absorption must be measured immediately at 405 nm against 620 nm as reference.

## 10. RESULTS

The following algorithms can be used alternatively to calculate the results. We recommend to use the "4-Parameter-algorithm".

1. 4-parameter-algorithm

It is recommended to use a linear ordinate for optical density and a logarithmic abscissa for concentration. When using a logarithmic abscissa, the zero calibrator must be specified with a value less than 1 (e. g. 0.001).

2. Point-to-point-calculation

We recommend a linear ordinate for optical density and a linear abscissa for concentration.

3. Spline-algorithm

We recommend a linear ordinate for optical density and a logarithmic abscissa for concentration. When using a logarithmic abscissa, the zero calibrator must be specified with a value less than 1 (e. g. 0.001).

The plausibility of the pairs of values should be examined before the automatic evaluation of the results. If this option is not available with the used program, a control of the paired values should be done manually.

### Saliva

For the calculation of the saliva values, the results from the microplate reader must be multiplied by **2.000**.

### Faeces

1. A **final dilution factor of 1:12500** results for a constant dilution factor 50 in dilution step 1:

Dilution step 1: 50

Dilution step 2: 250

Final dilution factor:  $50 \times 250 = 12500$

2. A **final dilution factor of 1: 15625** results for a sample weight of 80 mg and a buffer volume of 5ml in dilution step 1:

Weight: 80 mg (1ml stool = 1g) = 0,08 ml

Dilution step 1:  $5\text{ml} / 0,08\text{ml} = 62,5$

Dilution step 2: 250

Final dilution factor:  $62,5 \times 250 = 15625$

Multiply the result by **15625** to obtain the sample concentration. The dilution factor depends on the weight of the faeces.

## 11. LIMITATIONS

Samples with sIgA levels greater than the highest calibrator, should be diluted and re-assayed.

## 12. QUALITY CONTROL

Control samples should be analyzed with each run. Results, generated from the analysis of control samples, should be evaluated for acceptability using appropriate statistical methods. The results for the patient samples may not be valid, if within the same assay one or more values of the quality control sample are outside the acceptable limits.

### *Expected values*

#### **Normal range**

##### ***Secretory IgA in saliva***

|                      |                                  |
|----------------------|----------------------------------|
| Children (n=37)      | 18 - 237 µg/ml (mean 128 µg/ml)* |
| Age >16 years (n=33) | 102 - 471 µg/ml                  |

***Secretory IgA in faeces*** 510 - 2040 µg/ml (n = 76)

It is recommended for each laboratory to establish its own normal range.

\*Hofman LF, Le T (2002) Preliminary pediatric reference range for secretory IgA in saliva using an enzyme immunoassay. *Clinical Chemistry* **48 (6)**:A169, Suppl.

### 13. PERFORMANCE CHARACTERISTICS

#### *Precision and reproducibility*

The precision (intra-assay variation) of the Immundiagnostik anti-slgA ELISA test was calculated from 20 replicate determinations on each of one samples.

Intra-Assay CV n= 20

| Sample | slgA [ng/ml] | Intra-Assay CV [%] |
|--------|--------------|--------------------|
| 1      | 77.7         | 5                  |
| 2      | 92.5         | 9                  |

The total precision (inter-assay variation) of the Immundiagnostik anti-slgA ELISA test was calculated from data on 2 samples obtained in 20 different assays by three technicians on two different lots of reagents over a period of three months.

Inter-Assay CV n= 20

| Sample | slgA [ng/ml] | Inter-Assay CV [%] |
|--------|--------------|--------------------|
| 1      | 102.4        | 8                  |
| 2      | 1277.4       | 7.4                |

### Recovery

Two samples were spiked with slgA calibrator and measured with this assay.

Recovery n=2

| Sample [ng/ml] | Spike [ng/ml] | slgA Expected [ng/ml] | slgA Measured [ng/ml] |
|----------------|---------------|-----------------------|-----------------------|
| 103.7          | 150           | 253.7                 | 279.7                 |
| 103.7          | 75            | 178.7                 | 194.7                 |
| 103.7          | 50            | 153.7                 | 158.7                 |
| 103.7          | 25            | 128.7                 | 141.2                 |
| 100.3          | 150           | 250.3                 | 272.4                 |
| 100.3          | 75            | 175.3                 | 212.9                 |
| 100.3          | 50            | 150.3                 | 165.4                 |
| 100.3          | 25            | 125.3                 | 126.5                 |

### Sensitivity

n=20

| Sample | slgA Mean value [OD] | Standard variation | Detection limit [ng/ml] |
|--------|----------------------|--------------------|-------------------------|
| 1      | 0.171                | 0.099              | 13.4                    |

### Cross reactivity

No cross reactivity to other proteins in stool and saliva.

## Sample dilution

Linearity n= 2

Two patient samples were diluted with ELISA wash buffer. The results are shown below:

| Sample | Dilution  | Expected [ng/ml] | Measured [ng/ml] |
|--------|-----------|------------------|------------------|
| A      | undiluted | 126.8            | 126.8            |
|        | 1:2       | 63.4             | 65.5             |
|        | 1:4       | 31.7             | 35.1             |
|        | 1:8       | 15.9             | 25.6             |
| B      | undiluted | 184.9            | 184.9            |
|        | 1:2       | 92.5             | 93.7             |
|        | 1:4       | 46.2             | 52.1             |
|        | 1:8       | 23.1             | 21.9             |

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Brandtzaeg P, Johansen FE (2005) Mucosal B cells: phenotypic characteristics, transcriptional regulation, and homing properties. *Immunol Rev.* Aug; **206**: 32-63. Review

Burnett D, Crocker J, Stockley RA. (1987) Cells containing IgA subclasses in bronchi of subjects with and without chronic obstructive lung disease. *J Clin Pathol.* Oct; **40(10)**: 1217-20

Corthésy B and Spertini F (1999) Secretory immunoglobulin A: from mucosal protection to vaccine development. *Biol. Chem.* Nov; **380**: 1251–1262

- Hein M, Petersen A C, Helmig Rb, Uldbjerg N, Reinholdt J (2005) Immunoglobulin levels and phagocytes in the cervical mucus plug at term of pregnancy. *Acta Obstet Gynecol Scand.* **84**: 734–742
- Hocini H, Iscaki S, Bouvet JP, Kazatchkine MD, Bélec L (2000) An ELISA method to measure total and specific human secretory IgA subclasses based on selective degradation by IgA1-protease. *J Immunol Methods.* Feb 21; **235(1-2)**: 53-60
- Kitz R, Ahrens P, Zielen S (2000) Immunoglobulin levels in bronchoalveolar lavage fluid of children with chronic chest disease. *Pediatric Pulmonology* **29**:443–451
- Könönen E, Jousimies-Somer H, Bryk A, Kilpi T, Kilian M (2002) Establishment of streptococci in the upper respiratory tract: longitudinal changes in the mouth and nasopharynx up to 2 years of age. *J. Med. Microbiol.* **51**: 23-730
- Majkowska-Skrobek G, Augustyniak D, Jankowski A (2003) Assessment of IgA subclasses synthesis in children with selective and partial IgA deficiency. *Centr Eur J Immunol* **28(3)**: 110–118
- Mestecky J, Russell MW, Elson CO (1999) Intestinal IgA: novel views on its function in the defence of the largest mucosal surface intestinal. *Gut* **44**: 2-5
- Rose MA, Schubert R, Schmitt-Grohe S, Reichenbach J, Zielen S (2006) Immunoglobulins and inflammatory cytokines in nasal secretions in humoral immunodeficiencies. *Laryngoscope* **116**: 239–244
- Rüssmann H, Lissner R, Schmidt H, Karch H (1999) IgA/IgM and secretory immunity. *Sepsis* **3**: 219–224
- Publications based on Immundiagnostik's slgA ELISA**
- Hofman LF, Le T (2002) Preliminary pediatric reference range for secretory IgA in saliva using an enzyme immunoassay. *Clin Chem Vol* **48(6)** Suppl A169
- Noel N, Levinson U, Ho J, Lee S (2001) Measurement of secretory IgA in fecal samples preserved with PBS/0.02% thimerosal. Abstract P 155 *Clin Chem Vol* **47 (6)**, Suppl A47
- Martin M (1999) Diagnostik bei latent-entzündlichen Darmschleimhautveränderungen. *PRAXIS-telegramm* **1**: 17-21
- Martin M (1998) Diagnostik chronisch entzündlicher Erkrankungen des Darms. *PRAXIS-telegramm* **6**: 16-18
- Martin M (1998) Diagnostik zivilisationsbedingter Erkrankungen im Kindesalter. *PRAXIS-telegramm* **3-4**: 15-17
- Schütz B, Poschwatta-Rupp S, Rusch K, Zimmermann K (1998) Das KyberPlus-Konzept. Diagnostische Verfahren zur Abklärung unklarer intestinaler Beschwerden. *Biologische Medizin* **27(1)**: 31-36

## 15. GENERAL NOTES ON THE TEST AND TEST PROCEDURE

- This assay was produced and put on the market according to the IVD guidelines of 98/79/EC.
- For *in vitro* diagnostic use only.
- Quality control guidelines should be followed.
- Human materials used in kit components were tested and found to be negative for HIV, Hepatitis B and Hepatitis C. However, for safety reasons, all kit components should be treated as potentially infectious.
- Kit reagents contain sodium azide or thimerosal as bactericides. Sodium azide and thimerosal are toxic. Substrates for the enzymatic color reactions are toxic and carcinogenic. Avoid contact with skin or mucous membranes.
- Stop solution is composed of sulfuric acid, which is a strong acid. Even diluted, it still must be handled with care. It can cause acid burns and should be handled with gloves, eye protection, and appropriate protective clothing. Any spills should be wiped out immediately with copious quantities of water.
- Do not mix different lot numbers of any kit component.
- Reagents should not be used beyond the expiration date shown on the kit label.
- The assay should always be performed according the enclosed manual.
- Incubation time, incubation temperature and pipetting volumes of the components are defined by the producer. Any variation of the test procedure, which is not coordinated with the producer, may influence the results of the test. Immundiagnostik AG can therefore not be held responsible for any damage resulting from wrong use.
- Warranty claims and complaints in respect of deficiencies must be logged within 14 days after receipt of the product. The product should be send to Immundiagnostik AG along with a written complaint.

**Used symbols:**

Temperature limitation

Catalogue Number

In Vitro Diagnostic Medical Device

Contains sufficient for <n> tests

Manufacturer

Use by

Lot number