

YK140 Rat GLP-2 EIA Kit

I. Introduction

The proglucagon gene is expressed in both pancreatic A cell and intestinal L cell. Tissue-specific posttranslational processing of proglucagon by the prohormone convertase produced the different proglucagon derived peptides (PGDPs) in both pancreas and intestine. The most notable pancreatic PGDP is glucagon, whereas the L cell produces several structurally related peptides, including glucagon-like peptide 1 (GLP-1) and GLP-2, as well as glicentin and oxyntomodulin, which contain glucagon sequence in their molecules. Among PGDPs, GLP-2 has recently been found to show intestinal epithelial proliferation.

YK140 Rat GLP-2 EIA Kit	Contents
The assay kit can measure GLP-2 in the range of 0.137 - 100 ng/mL	1) Antibody coated plate
The assay completes within 16-18 hr. +1.5 hr.	2) Rat GLP-2 standard
With one assay kit, 40 samples can be measured in duplicate	3) Labeled antigen
Test sample: rat serum or plasma	4) Rat GLP-2 antibody
Sample volume: 25 μ L	5) SA-HRP solution
The 96-well plate in kit was consisted by 8-wells strips. The kit can be used separately.	6) Substrate buffer
Precision and reproducibility	7) OPD tablet
Intra-assay CV (%) serum 3.5 - 8.9	8) Stopping solution
Inter-assay CV (%) serum 7.6 - 13.0	9) Buffer solution
Intra-assay CV (%) plasma 3.1 - 7.2	10) Washing solution (concentrated)
Inter-assay CV (%) plasma 6.7 - 11.5	11) Adhesive foil
Stability and Storage	
Store all of the components at 2-8 .	
12 months from the date of manufacturing.	
The expiry date is described on the label of kit.	

II. Characteristics

. Characteristics

This EIA kit is used for quantitative determination of rat GLP-2 in serum or plasma samples. The kit is characterized for sensitive quantification, high specificity and no influences with other components in samples. Rat GLP-2 standard is highly purified synthetic product.

< Specificity >


The EIA kit has high specificity to rat GLP-2 and shows no cross reactivity with rat glucagon and rat GLP-1 even in the concentration of 300 pmol/mL.

< Test Principle >

This EIA kit for determination of rat GLP-2 in serum or plasma samples is based on a competitive enzyme immunoassay using combination of highly specific antibody to rat GLP-2 and biotin-avidin affinity system. The 96-wells plate is coated with goat anti rabbit IgG antibody. Rat GLP-2 standard or samples, labeled antigen and anti rat GLP-2 polyclonal antibody are added to the wells for competitive immunoreaction. After incubation and plate washing, HRP labeled streptoavidin (SA-HRP) are added to form HRP labeled streptoavidin-biotinylated rat GLP-2-antibody complex on the surface of the wells. Finally, HRP enzyme activity is determined by o-Phenylenediamine dihydrochloride (OPD) and the concentration of rat GLP-2 is calculated.

III. Performance Characteristics

Typical standard curve



Analytical recovery

< Rat serum >

Sample No.	Rat GLP-2 added (ng/mL)	Observed (ng/mL)	Expected (ng/mL)	Recovery (%)
1	0	2.55	-	-
2	1	3.59	3.29	109.0
3	5	8.11	7.29	111.2
4	20	27.08	22.29	122.8

< Rat plasma >

Sample No.	Rat GLP-2 added (ng/mL)	Observed (ng/mL)	Expected (ng/mL)	Recovery (%)
1	0	2.90	-	-
2	1	3.89	3.61	107.8
3	5	8.28	7.61	108.7
4	20	26.85	22.61	118.7

Precision and reproducibility

- Intra-assay/Rat serum CV (%) 3.5 ~ 8.9
- Inter-assay/Rat serum CV (%) 7.6 ~ 13.0
- Intra-assay/Rat plasma CV (%) 3.1 ~ 7.2
- Inter-assay/Rat plasma CV (%) 6.7 ~ 11.5

Assay range

0.137 ~ 100 ng/mL

IV. Stability and Storage

- < Storage > Store all of the components at 2-8 .
- < Shelf life > 12 months from the date of manufacturing
 The expiry date is described on the label of kit.
- < Package > For 96 tests per one kit including standards

V. References

1. Philippe, J.: Structure and pancreatic expression of the insulin and glucagon gene. *Endocr Rev* **12**: 252 - 271, 1991
2. Mojsov S. et al: Preproglucagon gene expression in pancreas and intestine diversifies the level of post-transcriptional processing. *J Biol Chem* **261**: pp11880 – 11889, 1986
3. Drucker, D. J. et al: Induction of intestinal epithelial proliferation by glucagon-like peptide 2. *Proc Natl Acad Sci USA* **93**: 7911 – 7916, 1996

<Manufacturer>

Yanaihara Institute Inc.

2480-1 Awakura, Fujinomiya-shi

Shizuoka, Japan 418-0011

TEL: +81-544-22-2771 FAX: +81-544-22-2770

Website: <http://www.yanaihara.co.jp> E-mail: ask@yanaihara.co.jp

Update at July 13, 2005